

Oct 27, 2020

🌐 Protein expression of hard-to-produce proteins in the cytoplasm of *Escherichia coli*

📁 In 1 collection

DOI


dx.doi.org/10.17504/protocols.io.bfbwjipe

Cristina Hernandez Rollan¹, Kristoffer Bach Falkenberg¹, Maja Rennig¹, Andreas Birk Bertelsen¹, Morten Norholm¹

¹Technical University of Denmark



Cristina Hernandez Rollan



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bfbwjipe>

Protocol Citation: Cristina Hernandez Rollan, Kristoffer Bach Falkenberg, Maja Rennig, Andreas Birk Bertelsen, Morten Norholm 2020. Protein expression of hard-to-produce proteins in the cytoplasm of *Escherichia coli*. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bfbwjipe>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: April 20, 2020

Last Modified: October 28, 2020

Protocol Integer ID: 35926

Keywords: E. coli, T7, Protein production, Cytoplasmic expression, sHuffle strain, disulfide bond, lytic polysaccharide monoxygenase, extraction of cytoplasmatic protein, cytoplasmatic protein, proteins in the cytoplasm, extraction of protein, protein expression, expressed protein, protein, cytoplasm, chemical lysis protocol, bacteria, based expression system, expression system, lpmo

Abstract

This protocol describes the expression and extraction of heterologously cytoplasmically expressed proteins in the bacterium *E.coli*, for T7 based expression systems.

Specifically, it describes how to produce "hard-to-produce", as for example, LPMOs.

It is based on a modified version of:

"[IPTG Induction and Extraction of Proteins from Bacteria](#)" by Swathi Arur and Sudhir Nayak, Schedl Lab. Washington University Genetics, St. Louis.

Hemsworth, G. R., Henrissat, B., Davies, G. J., and Walton, P. H. (2014) Discovery and characterization of a new family of lytic polysaccharide monoxygenases. Nat. Chem. Biol.10, 122–126.

The extraction of cytoplasmatic proteins is followed by a chemical lysis protocol, alternatively, sonication can also be employed.

Guidelines

Make sure to clone your favorite gene for expression in the appropriate expression vector.

Materials

MATERIALS

⊗ sHuffle T7 Competent E.coli - 6×0.05 ml **New England Biolabs Catalog #C3026H**

⊗ IPTG **Bio Basic Inc. Catalog #IB0168.SIZE.100g**

⊗ BL21(DE3) or BL21-Star(DE3) or Rosetta2(DE3) or etc for protein purification

⊗ Kanamycin **Research Products International Corp (RPI) Catalog #K22000-25.0**

⊗ LB **Research Products International Corp (RPI) Catalog #L24400-2000.0**

Troubleshooting



Safety warnings

! This protocol describes the construction of GMO classified organisms. Make sure that the local GMO and safety legislations are respected.

Before start

Prepare a fresh transformation of your expression vector in *E. coli* BL21 (DE3) cells or other appropriate strain (e.g. sHuffle strain, Rosetta, etc.)



Day 1 - Culture inoculation

- 1 Pick a fresh colony of your BL21 (DE3) strain previously transformed with your expression vector, and inoculate it in LB medium supplemented with relevant antibiotics. Grow the culture at 37 °C at 250 RPM shaking Overnight . The volume of the overnight culture depends on the volume of the expression culture and should be at least 1/100 of the expression culture

Day 2 - Induction

- 2 Dilute the overnight culture 1:100 in fresh LB supplemented with relevant antibiotics
- 3 Grow the culture at 37 °C with 250 RPM shaking until an OD₆₀₀ = 0.5 - 0.6
- 4 Move the culture into an incubator set to 18 °C with 180 RPM of shaking and allow the culture to reach OD₆₀₀ = 0.8 - 1.0
- 5 Induce the expression by adding IPTG to a final concentration of 1 millimolar (mM)
- 6 Grow the culture for 20:00:00 at 18 °C with 180 RPM shaking.

Note

The temperature can be lowered further allowing for a slower induction in order to increase the production of soluble product.

Day 3 - Collection and lysis

- 7 Collect induced samples by centrifugation at 8000 x g, 4°C, 00:20:00 and carefully discard the supernatant



Note

Pellets can be frozen here at -80 °C or -20 °C until protein extraction

- 8 Lyse the cells using a lysis protocol or a commercially available kit
- 9 Keep the extractions On ice when working with it and at 4 °C for storage