

May 07, 2019

## Protein Expression

DOI

[dx.doi.org/10.17504/protocols.io.sj2ecqe](https://dx.doi.org/10.17504/protocols.io.sj2ecqe)



Rengin Yamur Akbiyik<sup>1</sup>, Sevval Uysalcan<sup>1</sup>

<sup>1</sup>Istanbul Bilgi University



Rengin Yamur Akbiyik

Istanbul Bilgi University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.sj2ecqe](https://dx.doi.org/10.17504/protocols.io.sj2ecqe)

**Protocol Citation:** Rengin Yamur Akbiyik, Sevval Uysalcan 2019. Protein Expression. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.sj2ecqe>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 13, 2018

**Last Modified:** May 07, 2019

**Protocol Integer ID:** 14682

**Keywords:** Protein Expression



## Abstract

The protein expression is a process about synthesizing and modification for proteins in living organisms with using several laboratory techniques and devices. This technique is generally used for obtaining protein from the competent cell that is used for carrying plasmids that possess the protein that we want to express.

### **Protein expression involves two general steps that;**

- Gene expression; which is continuum that converting DNA to the mRNA that can be translate into protein.
- Protein manufacturing; which is translation process.

## Materials

### STEP MATERIALS

 LB Broth **Amresco Catalog #J106-2KG**

## Protocol materials

 LB Broth **Amresco Catalog #J106-2KG**

 LB Broth **Amresco Catalog #J106-2KG**



Transform BL21 cells and plate them on Amp (Kan).

1

Afer overnight incubation, take colonies and transfer into 6 ml LB with 6 µl Amp. Let them grow for 3-5 hours.

2

🌡️ 37 °C

🧪 LB Broth Amresco Catalog #J106-2KG

🧪 6 mL LB

🧪 6 µL Ampicillin

When incubation process is done, transfer it to 1 L ZYM-5052 Medium with 1 µl Amp.

3

Centrifuge at 6000 rpm for 10 minutes.

4

#### Equipment

SORVALL LYNX 4000

NAME

Centrifuge

TYPE

Thermo Scientific

BRAND

75006580


SKU

🕒 00:10:00 Centrifuge

🌡️ 15 °C




Discard liquid and repeat centrifuge step for 3 times total.

5  00:10:00 2nd Centrifuge

 00:10:00 3rd Centrifuge

Add Binding Buffer.

6  20 mL Binding Buffer


Transfer it to 50 ml falcon tube.

7 In this step 35  $\mu$ l of total sample should be obtained. Binding Buffer involves 5 ml Imidazole, 20 ml TBS, ddw up to 200 ml.

For the resuspension of pallet, tubes with Binding Buffer can be put into shaker.



8

Sonication will be occurred.

9 For sonication, plastic test tubes are usually used for sonication. This sonication process is done at 3 cycles for 7 minutes.  00:07:00 Sonication for 1 cycle.

Repeat this step for 3 times total.

Centrifuge at 13000 rpm for 30 minutes.

10  00:30:00 Centrifuge  6 °C