

Aug 14, 2023

Version 1

Protein Extraction from Dental Enamel V.1

DOI

dx.doi.org/10.17504/protocols.io.8epv5j8n4l1b/v1

Alexandra Burnett^{1,2}, simon.daled¹

¹ProGenTomics, Laboratory of Pharmaceutical Biotechnology, Faculty of Pharmaceutical Science, Ghent University;

²ArcheOs Laboratory for Biological Anthropology, Archaeology Department, Ghent University



Alexandra Burnett

Ghent University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5j8n4l1b/v1>

Protocol Citation: Alexandra Burnett, simon.daled 2023. Protein Extraction from Dental Enamel. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5j8n4l1b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working.

Created: April 21, 2023

Last Modified: August 21, 2023

Protocol Integer ID: 80877

Keywords: Enamel, archaeology, dental, proteomics, protein, protein extraction from dental enamel, powdered dental enamel for proteomic analysis, protein extraction, extracting protein, proteomic analysis, powdered dental enamel, dental enamel, extraction, protein, synthetic peptide, sample cleanup step, optional addition of synthetic peptide, sample

Funders Acknowledgements:

Ghent University Special Research Fund (BOF)

Grant ID: BOF.GOA.2022.0002.03

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol details a method for extracting proteins from ~5mg of powdered dental enamel for proteomic analysis by LC-MS/MS. This procedure includes the use of alkylating and reducing agents, a sample cleanup step, and features the optional addition of synthetic peptide heavy standards for absolute quantification and quality control of the analysis.

Total protocol duration is dependent on the number of samples to be processed. Protein extraction will typically take a day-and-a-half, after which time the samples can be safely stored prior to instrumental analysis.

Guidelines

Always process blank samples alongside your enamel samples to allow for the observation of laboratory contamination.



Materials

Consumables

1.5mL protein lo-bind Eppendorf® tubes
Strata™ -X 33 µm Polymeric Reversed Phase 30mg/1mL tubes
Pipettes
Pipette tips
Tube racks
MS sampling vials
Ice

Chemicals

⊗ Ultrapure Water

⊗ Ammonium Bicarbonate BioUltra ≥99.5% (T) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #09830-1KG**

⊗ Hydrochloric acid 37% a.r. 37+% HCl **Chem-Lab Analytical bvba Catalog #CL00.0310**

⊗ Trichloroacetic acid a.r. 99.5+% C₂HCl₃O₂ **Chem-Lab Analytical bvba Catalog #CL00.2037**

⊗ Urea 99.5-100.5% CH₄N₂O **Chem-Lab Analytical bvba Catalog #CL00.2101**

⊗ 2-propanol **Biosolve Catalog #162641**

⊗ Methanol for analysis EMSURE® ACSISOREag. Ph Eur **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1060092500**

⊗ Acetone p. 99+% C₃H₆O **Chem-Lab Analytical bvba Catalog #CL00.0102**

⊗ Formic acid 99% ULC/MS - CC/SFC **Biosolve Catalog #069141**

⊗ DTT for Biochemistry 99+% C₄H₁₀O₂S **Chem-Lab Analytical bvba Catalog #CL00.0481**

⊗ S-Methyl methanethiosulfonate 97% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #208795-10G**


⊗ Trypsin/Lys-C Mix, Mass Spec Grade, 100ug **Promega Catalog #V5072**

Equipment

Centrifuge (Eppendorf Centrifuge 5417R)
Centri-Vap SpeedVac Concentrator (Thermo Scientific Savant SPD111V SpeedVac Concentrator)
Sonicator (Elma Transsonic 460)
Vortex (Scientific Industries Vortex-Genie 2)
Drill (Dremel; Proxxon)
Diamond-tipped drill bits (Dremel; Proxxon)
Thermo-Shaker
Solid-Phase Extraction vacuum pump and chamber
Analytical balance
Tube racks for 1.5mL Eppendorf® tubes

Troubleshooting

Safety warnings

 This protocol includes the use of machinery and hazardous substances. Be sure to wear appropriate protective equipment, including a lab coat and gloves, at all times.



Sample acquisition

1

Note

Tip: label and weigh an Eppendorf tube for each sample prior to drilling. Use a fresh drill bit for each tooth and clean between uses to avoid transferral of enamel.

Optional: process additional, empty Eppendorfs alongside your samples to observe laboratory contaminants introduced throughout sample preparation (suggested 1:9).

Place your tooth on a fresh piece of weighing paper on a clean worktop.

Using a small drill bit, gently abrade the surface of the tooth crown to remove powdered enamel and collect it on the weighing paper.

- 2 Carefully tip the powdered enamel into a protein lo-bind Eppendorf tube. Weigh the tube and record the weight of the powdered enamel inside.

Protein Extraction

3

Safety information

This section includes the use of hazardous substances. Be sure to wear appropriate protective equipment including a lab coat and gloves.

Note

The following extraction protocol is predicated on archaeological enamel quantities of ~5mg. Smaller or larger solvent volumes (where listed) may be required for differing sample weights.

Tip: process extraction blanks alongside your samples as contamination controls.

4 **Demineralisation**

Add 200µL of 1.2M HCl to each sample tube. Vortex the tubes briefly to mix.

Incubate for two hours in a thermo-mixer, using moderate agitation at room temperature.

5 **Precipitation**

Add 100% TCA (in water) to a final concentration of 33% (V/V) (i.e. 100µL).



Incubate on ice for one hour.
Centrifuge for 10 minutes at 15,000G at 0°C.
Pipette off the supernatants and discard (or retain and freeze if desired).

6 **Wash the pellet**

Pipette 500µL of ice-cold acetone into each sample tube.
Centrifuge for 10 minutes at 15,000G.
Pipette off the supernatants and discard (or retain if desired).

7 **Resuspend in buffer**

Add 200µl of 8M urea in 1M ABC buffer (pH 8).
Vortex and sonicate until all pellets are resuspended.

8 **Reduction**

Add DTT to a final concentration of 5mM.
Incubate for 30 minutes at 37°C.

9 **Alkylation**

Create a fresh stock solution of MMTS in 2-propanol (200mM).
Add MMTS to a final concentration of 10mM in each sample.
Incubate at room temperature in the dark for 15 minutes.

Note

MMTS degrades in light and must be kept in darkness before use. Try to minimise exposure as much as possible and only prepare the alkylating solution immediately prior to use.

10 **Buffer dilution**

Dilute the samples to a final concentration of 2M urea with 1M ABC buffer.

11 **Digestion**

Prepare a solution of trypsin, or trypsin and lysine-C, in 1M ammonium bicarbonate.
Add enzyme to a 1:25 (w/w) enzyme/protein ratio (around 1-2µg per sample).
Incubate overnight at 37°C. If there is visible debris in the tubes following digestion, sonicate until this is broken up.

12h

12 **Sample Cleanup via Solid-Phase Extraction [Optional]**

Condition one Strata tube per sample with 1mL MeOH.
Equilibrate the Strata tubes twice with 1mL ultrapure water.

**Note**

Be careful not to let the filters inside the tubes dry out while conditioning, equilibrating, and loading the samples. It is best to stop draining the tubes ~1mm above the filter to avoid this. If a filter runs dry after conditioning, start over again.

Remove the waste MeOH and water and replace the waste collection tubes with fresh Eppendorf lo-bind tubes.

Slowly load 1mL of sample per tube.

Wash twice with 1 mL of 5% MeOH in ultrapure water. [The cartridge can be allowed to run dry at this stage.]

Elute the sample with 1 mL of 1% FA in MeOH; first let a small amount soak the cartridge for 2 minutes and then slowly elute the sample with the rest of the elution buffer.

13 Drying down

Vacuum-dry the samples. Drying at up to 40°C can accelerate the process without damaging the peptides.

Note

The samples can now be stored (refrigerated) in a dry state for transport or until it is time for LC-MS/MS. If storing for longer time periods (more than 2 days) it is recommended to freeze them at this point.

Resuspension for LC-MS/MS

- 14 [Optional] Prepare a solution containing amelogenin heavy labelled peptide standards in 0.1% formic acid. Add heavy peptides to a concentration of around 5fmol/ μ L, depending on the LOD and sensitivity of the LC-MS/MS system.

Add between 10-30 μ L of 0.1% formic acid [and heavy standards] to each sample tube.

Note

The volume of formic acid in which you resuspend your samples will be dependent on the desired injection volume, and the sensitivity of the analytical instrument to be used.

Sonicate and vortex the samples to ensure thorough solubilization.

Centrifuge at 16,000G for 10 minutes.

- 15 Transfer the samples to MS sampling vials ready for analysis.



Note

Tip: there may still be visible pellets in the bottom of the Eppendorf tubes. Avoid transferring any of this material into the MS vials. Should the pellet/s detach from the tube walls or begin to disintegrate into the supernatant, you may need to re-centrifuge the sample.

Protocol references

This protocol was used for analyses in the following publications:

Jackson, V., *et al.* (in preparation) 'Following the thread: understanding the complexities of medieval Ypres' spinning workforce using bioarchaeological and historical evidence'.