

Sep 26, 2019

Version 1

# • Protein gel sample preparation V.1

DO

dx.doi.org/10.17504/protocols.io.7p7hmrn

Cleo B.1

<sup>1</sup>Wageningen University



Cleo B.

# Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.7p7hmrn

Protocol Citation: Cleo B. 2019. Protein gel sample preparation . protocols.io

https://dx.doi.org/10.17504/protocols.io.7p7hmrn

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 26, 2019



Last Modified: September 26, 2019

**Protocol Integer ID: 28127** 

Keywords: sample, preparation, protein, SDS-page, gel, protein gel sample preparation, sample preparation for tricine sd, page gel, gel, tricine sd, sample preparation

#### **Abstract**

Sample preparation for tricine SDS-page gels and polyacrylamide SDS-page gels.

### **Materials**

#### **MATERIALS**

Beta-mercaptoethanol

MilliQ Water

tricine sample buffer Bio-Rad Laboratories Catalog #161-0739

**☎** 4x Laemmli Sample Buffer **Bio-Rad Laboratories Catalog #**161-0747

protein samples

# **Troubleshooting**



### cell cultures

OD correct cell cultures. For an OD of 0.6:

### cell cultures

1.1 Spin down 4 2 mL of culture at OD of 0.6 Remove supernatant Dissolve pellet in 4 100 µL MQ

## sample preparation tricine SDS page

2 If you are unsure about the amount of protein present in your sample, load multiple concentrations

Add beta-mercaptoethanol to tricine sample buffer until a final concentration of 2% Dilute 1 part sample with 1 part sample buffer

# sample preparation SDS page

3 If you are unsure about the amount of protein present in your sample, load multiple concentrations

Add  $\perp$  100  $\mu$ L of beta-mercaptoethanol per  $\perp$  900  $\mu$ L of 4x laemmli sample buffer Dilute 3 parts sample with 1 part sample buffer

4 Spin samples down ~ 10 sec

Place samples at \$\\$\\$\\$\\$\\$05 \cdots \text{for 5 min}\$

Load samples to gel

1m