

Aug 14, 2020

## Protein isolation from the extract of coelomic fluid of the *Lumbricus terrestris*.

DOI

[dx.doi.org/10.17504/protocols.io.bjsaknae](https://dx.doi.org/10.17504/protocols.io.bjsaknae)

Angel A Justiz-Vaillant<sup>1</sup>

<sup>1</sup>University of the West Indies St. Augustine

University of the West In...

[angel.vaillant@sta.uwi.e...](mailto:angel.vaillant@sta.uwi.edu)



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bjsaknae](https://dx.doi.org/10.17504/protocols.io.bjsaknae)

**Protocol Citation:** Angel A Justiz-Vaillant 2020. Protein isolation from the extract of coelomic fluid of the *Lumbricus terrestris*. [protocols.io <https://dx.doi.org/10.17504/protocols.io.bjsaknae>](https://dx.doi.org/10.17504/protocols.io.bjsaknae)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** August 14, 2020

**Last Modified:** August 14, 2020

**Protocol Integer ID:** 40482

- 1 Six earthworms, *Lumbricus terrestris*, after removal from the earth are placed in a 15 ml glass tube and washed thoroughly with PBS, pH 7.4.
- 2 The body of the annelid is pierced with an 18-gauge needle.
- 3 The coelomic fluid content of the annelid is added to 2 ml of PBS, pH 7.4.
- 4 After incubation at 4°C for 1h, 0.9 ml chloroform is added.
- 5 The mixture is centrifuged at 1000 x g for 5 min at 4°C.
- 6 The supernatant (1.1 ml) is collected and equal volume of cold ethanol is added drop wise to the preparation.
- 7 The mixture is then incubated at 4°C for 12 h.
- 8 Then, it is centrifuged at 4°C for 5 min.
- 9 The pellet (protein extract) is resuspended in 0.3 ml of PBS, pH 7.4.
- 10 The preparation is dialyzed against 1L of PBS, pH 7.4 for 12 h at 4°C.
- 11 The protein concentration is assessed by the Bradford method or any available method.