



Dec 03, 2024

## Protein Lysate and Immunoblotting

 [The Journal of Cell Biology](#)

DOI

[dx.doi.org/10.17504/protocols.io.5qpvo9bmdv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvo9bmdv4o/v1)

Amanda Bentley-DeSousa<sup>1,2,3,4,5</sup>, Shawn Ferguson<sup>1,2,3,4,5,6</sup>

<sup>1</sup>Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>2</sup>Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>3</sup>Program in Cellular Neuroscience, Neurodegeneration and Repair;

<sup>4</sup>Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>5</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

<sup>6</sup>Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Amanda Bentley-DeSousa

Yale University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo9bmdv4o/v1>

**Protocol Citation:** Amanda Bentley-DeSousa, Shawn Ferguson 2024. Protein Lysate and Immunoblotting. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpvo9bmdv4o/v1>

**Manuscript citation:**

<https://doi.org/10.1101/2023.10.31.564602>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 02, 2024

**Last Modified:** December 03, 2024

**Protocol Integer ID:** 113416

**Keywords:** RAW 264.7, protein lysate, immunoblotting, lysate, protein, protocol details step

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's Disease**

Grant ID: ASAP-000580

## Abstract

This protocol details steps taken to create protein lysates and perform immunoblotting in RAW 264.7 cells.

## Materials

DMEM (Thermo Fisher Scientific, 11965-092)

FBS (Thermo Fisher Scientific, 16140-071)

Penicillin/Streptomycin (Thermo Fisher Scientific, 15140122)

CellStripper (Corning, 356230)

PBS [1.1 mM KH<sub>2</sub>PO<sub>4</sub> (J.T. Baker, 3246-01), 155.2 mM NaCl (Sigma-Aldrich, 3624-05), and 3 mM Na<sub>2</sub>HPO<sub>4</sub> (J.T. Baker, 3828-05)]

Lysis Buffer [50 mM Tris Base (American Bio, AB02000-05000), 150 mM NaCl (Sigma-Aldrich, 3624-05), 1% Triton X-100 (Sigma-Aldrich, X100), 1 mM EDTA (Sigma-Aldrich, 03690) supplemented with cOmplete mini EDTA-free protease inhibitor (Roche, 11836170001) and PhosSTOP (Roche, 4906837001)

Coomassie Plus Protein Assay Reagent (ThermoFisher Scientific, 23236)

Laemmli Buffer [80 mM Tris-HCl pH 6.8 (American Bio, AB020000-05000 / HCl, Sigma-Aldrich, 3624-06), 25.3% Glycerol (American Bio, AB00751), 2.67% SDS (American Bioanalytical, AB01920-00500), and Bromphenol Blue (Sigma-Aldrich, B5525) supplemented with 6.187% fresh β-mercaptoethanol (Sigma-Aldrich, M3148)]

4-15% miniPROTEAN TGX Stain-Free pre-cast gels (BioRad, 4568084/4568085/4568086)

Electrophoresis Buffer [24.76 mM Tris Base (American Bio AB02000-05000), 191.87 mM Glycine (American Bio AB00730-05000), 10 mL 10% SDS (American Bio, AB01920-00500)]

0.45 μm Pore Nitrocellulose Membrane (Thermo Fisher Scientific, 1620115)

Transfer Buffer [24.76 mM Tris Base (American Bio AB02000-05000), 191.87 mM Glycine (American Bio AB00730-05000), 20% Methanol (Sigma-Aldrich, 179337-4I-pB)]

Non-fat dry milk omniblock (American Bio, AB 10109-01000)

TBST [10 mM Tris Base (American Bio AB02000-05000), 150 mM NaCl (Sigma-Aldrich, 3624-05), 0.1% Tween 20 (Sigma-Aldrich, P7949)]

BSA (Sigma-Aldrich, A9647)

SuperSignal West Pico PLUS Chemiluminescence Substrate (Thermo Fisher Scientific, 34580)

SuperSignal West Femto Maximum Sensitivity Substrate (Thermo Fisher Scientific, 34095)

## Troubleshooting








## Day 0

- 1 Plate  $1 \times 10^6$  RAW 264.7 cells per well in a 6-well dish.

## Day 1

1h 41m

- 2 Add the indicated compound at the provided concentration and time point in fresh DMEM media containing FBS and P/S.
- 2.1 Proceed directly to the next step if not treating cells.
- 3 Wash each dish 2X with cold PBS on ice.
- 4 Remove the PBS, scrape each well with  50  $\mu$ L of ice-cold lysis buffer, and transfer the lysates to Eppendorf tubes.
- 4.1 Lysis Buffer: 50 mM Tris Base, 150 mM NaCl, 1% Triton X-100, and 1 mM EDTA supplemented with cOmplete mini EDTA-free protease inhibitor and PhosSTOP.
- 5 Incubate on ice for  00:05:00 .
- 6 Centrifuge  14000 rpm, 4°C, 00:08:00 lysates. 8m
- 7 Measure the protein concentration of the supernatant using Coomassie Plus Protein Assay Reagent as per the manufacturer's instructions.
- 8 Mix the supernatant/lysate 1:1 with Laemmli Buffer and heat at  95 °C for  00:03:00 . 3m
- 8.1 Laemmli Buffer: 80 mM Tris-HCl pH 6.8, 25.4% Glycerol, 2.67% SDS, and Bromophenol Blue supplemented with 6.187% fresh B-mercaptoethanol.



9 Perform electrophoresis using between 20-30 uG of protein in 4-5% miniPROTEAN TGS Stain-Free pre-cast gels.

9.1 Electrophoresis Buffer: 24.76 mM Tris Base, 191.87 mM Glycine, and 1% SDS.

10 After electrophoresis, transfer gels onto 0.45 um pore nitrocellulose membrane at 100V for 01:00:00 .

1h

10.1 Transfer Buffer: 24.76 mM Tris Base, 191.87 mM Glycine, and 20% methanol.

11 Block membranes in 5% non-fat dry milk omniblock in TBST for 00:30:00 .

30m

11.1 TBST: 10 mM Tris Base, 150 mM NaCl, 0.1% Tween 20.

12 Add antibodies in 5% BSA TBST 4 °C Overnight at the indicated concentrations.

## Day 2

1h 20m

13 Wash membranes 2X 00:10:00 with TBST.

10m

14 Add secondary antibodies in either 5% milk TBST (phospho antibodies) or TBST for 01:00:00 at Room temperature .

1h

15 Wash membranes 3X 00:10:00 with TBST.

10m

16 Develop membranes via chemiluminescence (SuperSignal West Pico PLUS Chemiluminescence Substrate or SuperSignal West Femto Maximum Sensitivity Substrate) and image with a Biorad Chemidoc MP Imaging Station.