

May 07, 2019

Protein Purification (Column Cromotography)

DOI

dx.doi.org/10.17504/protocols.io.2ntgden



Sevval Uysalcan¹, Rengin Yamur Akbiyik²

¹İstanbul Bilgi University; ²Istanbul Bilgi University



Sevval Uysalcan

Istanbul Bilgi University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: https://dx.doi.org/10.17504/protocols.io.2ntgden

Protocol Citation: Sevval Uysalcan, Rengin Yamur Akbiyik 2019. Protein Purification (Column Cromotography). **protocols.io** https://dx.doi.org/10.17504/protocols.io.2ntgden

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: In development

We are still developing and optimizing this protocol

Created: May 07, 2019

Last Modified: May 07, 2019

Protocol Integer ID: 22963

Keywords: purification of protein, purification techniques appropriate for this protein, purification method, purification technique, more of the purification method, purification, proteins in pure form, protein mixture, resulting protein mixture, different protein, protein, concentration of this protein, desired protein, thousands of different protein, enzyme, many enzyme, amino acid composition, cell membrane, nucleus membrane, cell particle

Abstract

It is very difficult to purify a desired protein from a cell and tissue in a pure state. If the concentration of this protein is low, it is necessary to select the purification techniques appropriate for this protein to separate it from thousands of different proteins and to obtain it in pure form. The methods used in the purification of proteins are quite advanced.

Today, many enzymes and non-enzyme proteins have been isolated in pure and crystalline form. To study the amino acid composition and sequence, molecular weight, and other physical properties of a protein, it is first necessary to purify the protein. In order to isolate proteins from a group of tissues and cells, there are certain procedures that should be done first. In the first stage, the tissue or cell group where the protein we want to work at most is selected. The tissue and cells are homogenized by using buffer solutions to keep the protein we are interested in. The cell membrane and the nucleus membrane must be broken down. This homogenate is then subjected to centrifugation at certain rates to remove cell particles. The resulting protein mixture is then removed in the supernatant state to purify the proteins in pure form using one or more of the purification methods.


Guidelines

Binding buffer: 5 ml Immidiazole, 20 ml TBS, ddw up to 200 ml. Wash Buffer: 8 ml Immidiazole, 20 ml TBS, ddw up to 200 ml.

Troubleshooting



1 2 ml of resin (Ni-beads) should be put onto columns. Wash column with double distilled water before using resin. Δ 2000 μL Resin



- 2 1X Binding Buffer should be added until fill the whole column.
- When all Binding Buffer is leaked, samples should be put into columns.
- Wash columns and sample until all residues are gone. In the washing step take from sample and from bradford and put into 96 well plate.
 - $\stackrel{L}{=}$ 10 μ L sample

Δ 200 μL Bradford



5 As a last step, to keep the protein that is purified from samples, Elution Buffer is added and

transferred to new centrifuge tubes.

Δ 1000 μL Elution Buffer