

Sep 26, 2019

## Protein Purification strep-tag FPLC

DOI

[dx.doi.org/10.17504/protocols.io.7qchmsw](https://dx.doi.org/10.17504/protocols.io.7qchmsw)

Cleo B.<sup>1</sup>

<sup>1</sup>Wageningen University

iGEM Wageningen 2019



Cleo B.

---

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.7qchmsw](https://dx.doi.org/10.17504/protocols.io.7qchmsw)

**Protocol Citation:** Cleo B. 2019. Protein Purification strep-tag FPLC . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.7qchmsw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 26, 2019

**Last Modified:** September 26, 2019

**Protocol Integer ID:** 28132

## Abstract

Protein purification of proteins containing a strep-tag using FPLC.


Necessary:


Have performed Protein expression using E. coli strain BL21DE3

To continue see: Protein gel sample preparation V2.

## Materials

### MATERIALS

 MilliQ Water

 Roche Complete Protease Inhibitor EDTA-Free tablets **Merck MilliporeSigma (Sigma-Aldrich) Catalog #5056489001**



## Buffers

- 1 L Buffer W:
  - 100 mM Tris-HCl (pH 8.0)
  - 150 mM NaCl
  - 1 mM EDTA

Filter this buffer with a 22  $\mu$ m filter



- 100 mL Buffer E:
- 100 mM Tris-HCl (pH 8.0)
  - 150 mM NaCl
  - 1 mM EDTA
  - 2.5 mM desthiobiotin

1L 0.5 M NaOH (filtered with 22  $\mu$ m filter)


1L MQ (filtered with 22  $\mu$ m filter)


## Protein extraction

- 2 Place post-induction culture from Protein expression using E. coli strain BL21DE3 on ice. For FPLC have 2 L of cell culture.

- 3 Centrifuge culture at 5000xg for  00:10:00 at  4 °C

1m

- 4 The pellet is resuspended in  1 mL buffer W per 100 ml cell culture containing one crushed cComplete mini tablet for 2 L cell culture.


- 5 The cells are sonicated (VS70 T rod, 25% 1 sec on 2 sec off for  00:05:00 . On ice water). For small amounts use the MS72 rod.

From this point be very sure to keep the cells on ice as much as possible.

- 6 The cell extract is centrifuged for  00:45:00 at 30 000 g.

4m



7 The supernatant is collected and filtered first with 0.45  $\mu$ m filter and then with 0,22  $\mu$ m filter. If there are small sample volumes: add  3 mL buffer W.

8

## FPLC

9 Connect computer to system

10 Set a manual alarm specific for the column you're using  
All speeds are dependent on the machine/column that is used

11 Wash all pumps with MQ

12 Wash the system with MQ (~20 mL)

13 Set valve position to waste

14 Connect column wet (MQ) with a low flow rate (~1 mL/min)  
Wash column with MQ  
Increase flow rate slowly

15 Equilibrate column with buffer W  
Slowly increase flow rate  
Wait until the conductivity and UV have stabilised

16 Wash peristaltic pump with 15 mL MQ  
Introduce cell lysate in column using the peristaltic pump (connect wet)  
Story flow through and a bit of lysed cells

17 Wash the column with buffer W until the lines stabilize

18 Start elution with buffer E  
Collect all the fractions per 1/1.5 mL  
Wash until the protein peak has passed



Place all the fractions on ice immediately  
Fractions can be placed on analytical SDS-page gels

- 19 Wash the column with MQ (3 CV)  
Wash the column with NaOH (3 CV)  
Wash the column with ethanol (decrease speed before changing to ethanol) (3CV)  
Wash the system with ethanol