

Nov 22, 2023

Proteinase K digestion

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbx9qqlpk/v1>

Martin T. Henrich^{1,2}, Fanni F. Geibl^{1,2}

¹Northwestern University, Feinberg School of Medicine;

²Philipps University Marburg, Clinic for Psychiatry and Psychotherapy



Enrico Zampese

Surmeier

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbx9qqlpk/v1>

Protocol Citation: Martin T. Henrich, Fanni F. Geibl 2023. Proteinase K digestion. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbx9qqlpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: November 22, 2023

Last Modified: November 22, 2023

Protocol Integer ID: 91333

Keywords: ASAPCRN, Brain, Tissue digestion, Proteinase K, synuclein aggregate, immunofluorescence, subsequent visualization by immunofluorescence, insolubility of alpha, determining insolubility

Funders Acknowledgements:

ASAP

Grant ID: 020551

JPB Foundation, NIH

Grant ID: NS 121174

Stichting Parkinson Fonds

Von-Behring-Röntgen Stiftung

Abstract

This protocol describes the steps for performing proteinase K digestion and subsequent visualization by immunofluorescence, which can be helpful in determining insolubility of alpha-synuclein aggregates.

Attachments



PK digestion_2011202...

17KB

Materials

- Pre-cut brain sections (30-50 μm thick)
- 0.1 M PB
- Triton-X100 (detergent)
- Normal donkey serum (S30-100ML, Sigma-Aldrich)
- Proteinase K (e.g. Proteinase K, Invitrogen, #4333793)
- Primary antibodies (e.g. Anti-ChAT AB144P, Merck; anti-TH ab113, abcam; anti-p62, ab109012, abcam; anti-alpha synuclein(phosphoS129), ab51253, abcam)
- Secondary antibodies
- 12-well-plates with netwell inserts (e.g. Corning Costar Netwell)
- Hard-drying mounting medium (recommended ProLong Diamond, ThermoFisher Scientific)
- Aluminum foil
- Microscopy slides
- Glass coverslips
- Designated container for storage of microscopy slides
- Water bath with integrated shaker, temperature adjustable (e.g. GFL-1083)
- Orbital shaker (e.g. Heidolph Duomax 1030)
- Platform shaker (e.g. Heidolph Unimax 1010)
- Optional: DAPI, Sigma-Aldrich, D9542-5MG, 1:10,000 of 5 mg/ml



Safety warnings



Recommended PPE:

- Lab coat/disposable gown
- Safety goggles
- Examination gloves

Before start

- Prepare 0.1 M PB (PB)
- Prepare PBT (0.1 M PB with 0.3% Triton X-100)
- Prepare NDS solution (10% normal donkey serum diluted in PBT)
- Prepare proteinase K solution (12 µg/ml in 0.1 M PB)
- Set water bath to 65°C



Proteinase K digestion - Procedure

18h 8m

- 1 -Place individual sections in 12-well-plate with netwell inserts in PB.
- 2 -Wash sections for 3× 5 min on platform shaker in PB, change solution between washing steps. 15m
- 3 -Place sections in 12-well-plate with netwell inserts in PK solution. Use 4 ml solution per well.
- 4 -Perform PK digestion by immediately putting the filled 12-well-plate in the water bath (65 °C) followed by 10 min shaking. 65 °C 10m
- 5 -Take well-plate out of the water bath and wash sections 4× 5 min on platform shaker in 0.1 M PB at room temp., change solution between washing steps. 20m
Room temperature
- 6 -Block sections for 60 min at room temp. in NDS solution in 12-well-plates. We recommend 4 ml solution per well for good results. 1h
- 7 -For incubation with primary antibodies, transfer sections in a new 12-well-plate but without netwell inserts. This allows better shaking overnight. Solution for incubation with primary antibodies should contain NDS solution and the respective primary antibodies diluted according to manufacturer recommendation. Incubate sections on orbital shaker with gentle shaking at 4°C **overnight** (time may vary). We recommend at least 1 ml solution for each well. 4 °C 16h
- 8 -On the next day, transfer sections into 12-well-plate with netwell inserts and wash sections 4× 5 min in PBT on platform shaker at room temperature. Change washing solution after each washing step. 20m
- 9 -Next, transfer sections back in 12-well-plate without netwell inserts for incubation with secondary antibodies.



- 10 -Incubate sections with fluorophore-conjugated, species-specific secondary antibodies in NDS solution for 2 hours on orbital shaker at room temperature. For dilution of the secondary antibodies we generally use 1:500. We recommend at least 1 ml solution for each well. Cover the plate with aluminum foil during antibody incubation to avoid light exposure. 🌡 Room temperature
- 11 -Thereafter, transfer sections back to 12-well-plate with netwell inserts and wash sections for 5× 5min in PBT. For all washing steps, 12-well-plate should be covered with aluminum foil to avoid light exposure.

2h

STEP CASE

Stain with DAPI 8 steps

-In case you want to additionally stain with DAPI, perform only 2 washing steps with PBT and add DAPI (1:10.000 of 5 mg/ml) to the third PBT washing step and let sections wash for 10 min. Perform a final washing step without DAPI in PBT solution. For all washing steps, 12-well-plate should be covered with aluminum foil to avoid light exposure.

- 12 -After washing, mount sections on microscopy slides using a fine brush.
- 13 -Let sections dry a few minutes.
- 14 -Apply a small amount of hard-drying mounting medium sufficient to cover the sections. Carefully avoid the formation of air bubbles. Gently apply a coverslip over the sections and the mounting medium.
- 15 -Let cure overnight in the **dark** at room temperature. 🌡 Room temperature

3m

Proteinase K digestion - After the Procedure

- 16 -Dispose of waste and excess reagents/solution according to institutional guidelines.



- 17 -Clean tools/working station.
- 18 -Once the mounting medium is cured, slides are ready for observation.
- 19 -Fully mounted microscopy slides should be stored in a designated container at 4°C until time of observation. 🌡️ 4 °C