Sep 26, 2023

Protocol for abscisic acid (ABA) extraction from plant seeds

DOI

dx.doi.org/10.17504/protocols.io.81wgbxnnqlpk/v1



Jian You Wang^{1,2}, Lamis Berqdar¹, Salim Al-Babili^{1,2}

¹The BioActives Lab, Center for Desert Agriculture, King Abdullah University of Science and Technology, Saudi Arabia;

²The Plant Science Program, Biological and Environmental Science and Engineering Division, King Abdullah University of Science and Technology (KAUST), Saudi Arabia

Salim Al-Babili: *Correspondence: salim.babili@kaust.edu.sa;

Metabolomics Protocols & Workflows Tech. support email: bbmisraccb@gmail.com





DOI: dx.doi.org/10.17504/protocols.io.81wgbxnnqlpk/v1

Protocol Citation: Jian You Wang, Lamis Berqdar, Salim Al-Babili 2023. Protocol for abscisic acid (ABA) extraction from plant seeds. **protocols.io** <u>https://dx.doi.org/10.17504/protocols.io.81wgbxnnqlpk/v1</u>

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: September 22, 2023

Last Modified: September 26, 2023

Protocol Integer ID: 88357

Keywords: abscisic acid, plant seeds, LC-MS/MS, plant hormone

Funders Acknowledgements: King Abdullah University of Science and Technology Grant ID: CDA neodomestication

Disclaimer

The authors declare no competing interests.

Abstract

The plant hormone abscisic acid (ABA) regulates seed dormancy and germination. Here, we present a protocol for ABA extraction from plant seeds. We describe necessary steps required for preparation and extraction, followed by liquid chromatography-mass spectrometry (LC-MS/MS) analysis for ABA quantification.

Attachments



Guidelines

Author Contributions

S.A.-B. and J.Y.W. conceived the project. J.Y.W., and L.B. conducted experiments. L.B., J.Y.W., and S. A.-B. wrote, reviewed, and edited the protocol.

Competing interests

The authors declare no competing interests.

Materials

List of equipment

- Calibrated balance,
- Grinder (Mixer Mill MM 400, RETSCH (®),
- Sonicator (Branson 5510R-DTH Ultrasonic machine),
- ultra-centrifuge (Centrifuge 5424 R, EppendorfTM),
- speedVac (Concentrator Plus, EppendorfTM),
- 0.2 um filter (Non-Sterile Syringe Filter 13 mm; WhatmanTM UnifloTM),
- 0.1 mL Micro-insert clear glass tube (fisher scientific),
- 1.5 mL LCMS vial.

Before start

To prepare before to start:

Prepare $\Delta 2 \text{ mL}$ of working solution per sample: 100% Methanol spiked with $\Delta 1 \mu \text{L}$ D6-ABA stock solution.

Sam	ple preparation	2h 42m
1	Grind 7-8 seeds in a Safe-Lock 2.0 mL Eppendorf tube with 3-4 beads for 00:01:00 – frequency 25-26 Hz.	1m
2	Weight 4 100 mg fine powder of ground seeds.	
3	Add 🗸 1 mL of Standard solution in each tube.	0ª
4	Sonicate the samples for 00:15:00.	15m
5	Centrifuge for (14000 rpm, 00:05:00), then transfer the supernatant to a new 2 mL	5m
	Eppendorf tube.	•
6	Add another 📕 1 mL of Standard solution.	Ø
7	Sonicate 🕐 00:15:00 .	15m
8	Centrifuge for (14000 rpm, 00:05:00).	5m
9	Transfer the supernatant ($\boxed{2}$ 2 mL) to a 2 mL Eppendorf tube.	
10	Dry the supernatant under vacuum by speedVac for 02:00:00.	2h
11	After evaporating the solvent, either keep the extracted samples at 20 °C or proceed to the next step.	
12	Re-dissolve the final extract in $\boxed{2120 \ \mu L}$ of acetonitrile : water [25:75 (v/v)] followed by $\boxed{00:01:00}$ sonication.	1m

- 13 Filter the re-suspended solution through a $\rightarrow \leftarrow 0.22 \ \mu m$ filter into 0.1 mL micro-insert 29×5.7 mm clear glass tubes with inserted vials.
- 14 Tap the bottle to remove any bubbles.

Sample quantification

- 15 ABA quantification is performed by LC-MS/MS using a UHPLC-Triple-Stage Quadrupole Mass Spectrometer (Thermo Scientific Altis) machine.
- Chromatographic separation is achieved on the Hypersil GOLD C18 Selectivity HPLC Columns (150 × 4.6 mm; 3 μm; fisher scientific) with mobile phases consisting of water (A) and acetonitrile (B), both containing 0.1% formic acid, and the following linear gradient (flow rate, 0.5 mL/min): 0–10 min, 15%–100 % B, follow it by washing with 100 % B for 00:05:00 and equilibration with 15 % B for 00:02:00 .
- 17 Inject \underline{A} 10 μ L of sample, maintain the column temperature at \underline{B} 35 °C for each run.

18 Set the MS parameters of Thermo ScientificTM AltisTM as follows:

- negative mode,
- ion source of H-ESI,
- ion spray voltage of 3000 V,
- sheath gas of 40 arbitrary units,
- aux gas of 15 arbitrary units,
- sweep gas of 0 arbitrary units,
- ion transfer tube gas temperature of § 350 °C ,
- vaporizer temperature of § 350 °C ,
- collision energy of 20 eV,
- CID gas of 2 mTorr, and
- full width at half maximum (FWHM) 0.4 Da of Q1/Q3 mass.

Note

The characteristic Multiple Reaction Monitoring (MRM) transitions (precursor ion \rightarrow product ion) were the characteristic MRM transitions (precursor ion \rightarrow product ion) were 263.2 \rightarrow 153.1 for ABA; 269.2 \rightarrow 159.1 for D6-ABA.

7m

7m

Protocol references

References

1. Wang, J. Y., Alseekh, S., Xiao, T., et al. (2021). Multi-omics approaches explain the growth-promoting effect of the apocarotenoid growth regulator zaxinone in rice. Communications Biology, 4(1), 1-11. https://doi.org/10.1038/s42003-021-02740-8

2. Chen, G.-T.E., Wang, J.Y., Votta, C., et al. (2023) Disruption of the rice 4-DEOXYOROBANCHOL HYDROXYLASE unravels specific functions of canonical strigolactones. bioRxiv Article 535333, 10.1101/2023.04.03.535333.

3. Jamil, M., Alagoz, Y., Wang, J.Y., et al. (2023). Abscisic acid inhibits germination of Striga seeds and is released by them as a rhizospheric signal providing competitive advantage and supporting host infestation. bioRxiv Article 548005, 10.1101/2023.07.06.548005