

May 16, 2020

Version 1

Protocol for Dephosphorylation of 5' ends of DNA using Quick CIP (NEB #M0525) V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.7r2hm8e>



New England Biolabs¹

¹New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com



New England Biolabs

New England Biolabs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7r2hm8e>

External link: <https://neb.com/protocols/2019/06/04/protocol-for-dephosphorylation-of-5-ends-of-dna-using-quick-cip-neb-m0525>

Protocol Citation: New England Biolabs 2020. Protocol for Dephosphorylation of 5' ends of DNA using Quick CIP (NEB #M0525).
protocols.io <https://dx.doi.org/10.17504/protocols.io.7r2hm8e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 27, 2019

Last Modified: March 29, 2021

Protocol Integer ID: 28186

Keywords: phosphatase, intestinal alkaline phosphatase, need for multiple phosphatase, multiple phosphatase, protocol for dephosphorylation, phosphates from dna, snp analysis recombinant for purity, snp analysis recombinant, native enzyme faster reaction setup, active in any restriction enzyme buffer, restriction enzyme buffer, using quick cip, dephosphorylation, less enzyme, enzyme, shorter incubation time flexible reaction condition, dna sequencing, phosphate, quick cip

Abstract

Quick CIP is a heat-labile version of calf intestinal alkaline phosphatase (CIP) purified from a recombinant source.

- Rapid and irreversible heat inactivation eliminates unwanted activity
- Improved storage stability versus native enzyme
- Faster reaction setup (no supplemental additives like zinc required) and shorter incubation time
- Flexible reaction conditions (active in any restriction enzyme buffer, no clean-up required)
- Less enzyme required (high specific activity), resulting in a lower cost per reaction
- No need for multiple phosphatases (Quick CIP removes 5'- and 3'- phosphates from DNA, RNA and dNTPs)
- Active on unincorporated dNTPs in PCR products - improves DNA sequencing and SNP analysis
- Recombinant for purity, consistency and value

Materials

MATERIALS

 Quick CIP **New England Biolabs Catalog #M0525**

 CutSmart® Buffer **New England Biolabs Catalog #B7204S**

Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.



Before start

Dephosphorylation of DNA 5'-ends using Quick CIP in a Restriction Enzyme Reaction

- The phosphatase can be added directly into the digestion reaction during or after DNA digestion
- Add  1 μL of Quick CIP for every 1 pmol of DNA ends (about  1 μg of a 3 kb plasmid) and incubate at  37 °C for 10 minutes
- Quick CIP is active in all NEB restriction enzyme buffers
- The restriction enzyme should be heat inactivated at the same time as the phosphatase after digest and dephosphorylation
- If restriction enzyme(s) cannot be heat inactivated, DNA purification is required before ligation



1 Prepare a  20 μL reaction as follows:

DNA	1 pmol of DNA ends *
CutSmart® Buffer (10X)	2 μl
Quick CIP	1 μl
H ₂ O, purified	to 20 μl **

* Note: 1 pmol of DNA ends is about 1 μg of a 3 kb plasmid.

** Scale larger reaction volumes proportionally.

2 Incubate at  37 °C for  00:10:00 .



3 Stop reaction by heat-inactivation at  80 °C for  00:02:00 .