



Jul 14, 2023

Protocol for neurophysiological imaging in human iPSC-derived neurons

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbxp8ylpk/v1>

Natalia Garcia Perez^{1,2}, Yi-Ning Kang^{1,2}, Emma Moles-Garcia^{1,2}, Tom Venneman^{1,2}, Pieter Vanden Berghe^{1,2}

¹KU Leuven; ²LENS



Natalia Garcia Perez

KU Leuven, LENS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbxp8ylpk/v1>

Protocol Citation: Natalia Garcia Perez, Yi-Ning Kang, Emma Moles-Garcia, Tom Venneman, Pieter Vanden Berghe 2023. Protocol for neurophysiological imaging in human iPSC-derived neurons. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbxp8ylpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

The protocol works but we are still developing and optimizing this protocol.



Created: July 05, 2023

Last Modified: July 14, 2023

Protocol Integer ID: 84533

Keywords: neurophysiological imaging experiments on different type, neurophysiological imaging experiment, neurophysiological imaging, human ipsc, derived neuron, neurons the protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's

Abstract

The protocol describes how to perform neurophysiological imaging experiments on different types of human iPSC-derived neurons.

Materials

Solutions:

HEPES buffer consists of: (in mM) 140 NaCl, 5 KCl, 10 HEPES, 2 CaCl₂, 2 MgCl₂, and 10 d-glucose, adjusted to pH = 7.4 using NaOH ([Li et al. 2022](#)).

Electrodes:

- Focal electrodes were constructed by soldering a 50 µm tungsten wire to a gold connector inserted into a glass pipette that was pulled using a Sutter instruments pipette puller. The tip of the glass pipetted was sealed by melting the glass in the open flame of a lighter. This electrode was connected to a Grass stimulator and (30V) pulses were given at 1-20 Hz. A silver electrode was used as a reference electrode. ([Li et al. 2019](#))
- For field stimulation: two Pt wire electrodes (0.5 mm diameter) were flattened and cut to form two planar sides which were spaced 1 mm apart in a plexiglass rod. This rod was positioned using a Narishige 3D manipulator such that the 1 mm space covered the field of view that was under investigation. The electrode was gently put on the coverslip and current (5-10 mA) stimuli were given to depolarize the neuronal networks.

iPSC culture and maintenance

- 1 The iPSC cell lines (WTC11 and KOLF2.1J) are cultured in Essential 8 medium on Geltrex-coated plates. The cultures are passaged every 3 days with 0.5mM EDTA.

iPSC-derived neurons generation

- 2 Cells were differentiated from human iPSCs (WTC11, KOLF2.1 cell lines) according to published protocols, with some modifications to allow for easier live-cell imaging.
 - 2.1 To generate midbrain dopaminergic-like neurons, we followed the [Kriks et al. \(2011\)](#) protocol, with the following modifications:
 - Simplified basal medium:
Days 1-20: Neurobasal-A + 1% GlutaMAX + 2% B27 Supplement + 1% N2 Supplement + 1% PenStrep
Days 20+: Neurobasal-A + 1% GlutaMAX + 2% B27 Supplement + 1% PenStrep
 - Additional 1:2 replating step at Day 8 of differentiation to avoid cell detachment due to overgrowth.
 - For the replating step at Day 20, the cells are replated under low-density conditions (30×10^3 - 100×10^3 cells per cm^2) as opposed to the recommended high-density conditions (300×10^3 - 400×10^3 cells per cm^2)
 - If there are non-neuronal cells present in the culture after Day 20, a treatment with 2 μM Ara-C for 24-48h is carried out.
 - 2.2 To generate enteric nervous system-like neurons, we followed the [Barber et al. \(2019\)](#) protocol, with the following modifications:
 - Additional replating step at Day 20-22. We plate 60×10^3 - 100×10^3 cells per cm^2 .
 - If there are non-neuronal cells present in the culture after Day 20, a treatment with 2 μM Ara-C for 24-48h is carried out.

Live-cell imaging

- 3 Using calcium as a proxy for neuronal activity, we assess calcium dynamics during basal conditions, as well as during stimulations.

3.1 For midbrain dopaminergic-like neurons:

1h 15m

- The cells are loaded with the calcium dye Fluo-4-AM at a concentration of [M] 50-100 nanomolar (nM) for ⌚ 00:10:00 at 🌡️ 37 °C .
- The cells are washed 3x with HEPES and left to incubate at 🌡️ 37 °C for ⌚ 00:05:00 .
- The coverslips are mounted into an imaging holder and the cells covered with HEPES. As an easy alternative, glass bottom culture dishes can be used.
- The cells can be imaged for up to ⌚ 01:00:00 .

3.2 For ENS-like neurons:

1h 10m

- The cells are loaded with the calcium dye Fluo-4-AM at a concentration of [M] 1 micromolar (μM) for ⌚ 00:10:00 at 🌡️ Room temperature , on a shaker (<50 rpm).
- The cells are washed 1x with HEPES and mounted into an imaging holder. The cells are covered with HEPES. As an easy alternative, glass bottom dishes can be used.
- The cells can be imaged for up to ⌚ 01:00:00 .

4 Using mitochondrial dyes such as MitoTracker (MitoTracker Green/Red/Red CMXRos), we investigate neuritic mitochondrial transport in these cells (Van Steenberg et al. 2022).

1h 15m

- The cells are loaded in [M] 5-7.5 nanomolar (nM) (for dopaminergic-like neurons) or [M] 20 nanomolar (nM) (for ENS-like neurons) MitoTracker for ⌚ 00:10:00 at 🌡️ 37 °C .
- The cells are washed 3x with HEPES and left to incubate at 🌡️ 37 °C for ⌚ 00:05:00 .
- The coverslips are mounted onto an imaging holder and covered with HEPES. As an easy alternative, glass bottom culture dishes can be used.
- The cells can be imaged for up to ⌚ 01:00:00 .

Stimulation

5 Depending on the experiment, we use a combination of different stimuli to assess neuronal physiology.

5.1 Electrical stimulation

- Electrodes (custom field electrodes and focal electrodes; see Materials for details)



- Stimulator: Master8 (A.M.P.I) driving the A385 Stimulus Isolator (World Precision Instruments)
- Stimulation protocols:
 - Field electrode
 - 3 - 100 pulses at 20 Hz - 100 Hz
 - Current: 5 - 10 mA
 - Focal electrode
 - Train duration: 2×1000 ms
 - Rate: 2×10 pps
 - Delay: 1×1ms
 - Duration: 3×0.1 ms
 - Volts: 3×10 volts

5.2 Local perfusion (for ENS-like cells, as in **Boesmans et al. 2019**)

- High-K⁺ ([M] 75 millimolar (mM))
- Substance P ([M] 1 micromolar (μM))
- DMPP ([M] 10 micromolar (μM))

Microscopy

- 6 A variety of microscopes can be used. We use wide-field systems (Zeiss Axiovert, Nikon coupled to Till Vision image acquisition hard/software) and confocal systems (Zeiss LSM 780-NLO, Zeiss LSM 880 AiryScan). The image acquisition speed should be at least 1-2 Hz in order to capture the relevant physiological events.

Laser powers and exposure times (generally ~ 2 mWatt for widefield, 20-80ms exposure per 500 ms) should be adjusted to avoid bleaching and to preserve neuronal health.

Protocol references

Barber et al. 2019

Boesmans et al. 2019

Kriks et al. 2011

Li et al. 2019

Li et al. 2022

Van Steenbergen et al. 2022