

May 13, 2024

Protocol for performing PINK1 siRNA knockdown in mouse embryonic fibroblasts (MEFs)

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1>

Enrico Bagnoli¹, Miratul Muqit¹

¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK



Francesca Tonelli

MRC-PPU at The University of Dundee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1>

Protocol Citation: Enrico Bagnoli, Miratul Muqit 2024. Protocol for performing PINK1 siRNA knockdown in mouse embryonic fibroblasts (MEFs). protocols.io <https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 12, 2023

Last Modified: May 13, 2024

Protocol Integer ID: 92308

Keywords: ASAPCRN, pink1 sirna knockdown in mouse embryonic fibroblast, sirna knockdown in mouse embryonic fibroblast, mouse embryonic fibroblast, mef, sirna knockdown, pink1, performing pink1

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000463

Abstract

This protocol details the siRNA knockdown in mouse embryonic fibroblasts (MEFs) for PINK1, but applicable for any other target.

Attachments



[926-2384.docx](#)

17KB

Guidelines

72-Hr siRNA knockdown in MEFs in a 6-well plate

- If using Immortalised MEFs – seed 100,000 cells per well.
- If using Primary MEFs – seed 200,000 cells per well.



Materials

Key Reagents

- siRNA as purchased from Horizon Discovery as ON-TARGETplus siRNA Reagents (Dharmacon™ siRNA solutions).

Note

Protocol for 5 nmol quantities of purchased siRNA

-  Lipofectamine™ RNAiMAX Transfection Reagent **Thermo Fisher Catalog #13778150**
- Oligomycin
- Antimycin
- MLi2



Day 1 - Cell Seeding and siRNA preparation

- 1 For Primary MEFs seed 200,000 cells per well in a 6-well plate at a total volume of  2 mL .
- 2 For Immortalised MEFs seed 100,000 cells per well in a 6-well plate at a total volume of  2 mL .
- 3

Day 1 - Cell Seeding and siRNA preparation

- 4 Resuspend the dried siRNA – protocol for 5 nmol of purchased siRNA to reconstitute at [ 10 micromolar (μM)].
- 4.1 Add  400 μL of RNA-free water. 
- 4.2 Add  100 μL of 5x siRNA Buffer. 
- 4.3 Incubate in hood for  00:05:00 , vortex vigorously and store at  -20 °C . 

Note

The protocol can be amended if the nmol quantity of siRNA purchased is higher/lower by adjusting the resuspension to ensure a final [ 10 micromolar (μM)] concentration.

Day 2 - siRNA Knock-Down

5m



- 5 Prepare mastermix according to the number of wells. Four tubes have to be prepared, 2 with Optimem and lipofectamine and one each with the PINK1 and scramble siRNA respectively.

Note

For PINK1 knockdown a final concentration of 25 nM of siRNA in each well is used.

5.1 Tube1 PINK1 siRNA:

Per well -  5 μL of PINK1 siRNA at [10 μM 10 micromolar (μM)] diluted in  100 μL OPTI-MEMTube 3

5.2 Tube 2:

Per well -  10 μL of Lipofectamine diluted in  100 μL of siRNA-OPTI-MEM

5.3 Tube 3 scramble siRNA

Per well -  5 μL of scramble siRNA at [10 μM 10 micromolar (μM)] diluted in  100 μL OPTI-MEM

5.4 Tube 4 (same as tube 2)

Per well -  10 μL of Lipofectamine diluted in  100 μL of siRNA-OPTI-MEM

- 6 Vortex slowly and combine the content of Tube 1 with tube 2 (PINK1 siRNA) and similarly with tube 3 and 4.

- 7 Vortex again slowly and incubate for  00:05:00 at  Room temperature .

5m



- 8 Add drop by drop,  200 μL of PINK1 siRNA mix to each well for PINK1 KD and similarly of the scramble siRNA for controls.





Note

For an experiment with 4x 6 well plates (WT and mutant cells treated with scramble or PINK1)
the following volumes can be used:

- Tube1 and tube 3: 60 μL siRNA + 1200 μL OPTiMEM
- Tube 2 and 4: 120 μL of lipofectamine + 1200 μL of otpimem
- This gives 2580 μL of each mix (180 μL spare)

Day 4 - Mitochondrial depolarization

5m

9 Make a 500x of Oligomycin/antimycin solution. The final concentration in the well is:

- oligo: 1 micromolar (μM)
- antimycin: 10 micromolar (μM)

9.1 Add 4 μL of this 500x solution will be to each well.



Note

For a 100 μL of 500x solution, we use:

	A	B
	Compound	Quantity
	Antimycin A 50 mM solution	10 ul
	Oligomycin 6.4 mM solution	7.8 ul
	DMSO	82.2 ul



- 10 Treat cell with Oligomycin/antimycin A for 24:00:00 . OA should be added 48:00:00 after the addition of the siRNA. Include DMSO control.

3d



Day 5 - Cell lysis

1h 50m

- 11 Prepare working stock of MLi2 at a concentration of 100 micromolar (μM) in DMSO.

- 12 Treat cells with MLi2(10 nanomolar (nM)) for 01:30:00 (2 μl /well).

1h 30m



- 13 Lyse cells using 50 μL of Lysis buffer/well



- 14 Lysate can be precleared by centrifugation 17000 x g for 00:15:00 at 4 $^{\circ}\text{C}$.

15m



- 15 Perform protein estimation and subject lysate to immunoblotting.