



Mar 18, 2024

Version 1

# Protocol for Preparing Brain Samples for MUSIC V.1

DOI

[dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1)

Wenxin Zhao<sup>1</sup>, Sheng Zhong<sup>1</sup>

<sup>1</sup>Shu Chien-Gene Lay Department of Bioengineering, University of California San Diego

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)



Wenxin Zhao

Shu Chien-Gene Lay Department of Bioengineering, University ...

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1>

**Protocol Citation:** Wenxin Zhao, Sheng Zhong 2024. Protocol for Preparing Brain Samples for MUSIC. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 15, 2024

**Last Modified:** March 19, 2024

**Protocol Integer ID:** 96794

**Keywords:** brain samples for music study, brain samples for music, preparing brain sample, brain sample, music study, sample, detailed procedure, music, procedure, study

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

## Abstract

Here states the detailed procedure to prepare brain samples for MUSIC study.

## Troubleshooting



## Tissue pulverization and crosslinking

- 1 Cut a portion of post-mortem human brain frontal cortex sample on dry ice with heavy razor blades, and collect 50 mg of the sample in a 1.5 mL LoBind tube.
- 2 Thaw the 50 mg of brain sample on ice, and chop the tissue into smaller pieces by pestle. Store the rest of the sample at -80°C.
- 3 Incubate the sample with 10 mL of 2 mM disuccinimidyl glutarate (DSG) in 1X PBS in a 15 mL LoBind tube at room temperature for 45 min with gentle rotation.
- 4 Wash once with 10 mL of 1X PBS by centrifugation at 1,000 x g for 4 min.
- 5 Resuspend the sample in 15 mL of 1X PBS containing 3% formaldehyde, and incubate for 10 min with a gentle rotation.
- 6 Quench the crosslinking reaction by the addition of 5 mL of 1.25 M glycine followed by an incubation of 5 min with a rotation.
- 7 Centrifuge the sample at 1,000 x g for 4 min, and wash the sample twice with ice-cold 1X PBS containing 0.3% BSA (wt/vol).

## Nuclei isolation

- 8 Use Chromium Nuclei Isolation kit (10X genomics, 1000494) to isolate nuclei from crosslinked cortex samples.
- 9 Transfer 50 mg frozen tissue into pre-chilled sample dissociation tube.
- 10 Add 400 µL of Lysis Buffer to Sample Dissociation Tube. Dissociate tissue with plastic pestle until homogeneous.

	A	B
	Component	Volume (µL)



	A	B
	Lysis Reagent	1000
	Reducing Agent B	1
	Surfactant A	10
	<b>Total Volume</b>	1011

**Lysis Buffer**

- 11 Add 600  $\mu\text{L}$  of lysis buffer into the tube, and mix 10 times by pipetting. Incubate on ice for 10 min.
- 12 Equally load the solution into two nuclei isolation column, and centrifuge the tubes at 16,000 x g for 20 sec at 4°C.
- 13 Vortex the flowthrough in the collection tube that contains nuclei for 10 sec at 3,200 rpm or max speed to resuspend nuclei.
- 14 Centrifuge the collection tubes for 3 min at 500 x g at 4°C to pellet nuclei. Carefully discard the supernatant.
- 15 Resuspend the nuclei in 500  $\mu\text{L}$  of Debris Removal Buffer provided by the kit by pipetting 15 times.

	A	B
	<b>Component</b>	<b>Volume (<math>\mu\text{L}</math>)</b>
	Debris Removal Reagent	500
	Reducing Agent B	0.5
	<b>Total Volume</b>	500.5

**Debris Removal Buffer**

- 16 Centrifuge the nuclei at 700 x g for 10 min at 4°C. Carefully discard the supernatant.



- 17 Resuspend the nuclei in 1 mL of Wash and Resuspension Buffer.

A	B
Component	Volume (μL)
1X PBS	1750
10% BSA	200
RNase Inhibitor (40X)	50
<b>Total Volume</b>	2000

**Wash and Resuspension Buffer**

- 18 Centrifuge the nuclei at 500 x g for 5 min at 4°C. Carefully discard the supernatant.
- 19 Resuspend the nuclei again in 1 mL of Wash and Resuspension Buffer.
- 20 Centrifuge the nuclei at 500 x g for 5 min at 4°C. Carefully discard the supernatant as much as possible.
- 21 The nuclei are subjected to nuclei counting and the following procedures.