

Aug 03, 2023

Protocol for single cell and single nucleus prep for Cell Pellets/Multiple Myeloma (snRNA & snATAC)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l227oxl1y/v1>

Reyka Jayasinghe¹

¹Washington University in St. Louis



Reyka Jayasinghe

Washington University (St. Louis)

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l227oxl1y/v1>

Protocol Citation: Reyka Jayasinghe 2023. Protocol for single cell and single nucleus prep for Cell Pellets/Multiple Myeloma (snRNA & snATAC). protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l227oxl1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: August 03, 2023

Last Modified: August 03, 2023

Protocol Integer ID: 85943

Keywords: single nucleus prep for cell pellet, single nucleus prep, protocol for single cell, single cell, multiple myeloma, cell pellet, snatac, snrna

Abstract

Protocol for single cell and single nucleus prep for Cell Pellets/Multiple Myeloma (snRNA & snATAC)



Reagents and Protocols

1

All lysis buffers should be prepared fresh and stored at 4°C.

Buffers based on [Nuclei Isolation from Mouse Brain Tissue for Single Cell ATAC Sequencing Protocol from 10X](#)

Relevant 10X Protocols:

- [Nuclei Isolation for Single Cell ATAC Sequencing](#)
- [Nuclei Isolation from Mouse Brain Tissue for Single Cell ATAC Sequencing Protocol](#)
- [Chromium Single Cell ATAC Reagent Kits User Guide \(v1 Chemistry\)](#)
- [Chromium Single Cell 3' Reagent Kits User Guide \(v3.1 Chemistry\)](#)

Equipment

2

- Eppendorf DNA LoBind Tubes, 2.0 mL (022431048)
- Tips RT-LTS-A-1000µL-/L/S/W-768/8 (30389221)
- Flowmi® Cell Strainers, 40 µm, for 1000 µL Pipette Tips (BAH136800040)
- Sorvall™ ST 8 Small Benchtop Centrifuge (75007204 - refrigerated)
- Thermo Scientific™ Adapters for TX-150 Swinging Bucket Rotors (75-005-743)

Reagents to Order for Lysis and Wash Buffers:

3

- UltraPure™ 1 M Tris-HCl Buffer, pH 7.5 (Invitrogen, 15567027)
- NaCl (5 M), RNase-free (Invitrogen, AM9759)
- MgCl₂ (1 M) (Invitrogen, AM9530G)
- Nonidet P40 Substitute (Sigma, 74385-1L)
- MACS BSA Stock Solution (Miltenyi Biotec, 130-091-376)



- Phosphate Buffered Saline (1X) (Corning, 21-040-CM)
- SUPERase• In™ RNase Inhibitor (20 U/μL) (Invitrogen, AM2696)

Buffers

4 Stock Lysis buffer: Lysis Buffer

Reagent	Stock	Final	2 mL
Tris-HCl (pH 7.4)	1M	10 mM	20 mL
NaCl	5M	10 mM	4 mL
MgCl ₂	1M	3 mM	6 mL
Nonidet P40 Substitute	10%	0.1%	3 mL
Molecular Grade Water		-	1.974 mL

5 Lysis Dilution Buffer

Reagent	Stock	Final	10 mL
Tris-HCl (pH 7.4)	1M	10 mM	100 mL
NaCl	5M	10 mM	20 mL
MgCl ₂	1M	3 mM	30 mL
Molecular Grade Water		-	9.85 mL

6 0.1X Lysis Buffer

Reagent	Stock	Final	10 mL
1X Lysis Buffer	1X	0.1X	1 mL
Lysis Dilution Buffer	1X	0.9X	9 mL



7 Wash buffer: 1X PBS + 2% BSA + 0.2U/ul RNase inhibitor

Reagent	Stock	Final	10 mL
BSA	10%	2%	2 mL
PBS	1X	1X	8 mL
RNase Inhibitor	20 U/mL	0.2U/ul	1 mL

8 sn-ATAC-Seq submission

buffer: 1X Nuclei Dilution Buffer.

We store aliquots of 20X

Nuclei Dilution Buffer at -20°C in aliquots of no more than 25ul/each tube to minimize freeze/thaw. We make just enough fresh 1X Nuclei Dilution Buffer each time we use it. Add 425 uL Nuclease Free Water.

9 Trypan blue (2X) - filtered at 0.22 mm.

General Notes

- 10 • Keep everything on ice.
- Use wide-bore pipette tips for all steps when possible.
- Use RNase free reagents and consumables (Use filtered tips).
- For centrifugation steps, use swinging rotor bucket at 4°C.

Preparing sample for auto-MACS dead cell removal

- 11 1. Add 5 mL running buffer to tube
2. Thaw aliquot in 37°C water bath
3. Spray with 70% EtOH and wipe down
4. Add 1 mL thaw media to each aliquot and transfer to running buffer tube
5. Remove supernatant
6. Re-suspend each in 100 uL of beads and incubate at room temperature for 15 mins
7. Dilute 0.25 mL of 20x binding buffer with 4.75 mL of distilled water
8. Add 500 uL of 1x binding buffer to each tube

9. Run each through the DepleteS selection on the Automacs
10. Toss the positive fractions
11. Pull off 10 uL each from negative fractions to count
12. Proceed downstream with negative fractions, cells still suspended in buffer

Single Cell and Nucleus Prep for sc-RNA-Seq and sn-ATAC-Seq

- 12
 1. Spin down at 400 g for 5 min. Remove supernatant.
 2. Resuspend in 1000uL Wash Buffer.
 3. Spin down cells at 400 g for 5 min. Remove supernatant.
 4. Re-suspend in 200-1000 uL Wash buffer depending on size of the cell pellet with RNase Inhibitor. Count cells at this point – take 10uL of trypan blue and 10uL of sample and determine starting concentration of sample. If there is more than 1500 cells/uL then add more wash buffer and recount.
 5. Aliquot cell suspension to load approximately 10,000 cells/ul (refer to 10X 3' NEXT GEM protocol, Figure 2). Set aside another 50-100 uL of cell suspension for back up. Centrifuge remaining sample at 400g for 8' at 4°C.
 6. While remaining cell suspension is being spun down, load 10X CHIP G (single-cell). After loading chip you have 17 minutes to process remaining sample for sn-ATAC.
 7. After spin is finished, re-suspend cells in 500 µL of cold lysis buffer. Pipette gently for 25 times (until pellet is fully resuspended). Incubate on ice for 2 min.
 8. Add 400 uL Wash Buffer to same tube.
 9. Transfer the clear lysis-wash buffer suspension using the 40um flowmi cell strainer into a 2.0 mL tube. Centrifuge sample at 400g for 5' at 4°C.
 10. Once the sample is finished spinning, remove as much supernatant as possible. The pellet should be at the very bottom of tube if using ST8 Centrifuge. Add 7-10uL 1X Nucleus Dilution Buffer (sn-ATAC-Seq submission buffer). Gently resuspend until nuclei are completely re-suspended. Count nuclei at this point – Add 8 uL 1X Nucleus Dilution Buffer to 2 uL of sample and 10 uL of trypan blue. If concentration of nuclei is greater than 3,000 nuclei/uL then proceed with ATAC prep.
 11. Continue both 10X protocols accordingly.
 - a. Chromium Single Cell ATAC Reagent Kits User Guide (v1 Chemistry)
 - b. Chromium Single Cell 3' Reagent Kits User Guide (v3.1 Chemistry)