

Mar 18, 2025

Protocol for α -Synuclein Purification and Ionic Strength Modification Pivotal to High Yield and Reproducibility

DOI

dx.doi.org/10.17504/protocols.io.81wgbnd63gpk/v1

Chelva Janarthanam¹, Griffin Clabaugh¹, Manikandan Samidurai¹, Zerui Wang², Ileia Scheibe¹, Huajun Jin¹, Vellareddy Anantharam¹, Ramona J B Urbauer³, Jeffrey L Urbauer⁴, Arthi Kanthasamy¹, Xuemei Huang⁵, Vincenzo Donadio⁶, Wenquan Zou², Anumantha G Kanthasamy¹, Bradley Melvin⁷, Bradley R Melvin⁷

¹Isakson Center for Neurological Disease Research, Department of Physiology and Pharmacology, University of Georgia, Athens, GA 30602, USA.;

²Department of Pathology, Case Western Reserve University School of Medicine, Cleveland, OH 44106, USA.;

³Department of Chemistry, University of Georgia, Athens, GA 30602, USA.;

⁴Department of Chemistry, University of Georgia, Athens, GA 30602, USA;

⁵Department of Neurology, Pennsylvania State University College of Medicine, Hershey, PA 17033, USA;

⁶IRCCS Institute of Neurological Sciences of Bologna, Complex Operational Unit Clinica Neurologica, 40138 Bologna, Italy.;

⁷Carver College of Medicine



Piyush Padhi

University of Georgia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbnd63gpk/v1>

Protocol Citation: Chelva Janarthanam, Griffin Clabaugh, Manikandan Samidurai, Zerui Wang, Ileia Scheibe, Huajun Jin, Vellareddy Anantharam, Ramona J B Urbauer, Jeffrey L Urbauer, Arthi Kanthasamy, Xuemei Huang, Vincenzo Donadio, Wenquan Zou, Anumantha G Kanthasamy, Bradley Melvin, Bradley R Melvin 2025. Protocol for α -Synuclein Purification and Ionic Strength Modification Pivotal to High Yield and Reproducibility. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbnd63gpk/v1>

Manuscript citation:

Janarthanam, C.; Clabaugh, G.; Wang, Z.; Melvin, B.R.; Scheibe, I.; Jin, H.; Anantharam, V.; Urbauer, R.J.B.; Urbauer, J.L.; Ma, J.; et al. High-Yield α -Synuclein Purification and Ionic Strength Modification Pivotal to Seed Amplification Assay Performance and Reproducibility. *Int. J. Mol. Sci.* **2024**, *25*, 5988. <https://doi.org/10.3390/ijms25115988>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 18, 2025

Last Modified: March 18, 2025

Protocol Integer ID: 124533

Keywords: α -synuclein, Purification, SIZE-EXCLUSION, ION exchange, synuclein purification, synuclein seed amplification, promising diagnostic tools for parkinson, high yield α syn protein, protein purification protocol, fast protein liquid chromatography, isolated α syn, wild type α syn in bl21 escherichia, enhanced assay performance, ultra sensitivity of the assay afford, parkinson, high quantities of α syn, minute quantities of α syn, promising diagnostic tool, purification protocol, purification protocol for the efficient production, purification, wild type α syn, misfolded α syn, α syn, tissues for enhanced assay performance, separation between disease

Abstract

Here we describe a protocol for purifying a high yield α Syn protein purification protocol for the efficient production of monomers with a low propensity for self-aggregation. Alpha-synuclein seed amplification assays (α Syn-SAAs) have emerged as promising diagnostic tools for Parkinson's disease (PD) by detecting misfolded α Syn and amplifying the signal through cyclic shaking and resting in vitro. The ultra sensitivity of the assay affords the ability to detect minute quantities of α Syn in peripheral tissues, but it also presents various technical challenges in controlling batch-to-batch variability. To address the problem of variability, we expressed wild type α Syn in BL21 *Escherichia coli*, lysed the cells using osmotic shock, and isolated α Syn using acid precipitation and fast protein liquid chromatography (FPLC). Following purification, we optimized the ionic strength of the reaction buffer to better distinguish the fluorescence maximum (Fmax) separation between disease and healthy control tissues for enhanced assay performance. Our protein purification protocol yielded high quantities of α Syn (average: 68.7 mg/mL per 1 L of culture). Together, these methods are highly reproducible with high purity, stability and yield.

Attachments



[SYN purification PRO...](#)

2.2MB

Troubleshooting



1