



Feb 06, 2025

## Protocol for the production of RAB1A-10xHis in HEK293Gnti cells

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj9be5lk5/v1>

Annan Cook<sup>1</sup>

<sup>1</sup>University of California, Berkeley



Annan SI Cook

University of California, Berkeley

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.n2bvj9be5lk5/v1>

**Protocol Citation:** Annan Cook 2025. Protocol for the production of RAB1A-10xHis in HEK293Gnti cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.n2bvj9be5lk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** February 05, 2025

**Last Modified:** February 06, 2025

**Protocol Integer ID:** 119639

**Keywords:** ASAPCRN, hek293gnti cells protocol, protocol for the production, cell

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

## Abstract

Protocol for production of RAB1A-10xHis

## Materials

**Cell culture and transfection:**

HEK 293T GNTi cells

pCAG RAB1A-10xHis DNA

Polyethyleneimine (PEI) (Polysciences) 20,000 Da

Hybridoma media (Thermo)

Freestyle Media

**Buffers and chemicals:**

Phosphate-buffered saline (PBS)

Lysis buffer: 25 mM HEPES, pH 7.5, 200 mM NaCl, 2 mM MgCl<sub>2</sub>, 25 mM tris(2-carboxyethyl)phosphine (TCEP), 10% Glycerol

Protease inhibitor (EDTA-Free, Thermo Scientific)

Triton X-100

imidazole

Gel filtration buffer: 25 mM HEPES, pH 7.5, 150 mM NaCl, 1 mM MgCl<sub>2</sub>, 1 mM TCEP

**Equipment:**

Bucket Centrifuge

Ultracentrifuge

Pyrex Dounce homogenizer

Rocker

NiNTA resin (2 mL)

S75 10/300 column

Liquid nitrogen

## Cell growth and transfection

- 1 Grow 400 mL of HEK293T GNTi cells to a density of  $2.0\text{--}2.2 \times 10^6$  cells/mL.
- 2 Transfect the cells with 400  $\mu\text{g}$  of pCAG RAB1A-10xHis construct using a four fold molar excess of PEI in 80 mL Hybridoma

## Cell harvesting

- 3 After 48 hours, pellet the cells by centrifugation at 2200 rpm.
- 4 Wash and gently resuspend the pellet with PBS, centrifuge again at 500 x g for 10 minutes, and flash-freeze the pellet in liquid nitrogen. Store at  $-80^\circ\text{C}$ .

## Purification

- 5 Thaw the cell pellet at room temperature.
- 6 Resuspend the pellet in lysis buffer (25 mM HEPES (pH 7.5), 200 mM NaCl, 2 mM  $\text{MgCl}_2$ , 1 mM TCEP, 10% Glycerol) and add an EDTA-free protease inhibitor tablet (Thermo Fisher).
- 7 Transfer the resuspended pellet to a Pyrex Dounce homogenizer and homogenize the cells with 30 strokes.
- 8 Add Triton X-100 to the homogenate for a final concentration of 1%.
- 9 Gently mix and rock the mixture at  $4^\circ\text{C}$  for 1 hour.
- 10 Centrifuge the lysate at 17,000 rpm for 45 minutes at  $4^\circ\text{C}$ .



- 11 Apply the supernatant to 2 mL of NiNTA resin, repeating the application 3 times.
- 12 Wash the resin until the wash is free of protein, as determined by a Bradford assay.
- 13 Elute the protein using wash buffer supplemented with 300 mM imidazole.
- 14 Concentrate the protein and apply it to an S75 10/300 column equilibrated with buffer (25 mM HEPES, pH 7.5, 150 mM NaCl, 1 mM MgCl<sub>2</sub>, 1 mM TCEP).
- 15 Concentrate the protein to 50 μM.
- 16 Aliquot 50 μL of the purified protein and flash freeze in liquid nitrogen.