

Nov 01, 2024

Protocol of shotgun metagenomics library construction

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpzk9dlzp/v1>

Elizabeth Selene Gómez Acata¹, Braulio Ricardo Perez Alva², Yendi Navarro¹

¹Laboratorio de Interacciones Bióticas, Centro de Investigación en Ciencias Biológicas, Universidad Autónoma de Tlaxcala;

²Posgrado en Ciencias Biológicas, Universidad Autónoma de Tlaxcala



Elizabeth Selene Gómez Acata

Centro de Investigación en Ciencias Biológicas, Universidad ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpzk9dlzp/v1>

Protocol Citation: Elizabeth Selene Gómez Acata, Braulio Ricardo Perez Alva, Yendi Navarro 2024. Protocol of shotgun metagenomics library construction. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpzk9dlzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 18, 2024

Last Modified: November 01, 2024

Protocol Integer ID: 107913

Keywords: protocol of shotgun metagenomics library construction, shotgun metagenomics library construction, metagenomic library construction, optimization for metagenomic library construction, suitable for dna input, dna input, peat bog, dna, arable soil sample

Funders Acknowledgements:

CONAHCYT

Grant ID: 2649154

Abstract

This protocol is an optimization for metagenomic library construction from peat bog and arable soil samples. It is suitable for DNA inputs ranging from 20 pg to 10 ng, with a fragment size of 250 bp for Illumina sequencing.

Materials

QIAseq® FX DNA Library Core Kit Cat. No. 1120146

QIAseq® UDI Y-Adapter Kit B (96) Cat. No. 180314

AMPure XP Beckman Coulter Cat. No. A63880

FastGene Gel/PCR Extraction Kit Cat. No. FG-91302

Magnetic bead separation DynaMag™-2 Magnet Life technologies™ Cat. No. 12321D

Laptop cooler TrueNorth®

Quant-iT™ PicoGreen™ Cat. No. P7589

Ethanol (100% molecular grade) Sigma Aldrich Cat. No. E7023

QX200™ ddPCR™ EvaGreen Supermix Bio-Rad Cat. No. 186-4036

Tween 20 Sigma Aldrich Cat No. P1379-25ML

PCR tubes Axygen Cat. No. 321-02-051

qPCR tubes

Filter tips and pipettes

Equipment

Vortex Genie

Microcentrifuge Legend Thermo Scientific

Rotor-Gen Q Qiagen

Mastercycler Nexus Gradient Eppendorf

NanoDrop™ 3300 Fluorospectrometer

Before start

1. The workspace should be kept clean and sterile, using UV radiation for sterilization.
2. Library preparation needs to be carried out in a laminar flow hood.
3. Ensure pipettes are clean and calibrated.
4. A negative control should be included throughout the process to confirm the absence of reagent contamination.
5. To prevent batch effects, the same person should conduct the library construction.
6. Do not dissolve DNA in 1x TE or any solution containing EDTA.
7. Quantify DNA concentration using Quant-iT™ PicoGreen™ assay kit before starting the fragmentation process.
8. Each DNA sample must be normalized to 1 ng μL^{-1} using Tris 10 mM or water.

Thermocycler programs

- 1 Set up the thermocycler with the following programs:

Table 1. Fragmentation reaction program

	Step	Temperature	Time
	Pre-chilled	4°C	1 min
	Fragmentation	32°C	24 min*
	Enzyme inactivation	65°C	30 min
	Hold	4°C	Hold

*20 pg of input DNA, we used 14 min

Table 2. Adapter ligation reaction program

	Step	Temperature	Time
	Ligation	20°C	30 min
	Hold	4°C	Hold

Table 3. Library amplification program

	Step	Temperature	Time	Cycles
	Initial denaturation	98°C	2 min	1
	Denaturation	98°C	20 s	
	Annealing	60°C	30 s	10
	Extension	72°C	30 s	
	Final extension	72°C	1 min	1
	Hold	4°C		Hold

* 20 pg of input DNA, we used 16 cycles

Fragmentation



- 2
 1. Thaw the reagents on ice: Fx buffer 10x, FX Enzyme Mix and Fx Enhancer from QIAseq® FX DNA Library Core Kit.
 2. Keep the reactions on ice during setup.
 3. Prepare the following reaction mix for 10 ng or 20 pg of input DNA.

Table 4. Fragmentation reaction mix for samples containing 10 ng of input DNA

Component	Volume
Fx buffer 10x	5.0 µL
DNA (1 ng/µL)*	10.0 µL
Nuclease free-water	25.0 µL
Total	40.0 µL

*DNA dissolved in Tris 10 mM or water, do not contain EDTA

Table 5. Fragmentation reaction mix for samples containing 20 pg of input DNA

Component	Volume
Fx buffer 10x	5.0 µL
DNA (20 pg)*	20.0 µL
Nuclease free-water	12.5 µL
Fx Enhancer	2.5 µL
Total	40.0 µL

*DNA dissolved in Tris 10 mM or water, do not contain EDTA

4. Add 10 µL of FX Enzyme Mix to each reaction tube and mix thoroughly by pipetting up and down.
5. Place de PCR tubes into a pre-chilled thermocycler set to 4 °C.
6. Use the thermocycler program from Table 1, ensuring the fragmentation time matches the input DNA.
7. Proceed directly with adapter ligation.

Adapter ligation

- 3
 1. Prepare the adapter's dilutions according to Table 6.

Table 6. Adapter ligation factor

	DNA amount	Adapter dilution
	10 ng	1:15
	20 pg	1:300

2. Add 5 μL of diluted adapters to each fragmentation tube.
3. Add 45 μL of adapter ligation mix.

Table 7. Adapter ligation mix

	Component	Volume
	DNA ligase buffer 5x	20 μL
	DNA ligase	10 μL
	Nuclease free-water	15 μL
	Total	45 μL

4. Mix by pipetting up and down.
 5. Incubate 20 $^{\circ}\text{C}$ for 30 min (thermocycler program Table 2).
- Note:** do not use thermocycler with heated lid
6. Purify the adapter ligation tube immediately.

Adapter ligation cleanup

1. Use FastGene Gel/PCR Extraction Kit following the manufacturer instructions.
2. Elute the DNA with 40 μL of elution buffer.

Library amplification

- 4
 1. Thaw reagents on ice: HiFi PCR master mix 2x and Primer mix from QIAseq[®] FX DNA Library Core Kit.
 2. Use the following library amplification mix.

Table 8. Library amplification mix

	Component	Volume
	HiFi PCR master mix 2x	25.0 μL
	Primer mix (10 μM)	1.5 μL



Component	Volume
Library purified	23.5 μL
Total	50.0 μL

- Mix all the components by pipetting up and down.
- Use the library amplification program from table 3.

Library cleanup

Take 50 μL of the amplicon library and cleanup with FastGene Gel/PCR Extraction Kit following the manufacturer instructions.

Library quantification with qPCR

1. Prepare a 1:4000 dilution of the samples that used adapter 1:15, and a 1:500 dilution for those with adapter 1:300, using Tween 20 (0.1%).
Note: To prepare Tween 20 (0.1%), mix 10 mL of H_2O with 10 μL of Tween-20.
2. Thaw the reagents on ice: QX200™ ddPCR™ EvaGreen Supermix.
3. Use the following library amplification mix.

Table 9. Library quantification mix

Component	Volume
QX200™ ddPCR™ EvaGreen Supermix 2x	4.50 μL
Primer forward (10 μM)	0.18 μL
Primer reverse (10 μM)	0.18 μL
Diluted library	1.00 μL
Nuclease free-water	4.14 μL
Total	10.00 μL

- Mix all components thoroughly by pipetting up and down
- Open the lid of the Rotor-Gen Q, place the samples, and follow the instructions in the equipment's operation manual to start the quantification process.

Library pooling

1. Normalize each sample to 4 $\text{ng } \mu\text{L}^{-1}$ use the data obtained from qPCR.
2. Add 4 μL of each sample (4 $\text{ng } \mu\text{L}^{-1}$) to a tube and proceed with cleanup with AMPure XP Beckman Coulter beads.



Library pool cleanup

- 7
 1. Add 80 μ L of AMPure XP Beckman Coulter beads to each 75 μ L of pooled library and mix by pipetting.
 2. Incubate for 5 minutes at room temperature.
 3. Place the tube on a DynaMagTM-2 Magnetic Separator (Life TechnologiesTM) for a few minutes until a brown pellet forms.
 4. Carefully remove the supernatant with a micropipette, ensuring the magnetic bead pellet is not disturbed.
 5. Wash the beads with 200 μ L of 80% ethanol and remove the supernatant.
 6. Repeat the 80% ethanol wash once more.
 7. Incubate the beads on the DynaMagTM-2 for 5-10 min, or until dry. NOTE: Do not overdry the beads, as this reduces DNA recovery efficiency.
 8. Remove the tube from the DynaMagTM-2 . Add 37.5 μ L of elution buffer (1x TE) and mix using the micropipette.
 9. Place the tube back on the DynaMagTM-2 to separate the magnetic beads from the solution.
 10. Transfer the eluate to a new sterile tube.