

Feb 14, 2025

Protocol: RNeasy Plant Mini Kit for RNA Extraction

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk9yqdv5r/v1>

Gonzalo Villarino¹

¹Cornell University

Puregene



Gonzalo Villarino

Cornell University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk9yqdv5r/v1>

Protocol Citation: Gonzalo Villarino 2025. Protocol: RNeasy Plant Mini Kit for RNA Extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk9yqdv5r/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 14, 2025

Last Modified: February 14, 2025

Protocol Integer ID: 120297

Keywords: rneasy plant mini kit for rna extraction, rna extraction process from cannabis flower, rna extraction process, qiagen rneasy plant mini kit, rna extraction method, rna extraction, warmed rnase, rneasy plant mini kit, standardized rna extraction method, rna purification, rna elution efficiency, based rna purification, high rna yield, rna, cannabis flower, cannabis, liquid nitrogen, extraction, freezing in liquid nitrogen, purification

Abstract

This protocol describes the RNA extraction process from cannabis flowers using the Qiagen RNeasy Plant Mini Kit. Given the lack of a standardized RNA extraction method for cannabis, this protocol provides optimized steps to ensure high RNA yield and integrity by minimizing degradation. The workflow includes sample collection in the greenhouse (GH), flash-freezing in liquid nitrogen (LN2), mechanical disruption, and column-based RNA purification. Pre-warmed RNase-free water at 37°C is used to maximize RNA elution efficiency.

Attachments



HB-0435-

007_HB_RNYMi...

735KB

Guidelines

- Work quickly to prevent RNA degradation.
- Keep all tools and surfaces RNase-free by cleaning them with RNaseZap or 70% ethanol.
- Use fresh gloves and sterile tools for each sample.
- Flash-freeze all plant material immediately after liquid nitrogen or dry ice collection.
- Keep sample tubes and reagents cold during processing unless otherwise stated.
- Ensure ethanol (96-100%) is freshly prepared before use.



Safety warnings



- Do not allow tissue to thaw at any point before lysis.
- Buffer RLT and RW1 contain guanidine salts—avoid using with bleach or acid-based disinfectants.
- Ethanol is flammable—keep away from flames and heat sources.
- Wear protective equipment (gloves, lab coat, goggles) when handling β -mercaptoethanol.
- Always use an RNA-dedicated pipette and work in an RNase-free area.

Before start

Equipment:

- Portable analytical balance
- Liquid nitrogen (LN2) dewar
- Mortar & pestle (pre-chilled)
- Microcentrifuge (capable of $\geq 10,000$ rpm)
- Pipettes (2 μL – 1000 μL) and RNase-free pipette tips
- Heat block set at 37°C (for pre-warming elution buffer)

Reagents & Consumables:

- Qiagen RNeasy Plant Mini Kit
- RNase-free water
- Buffer RLT (with β -mercaptoethanol, if needed)
- Buffer RW1
- Buffer RPE (ensure ethanol is added before use)
- Ethanol (96–100%)
- Sterile scissors & tweezers
- RNase-free 50 mL Falcon tubes (for sample freezing)
- RNase-free 2 mL & 1.5 mL microcentrifuge tubes
- QIAshredder spin columns (lilac)
- RNeasy spin columns (pink)



Preparation Before Going to the Greenhouse (GH)

- 1 Set up a portable **analytical balance** in the GH
- 2 Label **one tube** per cultivar. Use 50 mL Falcon tubes
- 3 Bring **RNase-free** gloves and sterile tools (scissors, tweezers).

Go to the greenhouse and collect fresh cannabis flowers

- 4 Cut fresh flowers using **sterile scissors** (one sample at a time)
- 5 Place the **weighed ~100 mg** sample into a labeled 50 mL Falcon tube. Poke some holes in lid of Falcon tube. 
- 6 **Flash-freeze** the plant material (in 50 mL Falcon tube) immediately in liquid nitrogen

Go to lab

1m

- 7 Grind the frozen material **thoroughly** with a **mortar and pestle** while keeping it in liquid nitrogen. 
- 8 Transfer the finely powdered material to an **RNase-free, pre-chilled 2 mL** microcentrifuge tube (not supplied)

RNA Extraction Protocol

- 9 Add **450 µL Buffer RLT** to the 2 mL with the plant material. **Vortex** (1 minute) **vigorously** for complete lysis. Make sure the sample **looks homogenous** (no visible chunks or fibrous material)
- 10 Transfer **ALL** the lysate (**~450 µL**) from the 2 mL tube to a QIAshredder spin column (lilac) placed in a 2 ml collection tube. Centrifuge for **2 min at full speed**.



- 11 Transfer the supernatant (**~400 µL**) of the flow-through to a new **RNase-free 1.5 mL** microcentrifuge tube (not supplied) without disturbing the cell debris pellet in the collection tube. Use only this supernatant in subsequent steps. Discard only the QIAshredder spin column.
- 12 Add **~200 µL** (0.5 volume of 100% ethanol) to the cleared lysate and mix immediately by **pipetting (DO NOT VORTEX)**. Do not centrifuge. Proceed immediately to step 13.
- 13 Transfer the sample (**~600 µl**), including any precipitate that may have formed, to an RNeasy spin column (pink) placed in a 2 ml collection tube (supplied).
- 14 Close the lid and centrifuge for **15 seconds** at $\geq 8000 \times g$ ($\geq 10,000$ rpm). Discard the flow-through
- 15 Add **700 µl** Buffer **RW1** to the RNeasy spin column. Close the lid and centrifuge for **15 seconds** at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash the spin column membrane. Discard the flow-through.
- 16 Add **500 µl** Buffer **RPE** to the RNeasy spin column. Close the lid and centrifuge for **15 seconds** at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash the spin column membrane. Discard the flow-through.
- 17 Add **500 µl** Buffer **RPE** to the RNeasy spin column. Close the lid and centrifuge for **2 min** at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash the spin column membrane.
- 18 Place the RNeasy spin column in a **new 2 ml collection tube** (supplied), and discard the old collection tube with the flow-through. Close the lid and centrifuge at full speed for **1 min**.
- 19 Pre-warming the **RNase-free water** to **37°C** in heat block
- 20 Place the RNeasy spin column in a **new 1.5 ml collection tube** (supplied). Add **30 µl RNase-free water** directly to the spin column membrane. **Let it sit for 1 minute** at room temperature. Close the lid. Centrifuge for **1 min** at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to elute the RNA.



Protocol references

Qiagen RNeasy Plant Mini Kit Handbook

- Qiagen (2023). RNeasy Mini Handbook. Qiagen GmbH, Hilden, Germany.
- [Available Online](#)