Jan 01, 2025

Version 2

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<https://dx.doi.org/10.17504/protocols.io.q26g7m4r3gwz/v2>

Ricardo Borges¹, Juliana de Souza Moreira¹, Rômulo Pereira de Jesus¹


¹Federal University of Rio de Janeiro, Brazil.

LAABio-IPPN-UFRJ

Ricardo M. Borges

UFRJ

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7m4r3gwz/v2>

Protocol Citation: Ricardo Borges, Juliana de Souza Moreira, Rômulo Pereira de Jesus 2025. Protocolo: Extração Ácido-Base de Alcaloides. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7m4r3gwz/v2> Version created by **Ricardo M. Borges**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

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English:

The acid-base extraction procedure for alkaloids is a widely used technique for isolating and purifying nitrogen-containing compounds from natural materials. Traditionally performed in separatory funnels, this protocol has been modified to be executed in centrifuge tubes, aiming for greater practicality, reduced solvent use, and suitability for small-volume samples. This adaptation not only optimizes workflow but also makes the method more appropriate for laboratories with limited infrastructure or those handling small-scale samples. The methodology retains the fundamental principles of acid-base extraction, enabling efficient isolation of alkaloids from various matrices.

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