

Apr 18, 2024

Version 2

Pseudoislets for diabetes research: Purifying and selective re-aggregating human islet cells V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl85x7dl2w/v2>



Wei Liu¹, Craig Dorrell², Xiaojuan Chen³

¹Columbia Center for Translational Immunology, Department of Medicine, Columbia University Medical Center, New York, NY, USA;

²Oregon Stem Cell Center, Oregon Health & Science University, Portland, OR, USA;

³Columbia Center for Translational Immunology, Department of Surgery, Columbia University Medical Center, New York, NY, USA

Wei Liu: Current address: Department of Nephrology, The First Hospital of Jilin University, Changchun, Jilin, 130021, China;

Human Islet Research N...



Sandy Beshir

City of Hope

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl85x7dl2w/v2>

Protocol Citation: Wei Liu, Craig Dorrell, Xiaojuan Chen 2024. Pseudoislets for diabetes research: Purifying and selective re-aggregating human islet cells. protocols.io <https://dx.doi.org/10.17504/protocols.io.5jyl85x7dl2w/v2> Version created by **Sandy Beshir**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 18, 2024

Last Modified: April 18, 2024

Protocol Integer ID: 98375

Keywords: Pseudoislets, diabetes, human islet cells, HIRN, aggregating human islet cell, modeling of human islet cell, pseudoislets for diabetes research, human islet cell, various combinations of human islet cell, diabetes research, pseudoislet, cell

Funders Acknowledgements:

NIH

Grant ID: UC4 DK104207

Abstract

In vitro modeling of human islet cells for diabetes research utilizing purified and then selectively re-aggregated various combinations of human islet cells.

Note

Corresponding Author

Xiaojuan Chen, MD, PhD

Address: 650 West 168th Street, BB1701, New York, NY 10032, USA

Email: xc2248@cumc.columbia.edu

Tel: (212)342-3111

Fax: (212)342-6030



Guidelines

REFERENCES:

Citation

1. Dorrell, C. et al
. Isolation of major pancreatic cell types and long-term culture-initiating cells using novel human surface markers. Stem Cell Research 1, 183-194 (2008).

<https://doi.org/10.1016/j.scr.2008.04.001>

LINK

Citation

2. Dorrell, C. et al. Isolation of mouse pancreatic alpha, beta, duct and acinar populations with cell surface markers. Mol Cell Endocrinol 339, 144-150 (2011).

<https://dx.doi.org/10.1016%2Fj.mce.2011.04.008>

LINK

Citation

3. Liu, W. et al. Abnormal regulation of glucagon secretion by human islet alpha cells in the absence of beta cells. EBioMedicine 50, 306-316 (2019).

<https://doi.org/10.1016/j.ebiom.2019.11.018>

LINK



Citation

4. Ricordi, C. et al

. National Institutes of Health-Sponsored Clinical Islet Transplantation Consortium Phase 3 Trial: Manufacture of a Complex Cellular Product at Eight Processing Facilities.

Diabetes 65, 3418-3428.

<https://doi.org/10.2337/db16-0234>

LINK

Materials

- ✕ CMRL-1066 supplemented CIT medium **Corning Catalog #98-304-CV**
- ✕ Heparin **Sagent Pharmaceuticals Catalog #25021-400-10**
- ✕ IGF-1 (Cell SciencesCat. #CRI500 **Cell Sciences Catalog #CRI500**
- ✕ fetal bovine serum **GE Healthcare Catalog #SH30088.03HI**
- ✕ Trypsin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SM-2004-C**
- ✕ DPBS **Gibco - Thermo Fisher Scientific Catalog #70011044**
- ✕ Glucagon **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G2654**
- ✕ Insulin **Dako Catalog #A0564**
- ✕ Somatostatin (SST ZYMED) **Thermo Fisher Scientific Catalog #18-0078 (discontiued)**
- ✕ ghrelin **Abcam Catalog #ab209790**
- ✕ Bovine serum albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A6003**
- ✕ PicoGreen DNA assay kits **Invitrogen - Thermo Fisher Catalog #P7589**
- ✕ Flavin adenine dinucleotide (FAD) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F7378**
- ✕ Propidium Iodide (PI) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P4170**
- ✕ Human Insulin ELISA **Mercodia Catalog #10-1113-01**
- ✕ Human Glucagon ELISA **Mercodia Catalog #10-1271-01**
- ✕ somatostatin EIA kits **Phoenix Pharmaceuticals Catalog #EK06003**

Equipment			
Ultra-Low Attachment plates			NAME
Corning			BRAND
3473			SKU
https://ecatalog.corning.com/life-sciences/b2c/US/en/Microplates/Assay-Microplates/96-Well-Microplates/Costar%C2%AE-Multiple-Well-Cell-Culture-Plates/p/3473			LI NK
			

Equipment	
confocal laser-scanning microscope	NAME
LSM410	TYPE
Zeiss	BRAND
None	SKU
https://www.olympus-lifescience.com/en/laser-scanning/fv3000/ ^{LINK}	

Equipment	
cell culture inserts	NAME
EMD Millipore	BRAND
PITP01250	SKU
https://www.emdmillipore.com/US/en/product/Millicell-Cell-Culture-Insert-12mm-polycarbonate-3.0m,MM_NF-PITP01250?ReferrerURL=https%3A%2F%2Fwww.google.com%2F	LI N K
	

1 Human islet culture and single cell preparation

Islets isolated from non-diabetic deceased human donors within 2-4 days post-isolation will be used in experiments. Islets and the dispersed cells are cultured at  37 °C with 5% CO₂ in a CMRL-1066 supplemented CIT medium (Cellgro, Cat. #98-304-CV) with 10 IU/ml Heparin (Sagent Pharmaceuticals, Cat. #25021-400-10), 1 µg/ml IGF-1 (Cell Sciences, Cat. #CRI500), and 10% fetal bovine serum (GE Life Sciences, Cat. #SH30088.03HI). For single cell preparations[1-3], islets are incubated with 0.025-0.05% trypsin (Millipore Sigma, Cat. #SM-2004-C) for approximately 6 min at  37 °C . Undispersed material is excluded using a 40 µm cell strainer. The dispersed cells are re-suspended in CMRL-1066 supplemented CIT medium with 2% FBS.

2 Fluorescence-activated cell sorting (FACS)

The dispersed cells are incubated with two monoclonal antibodies, HIC1-2B4-APC (<https://apps.ohsu.edu/research/tech-portal/technology/view/269524>) and HIC3-2D12-PE (<https://apps.ohsu.edu/research/tech-portal/technology/view/3848191>) at a 1:100 dilution. After incubation on ice for 30-40 minutes, the cells are washed with cold DPBS (Gibco, Cat. #70011-044), re-suspended in CMRL with 2% FBS, and sorted using the BD Influx sorter at the Flow Cytometry Core of Columbia Center for Translational Immunology (CCTI). A two-step sorting process, using the enrich mode followed by the pure sort mode, is performed to eliminate contamination from other pancreatic cell types.

3 Sorted human islet cell purity and culture

Sorted α -, β -, δ , and/or ϵ -cells are counted and cultured in  500 µL of CMRL-1066 with 10% FBS and 5.6 mM glucose for a period of 3-7 days in Corning Ultra-Low Attachment plates (Corning/Costar, Ithaca, NY). Each well contains at least 20,000 sorted islet cells, calculated based on the number of cells determined by flow cytometry. The cells are given 3-7 days to aggregate in culture. Intact human islets are cultured under the same conditions for control purposes.

Samples of the sorted cells are concentrated on glass slides by cytocentrifugation, fixed in 2.5% paraformaldehyde for 10 mins, and then immunofluorescence stained for human glucagon (Sigma, Cat. #G2654), insulin (Dako, Cat. #A0564), Somatostatin (SST, ZYMED, Cat. #18-0078) and ghrelin (Abcam, Cat. #ab209790) to confirm purity. Secondary antibodies conjugated with Dylight-dyes (Jackson ImmunoResearch) are used for fluorescence detection. Digital images of fluorescence-labeled cells are

acquired using a Zeiss fluorescence axial microscope attached with a digital camera or with a Zeiss LSM410 confocal laser-scanning microscope.

4 **In vitro glucose challenge**

Duplicate or triplicate samples from each donor are tested for each treatment condition. Cells from each donor are used as their own controls when testing and comparing function under the various treatment conditions.

After the 3-7 day culture, the cell aggregates and control intact islets are subjected to glucose challenge using a static incubation approach described for intact islets[4] with modifications. The cell aggregates or islets are carefully transferred into individual cell culture inserts with a pore size of 3 μm (EDM Millipore, Cat. #PITP01250) in 24-well plates (Non-tissue culture treated, Corning). The cells are first equilibrated in 2.0 mM glucose-containing HEPES-balanced Krebs-Ringer bicarbonate buffer (KRB) with 2 mg/ml Bovine serum albumin (Sigma-Aldrich, Cat. #A6003,) at  37 °C , 5% CO₂ for 1 hr. Following this, the inserts are removed from the wells, drained of any residual media and transferred to new wells containing KRB with 2.0 mM glucose and incubated at  37 °C for one hour. Afterwards, the inserts are transferred to new wells containing KRB supplemented with various concentrations of glucose and growth factors or compounds (e.g., Glibenclamide) for another hour. The cells used in the high-to-low glucose challenge are first incubated in KRB with high concentration glucose for one hour and then transferred to KRB containing 2.0 mM glucose for another hour. The supernatants from each of the one hour-incubations of the two step-glucose challenge assay are collected for islet hormone measurements. The cells are lysed for cell number quantification at the end of the experiment using the PicoGreen DNA assay kits (Invitrogen, Cat. #P7589) according to the manufacturer's protocol. In addition, subgroups of cells are also evaluated for viability using FAD (Sigma, Cat. #F7378) and PI (Sigma, Cat. #P4170) staining, and for glucagon content assessment by acid ethanol extraction (1.5% HCL in 75% ethanol) followed by freezing and neutralization steps prior to glucagon concentration measurement.

5 **Islet endocrine hormone measurements**

The supernatant collected are analyzed for insulin and glucagon concentrations using the human insulin and glucagon ELISA kits (Mercodia, Cat. #10-1113-01 and #10-1271-01) and somatostatin EIA kits (Phoenix Pharmaceuticals, Cat. #EK06003) according to the manufacturers' instructions. The hormone released by each group of cells was expressed as the ratio of the hormone released during the second step incubation divided by that released in the first step incubation of the two-step glucose challenge assay, as described above. The ratio represents the function of cells in each sample in a



cell-number independent manner. In addition, in order to compare the absolute amount of hormone secretion by a given number of cells under different experimental conditions, glucagon or insulin release was standardized to cellular DNA content measured as described above.