



Nov 12, 2023

Version 1

Pull-down Assays V.1

DOI

dx.doi.org/10.17504/protocols.io.x54v9pw5pg3e/v1

Leonardo A Parra-Rivas¹

¹University of California, San Diego



Leonardo A Parra-Rivas

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9pw5pg3e/v1>

Protocol Citation: Leonardo A Parra-Rivas 2023. Pull-down Assays. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9pw5pg3e/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 12, 2023



Last Modified: November 12, 2023

Protocol Integer ID: 90822

Keywords: assay

Abstract

Pull-down Assays

Troubleshooting



- 1 To pull down proteins from brain lysates and LP2(P4) fractions, 1-2 mg of the sample was incubated with 25-50 µg of glutathione beads containing GST fusion proteins for 12-16 hrs.

- 2 The mixtures were washed three times with 1X PBS with 0.15% Triton X-100, and then resuspended in 20 µl of 1X NuPAGE LDS sample buffer for NuPAGE and immunoblotted analysis.

- 3 Alternatively, the samples were resuspended in 50µl of PBS for proteomics studies and stored at -20 °C before quality control and proteomics analysis.