



Jan 13, 2025

Pulldown of protein aggregates with a biotinylated peptide

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5wrwg47/v1>

Cole Sitron¹, Patricia Yuste Checa¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5wrwg47/v1>

Protocol Citation: Cole Sitron, Patricia Yuste Checa, Victoria A Trinkaus, F Ulrich Hartl 2025. Pulldown of protein aggregates with a biotinylated peptide. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge5wrwg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 10, 2024



Last Modified: January 13, 2025

Protocol Integer ID: 102730

Keywords: ASAPCRN, pulldown of protein aggregate, biotinylated peptide, peptide, protein aggregate, biotinylated peptide the purpose, a β 42 monomer, affinity reagent, protein, same pulldown protocol, aggregates as prey, fractions from the pulldown, streptavidin bead, pulldown, immunoblot, endpoint of this protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Aligning Science Across Parkinson's

Grant ID: ASAP- 024268

Abstract

The purpose of this protocol is to test whether monomeric or aggregated versions of a protein can bind to a peptide. This is tested by performing a pulldown using an affinity reagent that consists of a biotinylated peptide conjugated to streptavidin beads. This protocol describes how to perform this pulldown with α -syn and A β 42 monomers and aggregates as prey, but the same pulldown protocol could be extended to any prey of interest. The endpoint of this protocol is the production of input, flow-through, and eluate fractions from the pulldown, which can be analyzed by SDS-PAGE and/or immunoblot.

Materials

- Pierce Magnetic Streptavidin Beads
- Biotinylated peptide of interest
- Peptide Solubilization Buffer: 45% ethanol, 10% acetic acid, 1 mM DTT in ddH₂O.
- 330 μ M α -syn A53T monomers, prepared as in [dx.doi.org/10.17504/protocols.io.btyyypve](https://doi.org/10.17504/protocols.io.btyyypve).
- 330 μ M α -syn A53T aggregates (PFFs), prepared as in [dx.doi.org/10.17504/protocols.io.btyyypve](https://doi.org/10.17504/protocols.io.btyyypve).
- Lyophilized A β 42 protein, prepared exactly as in described in the “without urea” protocol in Linse 2020 (Section 3.6)
- 1. Alternatively, this protein can be bought, e.g. Thermo Fisher Scientific cat. no. 03-111
- Blocking Buffer:
 - ⊗ PBS, pH 7.2 **Thermo Fisher Scientific Catalog #20012068** , 1%
 - ⊗ Triton[®] X-100 Surfact-Amps[®] Detergent Solution **Thermo Fisher Catalog #85111** , additional 113 mM NaCl (final concentration 250 mM), 3% ⊗ BSA **Cell Signaling Technology Catalog #9998S**
- IP Buffer: PBS pH 7.2, 1% Triton X-100, additional 113 mM NaCl (final concentration 250 mM), 0.1% SDS, 1 mM DTT
- 500 mM NaCl in ddH₂O.
- GE Buffer: 6M GuHCl, 20 mM sodium phosphate, pH 8.5
- A β 42 Aggregation Buffer: 20mM sodium phosphate pH 8.5, 0.2 mM EDTA
- ⊗ NUPAGE LDS sample buffer (4x) **Thermo Fisher Scientific Catalog #NP0007** + 5%
 - ⊗ 2-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**
- HU Buffer: 8 M urea, 5% SDS, 200 mM Tris-HCl pH 6.8, 1 mM EDTA, 0.01% bromophenol blue, 5% β -mercaptoethanol
- Formic Acid
- Superdex 75 Column (GE HealthCare cat. no. GE17-5174-01)
- FPLC Setup
- Microcentrifuge
- Magnetic tube rack
- Tube rotator
- ⊗ Bioruptor Plus sonication system **Diagenode Catalog #B01020001** (or equivalent)
- Nanodrop
- Thermomixer
- SpeedVac
- ⊗ Protein LoBind 1.5mL microcentrifuge tubes **Eppendorf Catalog #0030108116**
- 96-well Half Area Black/Clear Flat Bottom Polystyrene (Corning cat. no. 3881)



Solubilization and Aggregation of A β 42

1d

- 1 Equilibrate a Superdex 75 column with A β 42 Aggregation Buffer.
- 2 Resuspend  500 μ g of lyophilized A β 42 in  700 μ L GE Buffer.
- 3 Use an FPLC setup to run the solubilized A β 42 through the Superdex 75 column, collecting the fractions in low bind tubes.
- 4 Check the fractions by  205 nanomolar (nM) absorbance on a Nanodrop and pool the A β 42-containing fractions.
- 5 Dilute the A β 42 to  10 micromolar (μ M) in A β 42 Aggregation Buffer.
- 6 Aliquot and snap freeze half of the sample to use later as A β 42 monomers.
- 7 Place  80 μ L in each well of a 96-well plate and incubate for  24:00:00 at  37 $^{\circ}$ C .
- 8 Aliquot the aggregated material, flash freeze, and store at  -80 $^{\circ}$ C .

1d



Conjugation of the biotinylated peptide to Streptavidin resin (bait)

2h

- 9 Dissolve in Peptide Solubilization Buffer to an appropriate concentration.
- 9.1 In our case, we resuspend our peptide to  2 μ L which was  390.6 micromolar (μ M) .
- 10 Add  25 μ L of Pierce Magnetic Streptavidin Beads to a low bind tube and place on the magnetic rack.



- 10.1 Be sure to pipette with a cut pipette tip to avoid damaging the resin.
- 10.2 Prepare as twice as many tubes as prey proteins that will be pulled down. One tube will be for the normal pulldown and the other will be a no peptide control to assess background binding.
- 11 Wash the beads once with  1 mL IP Buffer. 
- 12 Resuspend the beads in  400 μ L IP Buffer.
- 13 Add  3.2 μ L of the solubilized biotinylated peptide to the half of the tubes with beads. 
- 14 To the other half add the same volume of Peptide Solubilization Buffer. 
- 15 And rotate the tube for  01:00:00 at  Room temperature to allow the biotinylated peptide to complex with the streptavidin. 
- 16 Wash the beads twice with  500 μ L of  500 millimolar (mM) NaCl, then once with  1 mL IP buffer. 
- 17 Incubate the beads with rotation in  500 μ L Blocking Buffer for  01:00:00 . 

- 18 Wash the beads once with IP Buffer. 

Preparation of protein aggregates and monomeric controls (prey) for pulldown

- 19 Thaw the A β 42 aggregates, α -syn monomers, and α -syn PFFs  On ice .



20 Centrifuge the A β 42 and α -syn monomers at  20000 x g, 4°C, 00:30:00 to remove any aggregated material.

30m



21 Remove the supernatant and place in a new tube.

22 Sonicate the α -syn PFFs in the Bioruptor for 5 cycles of  00:00:30 on and  00:00:30 off.

1m

23 Sonicate the A β 42 aggregates for 5 cycles of  00:00:05 on and  00:00:05 off.

10s

24 Dilute each sample to  2.5 micromolar (μ M) in IP Buffer.

Pulldown

1h 0m 30s

25 Add  400 μ L of  2.5 micromolar (μ M) A β 42 monomers, A β 42 aggregates, α -syn monomers, and α -syn to appropriate tubes.



25.1 Each sample should go into one tube with peptide-conjugated beads and into one peptide-free control tube.

25.2 Reserve a small amount of the  2.5 micromolar (μ M) sample to run as an input fraction.

26 Rotate the tubes for  01:00:00 at  Room temperature, except for the A β 42 monomers, which need to be incubated at  4 °C to prevent aggregation.

1h



27 Gently centrifuge the tubes (~  50 x g, 00:00:30) to remove the resuspended beads from the cap area of the tube.

30s



28 Place the tube on the magnetic rack and set aside the supernatant to run as a flow-through fraction.



29 Using the magnetic rack, wash the beads three times with  500 μL IP Buffer.



29.1 On the last wash, transfer to a new tube.

Normal Elution (appropriate for α -syn)

10m

30 Replace the IP Buffer with  50 μL 4X NuPAGE LDS Sample Buffer + 5% β -mercaptoethanol.

31 Gently flick the tube to resuspend the resin.

32 Boil the tube for  00:10:00 at  95 $^{\circ}\text{C}$.

10m

33 Use the magnetic rack to remove the supernatant, which will be the eluate fraction.

34 The input and flow-through fractions can now be denatured by boiling with the NuPAGE LDS Sample Buffer + β -mercaptoethanol (diluted 1:4).

Harsh Elution (appropriate for $\text{A}\beta_{42}$)

1h 50m

35 Replace the IP Buffer with  100 μL formic acid.

36 Add  90 μL formic acid to  10 μL of the input and flow-through fractions.



37 Incubate at  19 $^{\circ}\text{C}$ for  00:40:00 with shaking at  500 rpm in a Thermomixer.

40m



38 Use the magnetic rack to remove the supernatant from the beads and place the supernatant in a new tube.



39 Evaporate the all fractions containing formic acid in a SpeedVac for  01:00:00 at  45 °C .

1h

40 Add  50 μ L HU Buffer + 5% β -mercaptoethanol to the tubes and incubate for  00:10:00 at  60 °C shaking at  500 rpm in a thermomixer.

10m



Protocol references

1. Linse, S. (2020). Expression and Purification of Intrinsically Disordered A β 42 Peptide and Setup of Reproducible Aggregation Kinetics Experiment. In: Kragelund, B.B., Skriver, K. (eds) Intrinsically Disordered Proteins. Methods in Molecular Biology, vol 2141. Humana, New York, NY. https://doi.org/10.1007/978-1-0716-0524-0_38