

Oct 04, 2025

Purification

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11752vr2/v1>

igemnthuth¹

¹iGEM_NTHU_Taiwan

iGEM NTHU



igemnthuth

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11752vr2/v1>

Protocol Citation: igemnthuth 2025. Purification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov11752vr2/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 03, 2025

Last Modified: October 04, 2025

Protocol Integer ID: 228969

Keywords: qiagen dna purification kit, purification pcr product purification, purification, pcr, dna

Abstract

PCR product purification using the QIAGEN DNA Purification Kit.

Materials

- QIAquick PCR Purification Kit
- PCR reaction sample (20 μ L)
- Benchtop centrifuge (capable of $\geq 13,000$ rpm)
- Microcentrifuge tubes (1.5 mL)
- Pipettes and tips



Purification

1 **Mix PCR product with binding buffer**

Add 20 μ L of PCR reaction mixture into a clean 1.5 mL microcentrifuge tube. Then add 100 μ L of Buffer PB (5 $\times$ the PCR volume). Mix thoroughly with a pipette. Do not centrifuge.

2 **Load the mixture onto the purification column**

Carefully transfer the entire mixture onto the center of a QIAquick spin column placed inside a 2 mL collection tube. Ensure all liquid is applied directly to the center of the membrane.

3 **First centrifugation and discard flow-through**

Centrifuge at 13,000 rpm for 1 minute. Discard the flow-through collected in the tube, but keep the spin column.

4 **Wash the silica membrane**

Add 750 μ L of Buffer PE (pre-diluted with 96–100% ethanol according to the manual) to the spin column. Centrifuge at 13,000 rpm for 1 minute.

5 **Remove residual ethanol**

Discard the flow-through and place the column back into the empty collection tube. Centrifuge again at 13,000 rpm for 1 minute to completely remove residual ethanol, which could interfere with downstream applications.

6 **Elute DNA**

Transfer the spin column to a clean 1.5 mL microcentrifuge tube. Add 30 μ L of Buffer EB directly to the center of the membrane (avoid touching the column walls). Let stand at room temperature for 1 minute to allow DNA release from the silica membrane.

7 **Final centrifugation and collection**

Centrifuge at 13,000 rpm for 1 minute to collect the eluate. The flow-through contains the purified PCR product.