



Feb 01, 2024

Purification of GST-LC3B

 Nature Structural & Molecular Biology

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4qnbjvo5/v1>

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qnbjvo5/v1>

External link: <https://doi.org/10.1038/s41594-024-01338-y>

Protocol Citation: Elias Adriaenssens 2024. Purification of GST-LC3B. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl4qnbjvo5/v1>

Manuscript citation:

Elias Adriaenssens, Thanh Ngoc Nguyen, Justyna Sawa-Makarska, Grace Khuu, Martina Schuschnig, Stephen Shoebridge, Marvin Skulsuppaisarn, Emily Maria Watts, Kitt Dora Csalyi, Benjamin Scott Padman, Michael Lazarou, Sascha Martens (2024) Control of mitophagy initiation and progression by the TBK1 adaptors NAP1 and SINTBAD. Nature Structural & Molecular Biology doi: [10.1038/s41594-024-01338-y](https://doi.org/10.1038/s41594-024-01338-y)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: January 31, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94565

Keywords: ASAPCRN, purification of gst, lc3b this protocol, purification, gst, lc3b, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

Abstract

This protocol describes the purification of GST-LC3B.

Attachments



[ppgyb32rp.pdf](#)

29KB



Materials

Materials:

- pGEX-4T1 vector (LC3B-encoding vector is available from Addgene)
- Glutathione Sepharose 4B beads (GE Healthcare)
- Amicon filter (Merck Millipore)
- Superdex 75 16/600 column (Cytiva)

Lysis buffer:

	A	B
	HEPES pH 7.5	50 mM
	NaCl	300 mM
	MgCl ₂	2 mM
	βmercaptoethanol	2 mM
	cOmplete EDTA-free protease inhibitors (Roche)	
	DNase (Sigma)	

Wash buffer:

	A	B
	HEPES pH 7.4	50 mM
	NaCl	300 mM
	DTT	1 mM

High-salt wash buffer:

	A	B
	HEPES pH 7.4	50 mM
	NaCl	700 mM
	DTT	1 mM

SEC buffer:

	A	B
	HEPES pH 7.4	25 mM



	A	B
	NaCl	150 mM
	DTT	1 mM

Equipment

- Beckman Ti45 rotor



Purification of GST-LC3B

8h 30m 30s

- 1 To purify GST-LC3B, insert human LC3B cDNA in a pGEX-4T1 vector. Note, the last five amino acids of LC3B are deleted to mimic the cleavage by ATG4.
- 2 After the transformation of the pGEX-4T1 vector encoding GST-LC3B in E.coli Rosetta (DE3) pLysS cells, grow the cells in LB medium at 37 °C until an OD₆₀₀ of 0.8-1 and induce the protein expression with 1 millimolar (mM) IPTG for 04:00:00 at 37 °C .
- 3 Collect cells by centrifugation and resuspend in lysis buffer.

4h



Lysis buffer:

A	B
HEPES pH 7.5	50 mM
NaCl	300 mM
MgCl ₂	2 mM
βmercaptoethanol	2 mM
cOmplete EDTA-free protease inhibitors (Roche)	
DNase (Sigma)	

- 4 Sonicate cell lysates.
- 4.1 Sonicate cell lysates twice for 00:00:30 (1/2).
- 4.2 Sonicate cell lysates twice for 00:00:30 (2/2).
- 5 Clear lysates by centrifugation at 140000 x g, 4°C, 00:30:00 in a Beckman Ti45 rotor.

1m

30s

30s

30m





6 Collect and incubate the supernatant with pre-equilibrated Glutathione Sepharose 4B beads (GE Healthcare) for  02:00:00 at  4 °C with gentle shaking to bind GST-LC3B.

2h



7 Centrifuge the samples to pellet the beads and remove the unbound lysate.



8 Wash beads twice with wash buffer, once with high salt wash buffer and two more times with wash buffer.



Wash buffer

A	B
HEPES pH 7.4	50 mM
NaCl	300 mM
DTT	1 mM

High-salt wash buffer

A	B
HEPES pH 7.4	50 mM
NaCl	700 mM
DTT	1 mM

9 Elute the proteins  Overnight with  20 millimolar (mM) reduced L-glutathione in  50 millimolar (mM) HEPES  7.4 ,  300 millimolar (mM) NaCl,  1 millimolar (mM) DTT buffer.

8h



10 Collect the supernatant, filter through a 0.45 µm syringe filter, concentrate using a 10 kDa cut-off Amicon filter (Merck Millipore), and load onto a pre-equilibrated Superdex 75 16/600 column (Cytiva).

11 Elute the proteins with SEC buffer.

SEC buffer

A	B
HEPES pH 7.4	25 mM



	A	B
	NaCl	150 mM
	DTT	1 mM

12 Analyze fractions by SDS-PAGE and Coomassie staining.



13 Pool fractions containing purified GST-LC3B protein.

14 After concentrating the purified protein, aliquot the protein and snap-frozen it in liquid nitrogen. Store proteins at  -80 °C .

