

Apr 25, 2024

PURIFICATION OF PROTEINS FROM PFA FIXED SAMPLES BAK_WITH BIOTIN PULLDOWN FOR LC-MS/MS_2023

DOI

dx.doi.org/10.17504/protocols.io.36wgq3xyylk5/v1

Bryan Killinger¹

¹Rush University



Bryan Killinger

Rush University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgq3xyylk5/v1>

Protocol Citation: Bryan Killinger 2024. PURIFICATION OF PROTEINS FROM PFA FIXED SAMPLES BAK_WITH BIOTIN PULLDOWN FOR LC-MS/MS_2023. protocols.io <https://dx.doi.org/10.17504/protocols.io.36wgq3xyylk5/v1>

Manuscript citation:

Killinger BA, Mercado G, Choi S, Tittle T, Chu Y, Brundin P, Kordower JH. Distribution of phosphorylated alpha-synuclein in non-diseased brain implicates olfactory bulb mitral cells in synucleinopathy pathogenesis. NPJ Parkinsons Dis. 2023 Mar 25;9(1):43. doi: 10.1038/s41531-023-00491-3. PMID: 36966145; PMCID: PMC10039879.

Choi SG, Tittle T, Garcia-Prada D, Kordower JH, Melki R, Killinger BA. Alpha-synuclein aggregates are phosphatase resistant. bioRxiv [Preprint]. 2024 Apr 9:2023.11.20.567854. doi: 10.1101/2023.11.20.567854. PMID: 38645137; PMCID: PMC11030248.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 27, 2024

Last Modified: April 25, 2024

Protocol Integer ID: 95966

Keywords: pfa fixed samples bak-with biotin pulldown for lc, bak-with biotin pulldown, extraction of protein, biotin pulldown, purification of protein, biotin, mass spectrometry, protein, pfa fixed sample, purification, ms-2023 this protocol, immunoblotting, extraction, pfa

Funders Acknowledgements:

NINDS

Grant ID: 1R01NS128467

Michael J. Fox Foundation

Grant ID: ASAP-024442

Abstract

This protocol details the purification of proteins from PFA fixed samples and extraction of proteins from formalin-fixed tissues. Also included, biotin pulldown prior to LC-MS/MS. Samples generated for this protocol have been used for mass spectrometry, immunoblotting, and pulldowns.

Attachments



PURIFICATION OF

PROT...

221KB

Image Attribution

Killinger BA, Mercado G, Choi S, Tittle T, Chu Y, Brundin P, Kordower JH. Distribution of phosphorylated alpha-synuclein in non-diseased brain implicates olfactory bulb mitral cells in synucleinopathy pathogenesis. NPJ Parkinsons Dis. 2023 Mar 25;9(1):43. doi: 10.1038/s41531-023-00491-3. PMID: 36966145; PMCID: PMC10039879.



Materials

Wash Buffer (1L):

	A	B
	10X TBS pH 7.6	100 mL
	Triton X-100	10 mL
	SDS	1 g
	0.5M EDTA pH 8.0	2 mL

Dissolve in milliq water. Up to final volume 1L. Store at room temperature for 6 months.

Reversal Buffer (RB) (500ML):

	A	B
	SDS	25g
	0.5M EDTA pH 8.0	2 mL
	Tris-Base	30.35 g
	NaCl	4.38 g

Dissolve in milliq water. pH to 7.6. Up to final volume 500mL. Store at room temperature for up to 1 year.

PMSF solution (100 mM). Acute toxicity, handle powder carefully.

1. Place 1.5mL tube onto scale, tare.
2. Add small amount of PMSF to tube, record weight.
3. Add calculated volume of isopropanol to PMSF. Mix well. To calculate: $\text{PMSF mass (mg)} / 17.4 = \text{volume of isopropanol}$.

Stable in isopropanol at RT for at least 6 months. Stable for ~15 min once added to aqueous solution.

Fixation solution:























	A	B
	Ethanol	50%
	Acetic acid	10%

Troubleshooting



Procedure Day 1 (Extract proteins and test concentration)

1h 35m

- 1 Wash sections in DM 3 X 10 min. 
- 1.1 Wash sections in DM for  00:10:00 (1/3) 
- 1.2 Wash sections in DM for  00:10:00 (2/3) 
- 1.3 Wash sections in DM for  00:10:00 (3/3) 
- 2 Place sections in 1.5mL Eppendorf tube.
- 3 Add  0.5 mL of reversal buffer. 
- 4 Briefly sonicate on low power to disperse tissue.
- 5 Add  5 μ L of  100 millimolar (mM) PMSF. Mix well. Quick spin. 
- 6 Heat on block for  98 $^{\circ}$ C for  00:30:00 with cap locks. 
- 7 Remove from heat block carefully (caps will pop if not careful).
- 8 After  00:05:00 of cooling, vortex well. 
- 9 Centrifuge at  22000 x g for  00:30:00  Room temperature . 

- 10 Collect S1 (extracted proteins).







Western blot (Optional)


- 11 Perform methanol/chloroform cleanup on $100\ \mu\text{L}$ of S1. *
- 12 Resuspend resulting pellet in $100\ \mu\text{L}$ 5% SDS. *
- 13 Perform BCA assay on $100\ \mu\text{L}$. Each BCA test takes $20\ \mu\text{L}$. The remaining volume can be used for western blot. Alternatively, more protein can be prepared from S1 if needed. *
- 14 Calculate volume required for $20\ \mu\text{g}$ protein. *
- 15 $20\text{ug} / \text{concentration of S1 ug/ul} = \text{ul of sample required}$ *

Procedure Day 2 (Capture biotinylated proteins and wash)


10h 3m

- 16 Add S1 to $10\ \text{mL}$ TBST in 15mL conical tube. Mix well.  
- 17 Add $40\ \mu\text{L}$ of prepared magnetic streptavidin beads to each tube. 
- 18 Nutation for 01:00:00 $^{\circ}$ Room temperature . 1h
- 19 Place tube on magnetic stand for 00:01:00 (until all beads have been drawn out of solution). 1m
- 20 Carefully remove supernatant no to disturb beads. Add $10\ \text{mL}$ wash buffer. 
- 21 Nutation for 00:30:00 . 30m



22 Place tube on magnetic stand for  00:01:00 (until all beads have been drawn out of solution).


1m

23 Add  10 mL wash buffer.



24 Nutation for  00:30:00 .

30m

25 Place tube on magnetic stand for  00:01:00 (until all beads have been drawn out of solution).

1m

26 Add  10 mL wash buffer




27 Nutation  Overnight at  4 °C .

8h



Day 3 (Elute captured proteins)

12m

28 Place tube on magnetic stand for  00:01:00 (until all beads have been drawn out of solution).

1m

29 Carefully remove supernatant no to disturb beads.

30 Add  1 mL of wash buffer.




31 Mix samples until beads are suspended in wash buffer.



32 Using a pipette with tip cut off, transfer suspended beads to a 1.5mL protein low bind tube.



33 Place tube on magnetic stand for  00:01:00 (until all beads have been drawn out of solution).

1m

34 Remove buffer and surface wash with milliq 2X.





35 Quick spin and remove liquid from tube (should only have beads left).

36 Add  80 μL 1X SDS-page sample buffer containing reducing agent.



37 Mix well and quick spin.





38 Place on heat block set to  98 °C for  00:10:00 .

10m

39 Mix well and quick spin.



40 Place on magnet.

41 Transfer eluent to two separate lo-bind tube and discard used beads.  35 μL in 1 tube (used for premass spec QC), the remaining  45 μL in another (used for mass spec).

42 Immediately place in  -80 °C for long-term storage.




















Day 4 (Prepare proteins for Mass Spec)

3h 8m

43 Prepare Bis-Tris 4-12% wedge gel 15 well.

44 Make  800 mL 1X MOPS running buffer.




- 45 Load  5 μL of MW standard in one side of gel. 
- 46 For MS load exactly  40 μL of the samples into wells. Careful not to spill into opposing wells. Only load samples into every other well to prevent cross contamination. 
- 47 Ensure all wells have a buffer. For empty wells fill with  40 μL 1X loading buffer. 
- 48 Run the gel at 150V for a few minutes. Watch carefully until the sample has completely entered gel.
- 49 Stop running, remove gel with clean gloves, and clean equipment.
- 50 Fix gel in 100mL fixation solution (refer materials section) for  01:00:00 
 Room temperature .
- 51 Wash gel in several changes of milliq, until gel has swollen to original size. 
- 52 Place in  100 mL colloidal Coomassie blue stain solution. Cover loosely, and heat in microwave for few minutes until solution just starts to boil.
- 53 Incubate gel in heated colloidal Coomassie blue stain solution for  00:08:00 . 

- 54 Remove gel from stain solution and wash several times with milliq water. Clear background is achieved with ~  02:00:00 of washing. 

- 55 Using a clean razor excise the entire sample from the gel.
- 56 Place each sample into a clean 1.5mL tube. Add  500 μL of milliq water. 



57 Samples can be stored at  4 °C until submitted to mass spec core.

Optional: Pre mass spec QC

45m


58 **Estimating capture concentration - Purpose:** To provide mass spec core with a rough estimate of the amount of protein in the capture sample. 



Note

This prevents potential overloading, or digestion problems in the mass spec core.


58.1 Prepare gel as described above (Day 4, Steps 43-45)

59 Load  20 µL of QC aliquot into wells. 

60 Load  2 µg BSA into 1 well.


61 Run the gel at 150V for approximately  00:45:00 or until the dye has reached end of gel. 

62 Stain gel with Coomassie as described above (Day 4, Steps 50-54)

63 Image on Odyssey using NIR. 

64 Forward image to NW mass spectrometry core.


Example of results; order = 2ug BSA, 3 captures, 2ug BSA, 3 captures, 2ug BSA, 3 captures



Dot blot to estimate target enrichment

8h 0m 30s

65 Cut PVDF membrane into ~1X3 inch piece.


66 Activate PVDF in 100% methanol for  00:00:30

30s

67 Equilibrate activated PVDF in diH₂O.

68 Just prior to blotting samples, place hydrated PVDF on wypall and quickly pat dry.

69 Place PVDF on new dry wypall.

70 Spot  1 μ L of each sample onto PVDF.



71 Allow PVDF to dry completely. (~1h or  Overnight)

8h



72 Reactivate spotted PVDF as described above and equilibrate in dH2O.

73 Blot can now be probed with antibodies using standard protocols.