

Apr 30, 2020

Version 1

## Purification of RNA from a DNA/RNA Extract V.1

DOI

[dx.doi.org/10.17504/protocols.io.qdqds5w](https://dx.doi.org/10.17504/protocols.io.qdqds5w)

**SOWA**

Roey Angel<sup>1</sup>, Eva Petrova<sup>1</sup>, Ana Lara-Rodriguez<sup>1</sup>

<sup>1</sup>Soil and Water Research Infrastructure

Anaerobic and Molecular Microbiology Lab, Biology Centre CAS

Tech. support email: [eva.petrova@bc.cas.cz](mailto:eva.petrova@bc.cas.cz)



**Roey Angel**

Soil and Water Research Infrastructure

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.qdqds5w>

**Protocol Citation:** Roey Angel, Eva Petrova, Ana Lara-Rodriguez 2020. Purification of RNA from a DNA/RNA Extract.  
**protocols.io** <https://dx.doi.org/10.17504/protocols.io.qdqds5w>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 24, 2018

**Last Modified:** July 30, 2020

**Protocol Integer ID:** 12432

**Keywords:** RNA, DNase, nucleic acid purification, purification of rna, total nucleic acids extraction from soil protocol, rna extract, rna cleanup, total nucleic acids extraction, purification, rna, extraction, soil protocol, dna, concentration micro kit, following protocol,

## Abstract

The following protocol is intended as a downstream application for our [Total Nucleic Acids Extraction from Soil](#) protocol. This protocol describes how to purify RNA from a DNA and RNA extract using [TURBO™ DNase](#) and [GeneJET RNA Cleanup and Concentration Micro Kit](#). This protocol is a simplified and condensed version of the full protocols provided by the manufacturers.

## Attachments



[4393900B.pdf](#)

281KB



[MAN0012671\\_GeneJET\\_](#)

[R...](#)

276KB



[THE RNA storage solu...](#)

205KB



## Materials

### MATERIALS

- ✕ TURBO™ DNase (2 U/μL) **Thermo Fisher Scientific Catalog #AM2238**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**
- ✕ RNase AWAY™ Surface Decontaminant **Carl Roth Catalog #A998.4**
- ✕ THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**
- ✕ RNaseOUT Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**

### STEP MATERIALS

- ✕ USB Dithiothreitol (DTT) 0.1M Solution **Thermo Fisher Scientific Catalog #707265ML**
- ✕ RNaseOUT™ Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**
- ✕ TURBO™ DNase (2 U/μL) **Thermo Fisher Scientific Catalog #AM2238**
- ✕ Nuclease-free autoclaved DEPC-treated water **Carl Roth Catalog #T143.1**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**

## Protocol materials

- ✕ USB Dithiothreitol (DTT) 0.1M Solution **Thermo Fisher Scientific Catalog #707265ML**
- ✕ TURBO™ DNase (2 U/μL) **Thermo Fisher Scientific Catalog #AM2238**
- ✕ TURBO™ DNase (2 U/μL) **Thermo Fisher Scientific Catalog #AM2238**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**
- ✕ THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**
- ✕ RNaseOUT Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**
- ✕ THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**
- ✕ RNase AWAY™ Surface Decontaminant **Carl Roth Catalog #A998.4**
- ✕ RNaseOUT™ Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**
- ✕ Nuclease-free autoclaved DEPC-treated water **Carl Roth Catalog #T143.1**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ RNaseOUT™ Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**
- ✕ USB Dithiothreitol (DTT) 0.1M Solution **Thermo Fisher Scientific Catalog #707265ML**
- ✕ Nuclease-free autoclaved DEPC-treated water **Carl Roth Catalog #T143.1**
- ✕ TURBO™ DNase (2 U/μL) **Thermo Fisher Scientific Catalog #AM2238**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**
- ✕ THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**

## Troubleshooting

### Before start

1. For each sample prepare one Gene JET RNA Purification Micro Column tube and two RNase-free collection tubes (1.5 ml).
2. Add the required amount of ethanol to Wash Buffer 1 and Wash Buffer 2 (amount is dependent on kit size).



## DNA digestion

45m

- 1 Prepare the following mixture in a 1.5 ml tube:
  1. 10  $\mu\text{L}$  to 42  $\mu\text{L}$  of TNA extract ( 1  $\mu\text{g}$  to 3  $\mu\text{g}$  of DNA).
  2. 5  $\mu\text{L}$  TURBO DNase buffer 10x
  3. 1  $\mu\text{L}$  RNaseOUT
  4. 1  $\mu\text{L}$  0,1M DTT
  5. 1  $\mu\text{L}$  Turbo DNase per up to 2  $\mu\text{g}$  DNA
  6. Complete to 50  $\mu\text{L}$  with RNase-free water

TURBO™ DNase (2 U/ $\mu\text{L}$ ) **Thermo Fisher Scientific Catalog #AM2238**

RNaseOUT™ Recombinant Ribonuclease Inhibitor **Thermo Fisher Scientific Catalog #10777019**

USB Dithiothreitol (DTT) 0.1M Solution **Thermo Fisher Scientific Catalog #707265ML**

Nuclease-free autoclaved DEPC-treated water **Carl Roth Catalog #T143.1**

- 2 Incubate at 37 °C for 00:30:00 .

30m



### STEP CASE

## Extended DNA digestion 10 steps

If this procedure still leaves out undigested DNA (for example due to the presence of inhibitors), increase the incubation time (to 40–60 min) and add another equal dose of DNase half-way through.


## RNA purification

7m

- 3 Add 250  $\mu\text{L}$  Binding Buffer .

GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**



4 Add  300  $\mu$ L absolute ethanol .



Ethanol, Absolute, Molecular Biology Grade **Thermo Fisher Scientific Catalog #BP2818500**


5 Transfer the mixture to the GeneJET RNA Purification Micro Column preassembled with a collection tube. Centrifuge the column for

1m



14000 x g, Room temperature, 00:01:00 . Discard the flow-through. Place the GeneJET RNA Purification Micro Column back into the collection tube.



6 Add  700  $\mu$ L Wash Buffer 1 (supplemented with ethanol) to the GeneJET RNA Purification Micro Column and centrifuge for


1m



14000 x g, Room temperature, 00:01:00 . Discard the flow-through and place the purification column back into the collection tube.



GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**

7 Add  700  $\mu$ L Wash Buffer 2 (supplemented with ethanol) to the GeneJET RNA Purification Micro Column and centrifuge for

1m



14000 x g, Room temperature, 00:01:00 . Discard the flow-through and place the purification column back into the collection tube.




GeneJET RNA Cleanup and Concentration Micro Kit **Thermo Fisher Scientific Catalog #K0841**

8 Repeat step 7.

1m

 [go to step #7](#)



9 Centrifuge the empty GeneJET RNA Purification Micro Column for an additional  14000 x g, Room temperature, 00:02:00 to completely remove residual Wash Buffer.

2m






#### Note


This step is essential to avoid residual ethanol in the purified RNA solution. The presence of ethanol in the RNA sample may inhibit downstream enzymatic reactions.

10 Transfer the GeneJET RNA Purification Micro Column into a clean 1.5 ml Collection Tube tube.





- 11 Add  10  $\mu\text{L}$  to  20  $\mu\text{L}$  RNA storage solution or nuclease-free water to the GeneJET RNA Purification Micro Column. Centrifuge for

 14000 rpm, Room temperature, 00:01:00 to elute the RNA.


 THE RNA Storage Solution **Thermo Fisher Scientific Catalog #AM7000**

1m



- 12 Discard the purification column. Use the purified RNA immediately in downstream applications or store at   $-20\text{ }^{\circ}\text{C}$  or   $-80\text{ }^{\circ}\text{C}$  until use.

#### Note

For prolonged storage (more than 1 month), storage at   $-80\text{ }^{\circ}\text{C}$  is recommended.