


Oct 24, 2020

pYCR cloning strategy

DOI

dx.doi.org/10.17504/protocols.io.bnvxme7n



Andreea S¹

¹University of Groningen

iGEM Groningen 2020



a.stan.6

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bnvxme7n

Protocol Citation: Andreea S 2020. pYCR cloning strategy . protocols.io <https://dx.doi.org/10.17504/protocols.io.bnvxme7n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Other

Protocol has been developed based on literature, but hasn't been tested yet.

Created: October 23, 2020

Last Modified: October 24, 2020

Protocol Integer ID: 43671

Abstract

Introducing NLP14a in the genome of *B. mycoides*: In order to create plasmid pYCR-gamyNLP, the backbone of Pycr will be digested with the PCR product of gamy_Fw and gamy_Rv to produce the pYCR_gamy (CRISPR vector containing the gRNA). The NLP14a sequence will be ordered from Twist Biosciences And PCR-amplified with the primers NLP_Fw and NLP_Rv. The resulting PCR product, as well as the pYCR_gamy vector will be digested with SfiI. The digested mixtures will be ligated using T4 ligase.

Kill switch cloning: The suggested kill switch mechanism is based on Trp auxotrophic strains as well as a toxin-antitoxin mechanism. The toxin will be introduced in the genome using crispr in place of trpE. The antitoxin will be maintained as a cytoplasmic plasmid (pAD-YqcF) and expressed only in the presence of solanine. Unfortunately, to the best of our knowledge, no operator that binds solanine of *B. mycoides* has been described in literature. If it is proven that no such operator is present, we suggest only keeping the tryptophan dependence strategy.











In order to create plasmid pYCR-gtrpEYqcG plasmid, the backbone of Pycr will be digested with BsaI and ligated with the annealing product of gtrpe_Fw and gtrpe_Rv to produce the pYCR_gtrp (CRISPR vector containing the gRNA). Synthetic dna for *yqcG* will be ordered from Twist Biosciences and PCR-amplified with the primers Yqcg_Fw and Yqcg_Rv. The resulting PCR product, as well as the pYCR_gtrpE vector will be digested with SfiI. The digested mixtures will be ligated using T4 ligase

Protocol was adapted after a paper of Yi et al., 2018

1. Yi, Y., Li, Z., Song, C. & Kuipers, O. P. Exploring plant-microbe interactions of the rhizobacteria *Bacillus subtilis* and *Bacillus mycoides* by use of the CRISPR-Cas9 system. *Environ. Microbiol.* 20, 4245–4260 (2018).



Clone sgRNA sequence into pYCR

- 1 Design sgRNA spacer sequence (~20 nts) using "Benchling" and choose *B. mycoides* M2E_15 genome
Choose the gRNA with the highest on-target score and the lowest off-target score.
- 2 Order the 2 complementary oligos flanked by overhang containing BsaI restriction site
- 3 Anneal oligos by mixing the following components:
 -  70 μ L Nuclease Free water
 -  10 μ L Oligo DNA annealing buffer 10x ($[M]$ 100 millimolar (mM) Tris  8 ,
 $[M]$ 500 millimolar (mM) NaCl, $[M]$ 10 millimolar (mM) EDTA)
 -  10 μ L DNA oligos A and B ($[M]$ 10 micromolar (μ M) each)
- 4 Mix well and incubate for  00:05:00 at  90 °C 5m
- 5 Cool slowly to room temperature (aprox.  01:00:00) by removing the heat block from the apparatus 1h
- 6 Run agarose gel (1%) to check the annealing product
- 7 Digest pYCR with BsaI and inactivate the enzyme by PCR clean up
- 8 Ligate the annealing product with the pre-cut pYCR:
Mix:
 -  10 μ L μ L of annealing product
 -  60 ng of digested pYCR
 -  1 μ L T4 ligase
- 9 Transform the ligation mixture into *E. coli* MC1061 using the heat-shock method.



Grow the cells on LB agar plates + Erithromycin [M] 0.1 mg/mL

- 10 To select the transformants perform a colony PCR with Fw primer binding upstream of PvanP and Rv primer downstream Sfil cloning sites.

Expected result



Succesful transformation will have a size of 427 pb

- 11 Miniprep and sequence the constructs

Clone flanking regions into pYCR containing sgRNA

- 12 Select the repair fragment of choice:
-for knock-out: PCR amplify aprox. 1 kpb upstream and downstream of the sgRNA targeted region
-for knock-in: sequence of the desired inserction sequence flanked by Sfil restriction sites compatible with the backbone of pYCR

- 13 Digest the pYCR(+sgRNA) and the repair fragments with Sfil
Inactivate Sfil by PCR clean up

- 14 Ligate the components using T4 ligase
Heat inactivate the ligase (incubate  00:20:00  60 °C)

20m

- 15 Perform Smal digestion to get rid of the original pYCR

- 16 Transform the ligation mixture into E. coli MC1061 using the heat-shock methos.
Grow the cells on LB agar plates + erythromycin

- 17 To select the transformants perform a colony PCR with Fw primer binding upstream of PvanP and Rv primer downstream Sfil cloning sites.

- 18 Miniprep and sequence the constructs

Transform final construct in *B. mycoides* M2E_15

19 Transform using an electroporation protocol using  5 μ L plasmid DNA

Protocol





NAME

Cloning of *Bacillus mycoides*

CREATED BY

a.stan.6

PREVIEW

- 19.1 Pick 1 colony of *B. mycoides* M2E_15 and inoculate it in BHIS
- 19.2 When OD600nm reached 0.85 add 2% glycine and 2% threonine in order to weaken the cell wall.
- 19.3 Grow the cells overnight at  200 rpm, 30°C
- 19.4 Dilute the overnight culture 50 times in LBSP medium until the OD600nm reaches 0.65
- 19.5 Collect the cells by  4000 x g, 4°C, 00:10:00 . Discard supernatant
- 19.6 Wash cell pellet with pre-chilled electroporation buffer (10% glycerol, 0.25 M sorbitol, 0.25M trehalose) x4
- 19.7 Suspended in 1 ml electroporation buffer (10% glycerol, 0.25 M sorbitol, 0.25M trehalose).



19.8 Snap-freeze the electrocompetent cells in liquid nitrogen and store at -80 °C

19.9 For the PCR mix pipette:

- 5 µL primers mixture Fw + Rv (5 uM) (see table "Primers")
- 10 µL HF buffer
- 1 µL dNTP (10mM)
- DNA template (0.1 ng)
- 0.25 µL phusion polymerase
- MiliQ to 50 µL

19.10 PCR reaction is performed with the following protocol:













Temperature (*C)	Time (mm:ss)
98	05:00
98	10:00
Variable	00:30
72	00:30
repeat above steps	30x
72	10:00
11	Infinit e

19.11 Load 5 µL of sample plus 1 µL of staining solution 6x. Check the length of the band (bp).

Purify sample : PCR cleanup/ gel extraction according to manufacturer.

19.12 Digest the polymerized insert and the desired vector with 2U of appropriate restriction enzyme. See table "Primers" in order to choose the enzyme.



- 19.13 Digestion mixture ( 20 μL):
- Sfil 2 U
 - Plasmid ( 50 ng)
 - Insert ( 20 ng)
 -  2 μL CutSmart 10x (or appropriate buffer)
 - MiliQ to  20 μL
- 19.14 Incubate  01:00:00 at the temperature recommended by the manufacturer of the restriction enzyme
- 19.15 Inactivate the restriction enzyme by incubating 20 min  65 °C
- **Sfil can't be heat inactivated so in order to inactivate it a PCR clean up is necessary
- 19.16 Ligate the digested insert and vector.
Ligation mixture :
- 10 ul digestion product
 - 1 ul T4
 - 2 ul T4 ligase buffer
 - 7 ul MQ
- 19.17 Incubate  02:00:00  Room temperature
- 19.18 Thaw on ice  100 μL of electrocompetent cells
- 19.19 Add  2 μg of plasmid and the aliquot of electrocompetent cells to ice-cold electroporation cuvettes
- 19.20 Electroporate 25 uF, 10 kV/cm, 200 Ohms
- 19.21 Add 1 ml of BHIS and incubate  100 rpm, 30°C, 05:00:00 for recovery
- 19.22 Plate on LB + Agar (1.5%) + 100 $\mu\text{g}/\text{ml}$ spectinomycin.
For crispr cloning, to activate cas9 expression, add 0.2% mannose.



- 19.23 Incubate at 30 °C Overnight
- 19.24 Randomly pick potential mutants and test by colony PCR (for primers use the ones that flank the region of insertion sites).
- 19.25 After selecting colonies that show the expected band size, purify the desired plasmid (miniprep) and transform into *B. mycoides*.
- 20 Plate on LB agar containing erythromycin
Grow at 28 °C overnight
- 21 Pick colonies and test them by colony PCR
If the expected band size is present, purify the plasmid and sequence
- 22 Grow the colonies with the expected sequence on liquid LB + erythromycin + manose
Incubate 200 rpm, 28°C, 12:00:00
- 23 Next day pick colonies in BHI medium (without antibiotic) and grow at 37 °C overnight
- 24 Dilute the liquid cultures and plate them on LN (no antibiotic)
grow overnight at 37 °C