

Apr 09, 2020 Version 3

# Q5® Site-Directed Mutagenesis (E0552) V.3

DOI

[dx.doi.org/10.17504/protocols.io.beurjev6](https://dx.doi.org/10.17504/protocols.io.beurjev6)



New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs



**New England Biolabs**

New England Biolabs

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.beurjev6](https://dx.doi.org/10.17504/protocols.io.beurjev6)

External link: <https://www.neb.com/protocols/2014/03/21/q5-site-directed-mutagenesis-kit-without-competent-cells-protocol-e0552>

**Protocol Citation:** New England Biolabs 2020. Q5® Site-Directed Mutagenesis (E0552). [protocols.io](https://www.protocols.io)  
<https://dx.doi.org/10.17504/protocols.io.beurjev6>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 09, 2020

**Last Modified:** April 09, 2020


**Protocol Integer ID:** 35441

## Abstract


This protocol defines methods for the Q5® Site-Directed Mutagenesis Kit without competent cells.

## Materials

### MATERIALS

 Q5 Site-Directed Mutagenesis Kit (Without Competent Cells) - 10 rxns **New England**  
Biolabs Catalog #E0552S

## Safety warnings

 Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

## Exponential Amplification (PCR)

1 Assemble the following reagents in a thin-walled PCR tube:



	<b>25 <math>\mu</math>l RXN</b>	<b>FINAL CON C.</b>
Q5 Hot Start High-Fidelity 2X Master Mix	12.5 $\mu$ l	1X
10 $\mu$ M Forward Primer	1.25 $\mu$ l	0.5 $\mu$ M
10 $\mu$ M Reverse Primer	1.25 $\mu$ l	0.5 $\mu$ M
Template DNA (1–25 ng/ $\mu$ l)	1 $\mu$ l	1–25 ng
Nuclease-free water	9.0 $\mu$ l	

2 Mix reagents completely.



3 Transfer to a thermocycler and perform the following cycling conditions:  
Thermocycling Conditions for a Routine PCR:



<b>STEP</b>	<b>TEMP</b>	<b>TIME</b>
Initial Denaturation	98°C	30 seconds
25 Cycle	98°C	10 seconds

	s		nds
		50–72°C*	10–30 seconds
		72°C	20–30 seconds/kb
	Final Extension	72°C	2 minutes
	Hold	4–10°C	

Note

\* For a Q5-optimized annealing temperature of mutagenic primers, please use **NEBaseChanger™**, the online NEB primer design software. For pre-designed, back-to-back primer sets, a  $T_a = T_m + 3$  rule can be applied, but optimization may be necessary.

## Kinase, Ligase & DpnI (KLD) Treatment

4 Assemble the following reagents:



	VOLUME	FINAL CONC.
PCR Product	1 µl	
2X KLD Reaction Buffer	5 µl	1X
10X KLD Enzyme Mix	1 µl	1X
Nuclease-free Water	3 µl	



5 Mix well by pipetting up and down.



6 Incubate at  for .



## Transformation

7 Thaw  .

### Note

NEB 5-alpha Competent E. coli (High Efficiency), **NEB #C2987**, are recommended

8 Add  from the "KLD Section" above to the tube of thawed cells.



9 Carefully flick the tube 4-5 times to mix. **Do not vortex.**



10 Place the mixture  for .



11 Heat shock at  for .

12 Place  for .

13 Pipette  into the mixture.







14 Incubate at  for  with shaking (250 rpm).



15 Mix the cells thoroughly by flicking the tube and inverting.



16 Spread 50-100  $\mu\text{l}$  onto a selection plate. 

17 Incubate  Overnight at  37 °C . 

#### Note

It may be necessary (particularly for simple substitution and deletion experiments) to make a 10- to 40-fold dilution of the transformation mix in SOC prior to plating, to avoid a lawn of colonies.