


Aug 06, 2018

QA-Bio PNGase F Denatured Glycoprotein Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.sd6ea9e>



Mike Gibson¹

¹QA-Bio Inc



Mike Gibson

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.sd6ea9e>

External link: <https://www.qa-bio.com/product/pngase-f/>

Protocol Citation: Mike Gibson 2018. QA-Bio PNGase F Denatured Glycoprotein Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.sd6ea9e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 06, 2018

Last Modified: August 06, 2018

Protocol Integer ID: 14494

Keywords: PNGase F, oligosaccharides from glycoprotein, oligosaccharide, plant glycoprotein, deglycosylation, glycoprotein, high mannose glycan, truncated sugar molecule, alternative enzyme, glycan, enzyme, number of alternative enzyme, keeping protein, most native protein, proteins in solution, aspartic acid, charged sugar

Abstract

PNGase F cleaves N-linked (asparagine-linked) oligosaccharides from glycoproteins. The enzyme deaminates asparagine to aspartic acid, leaving the oligosaccharides intact. Denaturation increases the rate of cleavage. Most native proteins can still be completely N-deglycosylated but incubation time must be increased. The enzyme will remain fully active under reaction conditions (37°C) for at least 96 hours. PNGase F will not remove oligosaccharides containing Alpha-(1,3)-linked core fucose commonly found on plant glycoproteins; for this purpose, use peptide N-glycosidase A.

There are a number of alternative enzymes which can be used to remove N-glycans, most especially the Endo F family of enzymes and Endo H. These enzymes cleave between the two N-acetylglucosamine residues in the core of the oligosaccharide, generating a truncated sugar molecule with one N-acetylglucosamine residue remaining on the asparagine. This leaves a charged sugar which can assist in keeping proteins in solution that precipitate after deglycosylation with PNGase F which removes the oligosaccharide intact. Endo F1 cleaves high mannose and some hybrid type N-glycans. Endo F2 will removes biantennary and high mannose (at a 40X reduced rate). Endo F3 releases of triantennary and fucosylated biantennary N-glycans. Endo H removes hybrid or high mannose glycans.

Attachments



QA-Bio.E-

PNG01.specs...

271KB

Materials

MATERIALS

 PNGase F QA-Bio Inc Catalog #E-PNG01



- 1 Add up to 200 µg of glycoprotein to an Eppendorf tube.
- 2 Adjust to 35 µl final volume with de-ionized water.
- 3 Add 10 µl 5x Reaction Buffer 7.5 and 2.5 µl of Denaturation Solution.
- 4 Heat at 100°C for 5 minutes.
- 5 Cool.
- 6 Add 2.5 µl of Triton X-100 and mix.
- 7 Add 2.0 µl of enzyme to the reaction.
- 8 Incubate 3 hours at 37°C.