

Sep 04, 2020

Qbiotix protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.bks3kwgn>

vito¹

¹Qbiotix



vito

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bks3kwgn>

Protocol Citation: vito 2020. Qbiotix protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.bks3kwgn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol.

Created: September 04, 2020

Last Modified: September 04, 2020

Protocol Integer ID: 41531

Keywords: qbiotix protocol, protocol



- 1 1- collect saliva vial  1 mL
- 2 2- open vial and scan barcode
- 3 3- dispense  250 μ L sample to single well in plate with a pipette
- 4 4- dispense  250 μ L of magnetic bead solution in well with a pipette
- 5 5- incubate at controlled temperature on the plate reader  00:10:00
- 6 6- place plate on magnetic separator  00:00:30
- 7 7- insert plate on analyzer (reading speed depends on chosen plate reader instrument)
- 8 6- read results