

Jan 16, 2026

🌐 QCSOP25-2: Identification of Anti-Rhesus IgA by Direct and Indirect ELISA

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl41m2zlo5/v1>

Yasmine IACONE¹, Vimala Gollamudi¹, Diogo Magnani¹


¹NHPRR, UMass Chan Medical School, Worcester, MA (USA)

NHPRR



Yasmine IACONE

NHPRR - UMass Chan Medical School



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl41m2zlo5/v1>

External link: <http://nhpreagents.org>

Protocol Citation: Yasmine IACONE, Vimala Gollamudi, Diogo Magnani 2026. QCSOP25-2: Identification of Anti-Rhesus IgA by Direct and Indirect ELISA. [protocols.io https://dx.doi.org/10.17504/protocols.io.3byl41m2zlo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl41m2zlo5/v1)



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 21, 2025

Last Modified: January 16, 2026

Protocol Integer ID: 220641

Keywords: ELISA, Anti-IgA antibody, Rhesus Macaque, Nonhuman Primate (NHP), Immunoglobulin A (IgA), Mucosal Immunity, Monoclonal antibodies, Polyclonal Antibodies, Monomeric IgA, Dimeric IgA, NHPRR, Nonhuman Primate Reagent Resource, performing antibody identity testing, antibody identity testing for research, indirect elisa this elisa, standardized elisa method, antibody reagent, antibody specificity, antibody, indirect elisa, direct elisa, antigen specificity, immunoassay, nhprp reagent, identification of anti, reliable use of nhprp reagent, polyclonal goat, elisa, nhprp laboratory staff, nonhuman primate reagent resource,

Funders Acknowledgements:

NIH

Grant ID: AI126683

Abstract

This ELISA-based protocol confirms the identity and reactivity of **Nonhuman Primate Reagent Resource (NHPRR)** polyclonal **goat anti-rhesus IgA** and monoclonal **mouse anti-rhesus IgA** antibody reagents. Both **Indirect and Direct ELISA** formats are used to evaluate antibody binding specificity and lot-to-lot consistency. The procedure involves plate coating, blocking, sample incubation, detection, and data analysis, all performed according to defined acceptance criteria.

Objective:

To provide a standardized ELISA method for verifying the antigen specificity and identity of NHPRR anti-IgA antibodies.

Scope:

Applies to NHPRR laboratory staff performing antibody identity testing for research and quality control purposes.

Background:

Accurate identification of anti-IgA antibodies ensures reliable use of NHPRR reagents in immunoassays. The **Indirect ELISA** detects antibody binding through a secondary conjugate, while the **Direct ELISA** uses labeled anti-IgA for direct detection. Together, they confirm antibody specificity and support NHPRR's quality testing program.

Image Attribution

Created in BioRender. de Vasconcellos Castro, J. (2026) <https://BioRender.com/m2svlm2>

Materials

Standard ELISA materials and reagents:

Item	Purpose	Vendor	Catalog Number
TMB Substrate	Substrate for HRP	1-Step Ultra TMB ELISA Substrate Solution, Thermo Fisher Scientific	34029
Corning® 96-well Clear Flat Bottom Polystyrene High Bind Microplate, Nonsterile	High-binding plates for ELISA assay	Corning	9018
Gibco™ [1X] DPBS, calcium, magnesium, glucose, pyruvate	Diluent for reagents used for coating the working plates	Gibco™	14-287-080
Gibco™ PBS (10X), pH 7.2	Used in the preparation of the Washing Buffer	Gibco™	70-013-032
VWR® High-Performance [15 mL] Centrifuge Tubes with Flat or Plug Caps, PP	Used for the preparation/dilution of the test sample and other reagents	VWR or equivalent	89039-668
ELISA Stop Solution	Stops the enzymatic reaction	Invitrogen™	SS04
Non-fat Dried Milk	Used for the Blocking Buffer solution	LabScientific Biotin-free non-fat dry milk	M0841
Low-binding plates	Dilution plate	Of your choice	Of your choice
1-Step™ TMB ELISA Substrate Solutions	Acts as the substrate for the HRP enzyme	Thermo Fisher Scientific	34029
TWEEN® 20	Used in the preparation of the Washing Buffer	Sigma Aldrich	P6585-100ML
Fisherbrand™ Hexagonal Polystyrene Weighing Dishes	Used for weighing dry milk during Blocking Buffer preparation	Fisher Scientific	02-202-103
Timer	To keep track of incubation and developing time	of your choice	N/A
Rainin Classic™ Manual Single-Channel Pipettes and Allocated Pipette Tips	Used for the preparation, dilution, and transfer of the test	Rainin or equivalent	N/A



Item	Purpose	Vendor	Catalog Number
	sample and other reagents		
Pipette Tips SR LTS 200µL F 960A/5	Used for transferring the test sample and other reagents from the dilution plate to the working plate	Rainin	17005859
PPE (Personal Protective Equipment)	For analyst protection: lab coat, nitrile gloves, and safety goggles	of your choice	N/A
Plate Sealers	To cover ELISA plates for incubation	BioLegend	423601
Elisa Plate Reader Software	Analyze the OD values	Molecular Devices	Model: SoftMax® Pro
Nalgene™ Rapid-Flow™ Sterile Filter Storage Bottles	Used for the preparation of the Blocking Buffer	Thermo Fisher Scientific	455-0250
Pipet-Lite™ XLS+ Manual 12-Channel Pipette, 20-200 µL	Used for transferring the test sample and other reagents from the dilution plate to the working plate	Rainin	17013810
Sterile GRIPTIPS® INTEGRA 1250 µL Pipette Tips	Used for the preparation, dilution, and transfer of the test sample and other reagents	INTEGRA	6544 (non-filtered) 6545 (filtered)
ELISA Plate Washer	Wash ELISA plates	BioTek	Model #: 405LSRS
ELISA Microplate Reader	Measure the OD values	Molecular Devices	Model: SpectraMax ABS
INTEGRA VIAFLO 12-Channel 50-1250 µL Electronic Pipette	Used for the preparation, dilution, and transfer of the test sample and other reagents	INTEGRA	Model: 12-channel VIAFLO
[Serological] Pipette, 10 mL and 5 mL, Graduated 1/10 mL,	Used for the preparation/dilutio	Greiner Bio-One	606180 (5 ml) 607107 (10 mL)



Item	Purpose	Vendor	Catalog Number
Sterile	n of the test sample and other reagents		
Falcon® 96-well Polystyrene [Dilution] Microplates	Low-binding plates for preparation of the test sample and other reagents	Corning	353936
Calculator	Calculate the dilution volumes	of your choice	N/A
Biohazard Waste Bags	Collect used pipette tips and other waste from the experiment	of your choice	N/A
Analytical Balance	For weighing dry milk (part of the Blocking Buffer preparation)	of your choice	N/A
VWR® [50 mL] Disposable Pipetting Reservoirs	For holding the test sample and other reagent solutions before transfer	VWR or equivalent	89094-684

Assay-Specific ELISA reagents:

 Reference rhesus monomeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00377**

 Reference rhesus dimeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00378**

 Goat Anti-Rhesus IgA [GARI] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00366**

 Anti-IgA [NHPMab001] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00423**

 Anti-IgA [NHPMab001]-HRP **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00392**

 Goat Anti-Rhesus IgA [GARI]-HRP **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00365**

 Mouse Anti-Goat IgG Fc-HRP (SB115d) **Southern Biotech Catalog #6158-05**

 Peroxidase AffiniPure® F(ab)₂ Fragment Goat Anti-Mouse IgG, Fcγ fragment specific **Jackson ImmunoResearch Laboratories, Inc. Catalog #115-036-071**



	Antibody	Isotype	Species	Vendor	Catalog Number	Purpose	RRID
	Reference rhesus monomeric IgA [CC6.29]	IgA	rhesus	NHPRR	CAT-00377	Coating	AB_3714921
	Reference rhesus dimeric IgA [CC6.29]	IgA	rhesus	NHPRR	CAT-00378	Coating	AB_3714922
	Goat Anti-Rhesus IgA [GARI]	IgG	goat	NHPRR	CAT-00366	Primary Antibody	AB_3086850
	Anti-IgA [NHPMab001]	IgG	mouse	NHPRR	CAT-00423	Primary Antibody	AB_3717991
	Mouse Anti-Goat Fc-HRP (SB115d)	IgG	mouse	Southern Biotech	6158-05	Secondary Antibody	AB_2796222
	Peroxidase AffiniPure® F(ab') ₂ Fragment Goat Anti-Mouse IgG, Fcy fragment specific	IgG	goat	Jackson ImmunoResearch	115-036-071	Secondary Antibody	AB_2338524
	Goat Anti-Rhesus IgA [GARI]-HRP	IgG	goat	NHPRR	CAT-00365	Primary Antibody already conjugated	AB_3086849
	Anti-IgA [NHPMab001]-HRP	IgG	mouse	NHPRR	CAT-00392	Primary Antibody already conjugated	AB_3105907



Protocol materials

⊗ Mouse Anti-Goat IgG Fc-HRP (SB115d) **Southern Biotech Catalog #6158-05**

⊗ Peroxidase AffiniPure® F(ab)₂ Fragment Goat Anti-Mouse IgG, Fcγ fragment specific **Jackson ImmunoResearch Laboratories, Inc. Catalog #115-036-071**

⊗ Reference rhesus dimeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00378**

⊗ Anti-IgA [NHPMab001] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00423**

⊗ Goat Anti-Rhesus IgA [GARI]-HRP **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00365**

⊗ Reference rhesus monomeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00377**

⊗ Goat Anti-Rhesus IgA [GARI] **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00366**

⊗ Anti-IgA [NHPMab001]-HRP **Nonhuman Primate Reagent Resource (NHPRR) Catalog #CAT-00392**



Troubleshooting

Problem

Blank Wells Show High Signal

Solution

The issue may be caused by incomplete washing, insufficient blocking, incorrect buffer volume, contamination of the washing buffer, or an inadequate number of wash cycles. To address this, repeat the assay using freshly prepared blocking reagent and ensure reagent reservoirs are clean or replaced. Confirm that all wash steps are performed using the correct buffer and that full plate coverage is achieved during the blocking step. Use new pipette tips for all reagent transfers to prevent cross-contamination.

Problem

All Wells Show Low Signal

Solution

The issue may be due to inactive streptavidin-HRP or TMB substrate, incorrect incubation temperature, or systematic pipetting errors. To address this, verify reagent activity by testing streptavidin-HRP and TMB substrate using a known positive control. Additionally, confirm that all incubation steps were performed at room temperature or in accordance with SOP requirements. Inspect the performance of the plate reader to ensure proper function. If necessary, repeat the assay using a new high-binding ELISA plate and freshly prepared signal amplification reagents from different lots or new aliquots.

Problem

High Variability Between Duplicates

Solution

The issue may be caused by inconsistent pipetting, incomplete mixing of reagents, or non-uniform loading of samples and reagents. To address this, use calibrated pipettes and maintain a consistent pipetting technique; the use of a multichannel pipette is recommended for uniform application. Additionally, ensure that all samples and reagents are thoroughly mixed before use.



Before start

Note on Concentration Units: All antigen and antibody dilutions in this protocol are prepared using molar concentrations (e.g., nM) to ensure accurate comparison across immunoglobulin types with different molecular weights. When applicable, molecular weight values are used to convert between ng/mL and nM to reflect molar equivalency (e.g., IgG and monomeric IgA ~150 kDa, dimeric IgA ~300 kDa, IgM ~950 kDa).

Concentrations are lot-dependent; always make sure to check the lot even if the catalog number is the same.

Example conversion: To convert 1 mg/mL of monomeric IgA (~150 kDa) to molarity (mM), ensure that the units are consistent: mass in **mg/mL** and molecular weight in **kDa**. If using µg/mL, your result will be in µM instead of mM.

$$(1 \text{ mg/mL}) / (150 \text{ kDa}) = 0.00667 \text{ mM} = 6.67 \text{ µM}$$

$$\text{Alternatively, } [(1 \text{ mg/mL}) / (150 \text{ kDa})] \times 1000 = 6.67 \text{ µM}$$

Prepare all the calculations beforehand.

Definitions and Acronyms

1 Definitions

Analyte / Antigen (recombinant rhesus IgA, monomeric and dimeric)

Recombinant rhesus IgA, present in either monomeric or dimeric form, is coated onto the ELISA plate and serves as the target analyte. Anti-rhesus IgA antibodies in the sample bind specifically to the immobilized IgA. The assay signal generated reflects the extent of antibody binding to the coated antigen.

Test Articles (TAs), TA-pAb, and TA-mAb

Test Articles (TAs): Anti-rhesus IgA antibody reagents that are being evaluated for identity and specificity.

TA-pAb: Goat anti-rhesus IgA polyclonal antibody Test Article.

TA-mAb: Mouse anti-rhesus IgA monoclonal antibody Test Article.

These reagents are tested in parallel on plates coated with monomeric and dimeric IgA to assess binding profiles and confirm identity.

Reference Control (Positive Control) Antibody

A previously qualified anti-rhesus IgA antibody lot which is included on each plate (typically as a serial dilution/titration curve) to confirm that the assay produces the expected binding behavior (signal range and curve shape) and to support run-to-run trending and lot-to-lot comparability of reagents. (Abcam, 2024). In this ELISA, the positive control is an anti-rhesus IgA antibody sample that is known to specifically bind both monomeric and dimeric recombinant rhesus IgA. Its binding curve serves as a benchmark for the expected signal and plate validity.

Note

If the positive control fails to generate the expected signal, the assay may be invalid, and the run should be repeated (BosterBio, n.d.).

Negative Control Antibody

An antibody of the same species, isotype, and format as the Test Articles (TAs), but known **not** to bind rhesus IgA. The negative control antibody should produce only baseline signal. It is used to verify assay specificity and to help identify non-specific binding or background issues.

Blank

A **buffer-only well** that is used to establish the background signal and verify the absence of non-specific binding. Unexpected signal indicates contamination or assay

error.

Note

If a blank well exhibits an unexpected signal, this may indicate contamination or non-specific binding, and the assay should be considered invalid (Jackson ImmunoResearch, n.d.).

Monoclonal Antibody (mAb)

An antibody preparation derived from a single B-cell clone, recognizing one specific epitope on an antigen and therefore having a single binding specificity. In this SOP, the monoclonal antibody Test Article (TA-mAb) is a mouse anti-rhesus IgA antibody evaluated for its binding to recombinant monomeric and dimeric rhesus IgA.

Polyclonal Antibody (pAb)

An antibody preparation containing a mixture of immunoglobulin molecules produced by multiple B-cell clones, recognizing multiple epitopes on the same antigen. In this SOP, the polyclonal antibody Test Article (TA-pAb) is a goat anti-rhesus IgA antibody evaluated for its binding profile to recombinant monomeric and dimeric rhesus IgA.

Monomeric IgA

A single IgA immunoglobulin unit composed of two heavy and two light chains (one “Y-shaped” IgA molecule). In this SOP, recombinant monomeric rhesus IgA is used as a coated antigen on ELISA plates to assess antibody binding to the monomeric form of IgA.

Dimeric IgA

A form of IgA composed of two IgA monomers covalently linked, typically via a joining (J) chain, representing the main configuration of secretory IgA at mucosal surfaces. In this SOP, recombinant dimeric rhesus IgA is used as a coated antigen on ELISA plates to assess antibody binding to the dimeric form of IgA.

Direct ELISA

A plate-based immunoassay format in which the antigen is immobilized on the plate and detected using a single enzyme-labeled (e.g. HRP-conjugated) primary antibody that binds directly to the antigen. Signal is generated when the enzyme acts on a chromogenic substrate (e.g. TMB), and the measured absorbance is proportional to the amount of antigen–antibody complex formed. Direct ELISA uses fewer steps and reagents than indirect ELISA but offers less signal amplification because only one labeled antibody binds per antigen site. (*ELISA Guide; Part 1: Introduction to ELISA, Formats and Signal Amplification*, 2023; NCBI, 2020)

Indirect ELISA

A plate-based immunoassay format in which the antigen is immobilized on the plate,

followed by incubation with an unlabeled primary antibody specific for the antigen, and then a labeled secondary antibody (e.g. HRP-conjugated anti-species antibody) that recognizes the primary antibody. The enzyme on the secondary antibody converts a chromogenic substrate (e.g. TMB) to a colored product measured by absorbance. Indirect ELISA provides signal amplification because multiple secondary antibodies can bind each primary antibody, increasing assay sensitivity compared with direct ELISA. (*ELISA Guide; Part 1: Introduction to ELISA, Formats and Signal Amplification*, 2023; NCBI, 2020)

Dilution Plate

A low-binding microplate used to prepare antibody dilutions and perform vertical serial dilutions of TAs, controls, and HRP-conjugated antibodies before transfer to the working plate. The dilution plate mirrors the layout of the working plate but is not read in the plate reader.

Working Plate

The high-binding ELISA plate on which antigen coating, blocking, antibody incubation, washing, substrate addition, and OD reading at 450 nm are performed. The working plate contains the bound analyte and is the plate that generates the analytical signal.

Blocking Buffer

In this SOP, the blocking buffer is 5% non-fat dry milk in PBS with 0.05% Tween-20. It is applied after the coating step to occupy unbound sites on the plate surface and is also used as a diluent for antibodies. Blocking buffer minimizes non-specific binding and reduces background signal. (Excedr, 2023)

Washing Buffer

PBS containing 0.05% Tween-20 (1× DPBS or PBS base, as prepared under “Wash Buffer Preparation”). It is used between assay steps to remove unbound reagents and reduce non-specific background.

Detection Reagent (HRP-conjugated antibody)

For the indirect ELISA, the detection reagents are HRP-conjugated secondary antibodies (mouse anti-goat HRP and goat anti-mouse HRP) that bind the species-specific primary anti-IgA antibodies. For the direct ELISA, the detection reagents are HRP-conjugated anti-rhesus IgA antibodies themselves. In both formats, the detection reagent binds to the antibody-antigen complex and enables signal generation via HRP-mediated TMB conversion. (*ELISA Guide; Part 1: Introduction to ELISA, Formats and Signal Amplification*, 2023; NCBI, 2020)

TMB Substrate and Stop Solution

TMB (3,3',5,5'-Tetramethylbenzidine) is a chromogenic substrate that is oxidized by HRP, producing a blue color proportional to bound HRP. The Stop solution (acid) halts

the HRP reaction and converts the color to yellow, allowing stable measurement at 450 nm.

Optical Density at 450 nm (OD_{450})

The absorbance measured at 450 nm using an ELISA plate reader after TMB development and stopping. OD_{450} values reflect the amount of HRP-mediated color change and are used to construct binding curves for each antibody.

Acronyms:

CV – Coefficient of Variation.

DPBS – Dulbecco's Phosphate-Buffered Saline.

EC₅₀ – Half-maximal effective concentration.

ELISA – Enzyme-Linked Immunosorbent Assay.

HRP – Horseradish Peroxidase.

IgA – Immunoglobulin A.

kDa – KiloDalton (molecular mass unit).

LoD – Limit of Detection.

NHP – Nonhuman Primate.

NHPRR – Nonhuman Primate Reagent Resource.

OD – Optical Density (absorbance).

O/N – Overnight (incubation at 4 °C unless otherwise stated).

PBS – Phosphate-Buffered Saline.

RT – Room Temperature.

TA – Test Article.

TA-pAb – Test Article polyclonal anti-rhesus IgA antibody (goat).

TA-mAb – Test Article monoclonal anti-rhesus IgA antibody (mouse).

TMB – 3,3',5,5'-Tetramethylbenzidine.

Principle

- 2 **Analyte/Antigen:** recombinant **rhesus IgA**, monomeric and dimeric forms (two plates processed in parallel).

Test Articles (TA):

TA-pAb: goat anti-rhesus IgA (polyclonal)

TA-mAb: mouse anti-rhesus IgA (monoclonal)

Controls:

Negative Control Ab: same species/format as TA, known non-binder to rhesus IgA

Reference Positive Control Ab: previously qualified anti-IgA (optional but recommended for plate validity)

Identity is confirmed by a specific binding signature to both monomeric and dimeric IgA with acceptable curve fit and specificity, and by matching expected relative response with defined limits.

See Acceptance Criteria Section below.




Figure 1: Workflow for Identification of Anti-Rhesus IgA by ELISA

Publication citation: Created in BioRender. de Vasconcellos Castro, J. (2026)

<https://BioRender.com/m2svlm2>

Concentration & Dilution Notes

3 Please read the guidelines on concentration units under Guidelines & Warnings.

Use the dilution formula:

$$C_1V_1 = C_2V_2$$

Where:

- C_1 = stock concentration
- V_1 = volume of stock needed
- C_2 = desired final concentration
- V_2 = final total volume

Example: To prepare 1 mL of antigen at 1 μ M from a 1 mM stock:

- C_1 = 1 mM = 1000 μ M
- C_2 = 1 μ M



- $V_2 = 1 \text{ mL} = 1000 \text{ }\mu\text{L}$
- $V_1 = (C_2 \times V_2) / C_1 = 1 \text{ }\mu\text{M} \times 1000 \text{ }\mu\text{L} / 1000 \text{ }\mu\text{M} = 1 \text{ }\mu\text{L}$

So: **Add 1 μL of stock + 999 μL of DPBS** to reach 1 μM in 1 mL total volume.

Wash Buffer Preparation (1x DPBS + 0.05% Tween-20)

4 Add  18 L of deionized water in a container.

4.1 Add  2 L of 10x DPBS and mix gently.

4.2 Add  50 mL of Tween-20 and mix gently.

Note

Scale up or down according to the need.

Day 1 - Plate Coating (~30 min preparation + O/N incubation)

1d

5 Prepare  10 nanomolar (nM) (C_2) of IgA in 1x DPBS.

Two working plates in parallel:

- Plate 1: **Monomeric IgA**
- Plate 2: **Dimeric IgA**

*Plate Layout (See **Figure 2**):*

- Columns 1-9: Coating with IgA (as per plate type);
- A10-A11: Coat only (IgA);
- B10-F11: 1x DPBS only controls;
- G10-G11: Blank (1x DPBS)




Figure 2: Coating Plate Layout

Publication citation: Created in BioRender. de Vasconcellos Castro, J. (2026)
<https://BioRender.com/rh7lu1d>

Calculation Volumes (in Table 1)

$V_2 = 74 \text{ wells/plate} \times 100 \text{ }\mu\text{L} = 7400 \text{ }\mu\text{L}$; prepare 7800 μL /plate to allow excess.

 Reference rhesus monomeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR)** Catalog #CAT-00377

 Reference rhesus dimeric IgA [CC6.29] **Nonhuman Primate Reagent Resource (NHPRR)** Catalog #CAT-00378

	Antibody	Initial Concentration (C1)	Final Concentration (C2)	Final Volume (V2)	Initial Volume (V1)
	Referenc e rhesus monome ric IgA	2.69 mg/mL, 146.6 kDa, 18.3 μM	10 nM	7800 μL	4.3 μL
	Referenc e rhesus dimeric IgA	0.56 mg/mL, 308.7 kDa, 1.8 μM	10 nM	7800 μL	43.3 μL

Table 1: Calculation for the coating solution.

**Note**

For accuracy, remove the equivalent DPBS volume and add the same volume of stock.

5.1 Add  100 μ L /well.

5.2 Seal the plate/s and incubate  Overnight at  4 °C .

Day 2 - Blocking (~15 min preparation + 1 h incubation)**1h 15m**

6 Wash the working plates twice with wash buffer.

2m

6.1 Prepare the **Blocking buffer** with  5 % volume milk in wash buffer:

9m

For  150 mL : add  7.5 g milk to  150 mL PBS +
 0.05 % volume Tween-20 .

6.2 Add  200 μ L /well.

3m

6.3 Seal the plate/s.

1m

6.4 Incubate  01:00:00 at  Room temperature

1h

Note

Potential STOP point: store up to 1-2 days at  4 °C in blocking buffer

OPTION A - Indirect ELISA: Primary antibody preparation and incubation (~30 min preparation + 1 hour incubation)

1h 30m

- 7 Prepare [IM] 40 nanomolar (nM) (C₂) working stocks of **TA-pAb**, **TA-mAb**, and **Negative Control Ab** in blocking buffer.

Plate Layout (Figure 3):

- Columns 1-3: TA-pAb (goat anti-rhesus IgA, pAb);
- Columns 4-6: TA-mAb (mouse anti-rhesus IgA, mAb);
- Columns 7-9: Negative control Ab (non-binder);
- Row A10-11: Blocking buffer (coated wells only);
- Row B10-11: TA-pAb controls;
- Row C10-11: TA-mAb controls;
- Row D10-11: Negative control Ab controls;
- Row E10-F11: Reserved for secondary-only controls;
- Row G10-11: Blank.




Figure 3: Primary Antibodies Plate Layout

Publication citation: Created in BioRender. de Vasconcellos Castro, J. (2026)

<https://BioRender.com/kvj3gee>

Volume Calculations V_2 (in Table 2):

5 wells/plate/each antibody (3 wells for the triplicate + 2 wells for the control)

2 plates total

200 μL of primary antibody solution/well

$V_2 = 5 \times 2 \times 200 = 2000 \mu\text{L}$ - Make it in excess

$V_2 = 2500 \mu\text{L}$

 Goat Anti-Rhesus IgA [GARI] **Nonhuman Primate Reagent Resource**
(NHPRR) Catalog #CAT-00366

 Anti-IgA [NHPMab001] **Nonhuman Primate Reagent Resource**
(NHPRR) Catalog #CAT-00423

	Antibody	Initial Concentration (C1)	Final Concentration (C2)	Final Volume (V2)	Initial Volume (V1)
	Goat anti-rhesus IgA [GARI]	6.64 mg/mL, 150 kDa, 44 μM	40 nM	2500 μL	2.272 μL
	Mouse anti-rhesus IgA [NHPMab001]	3.78 mg/mL, 150 kDa, 25 μM	40 nM	2500 μL	4 μL

Table 2: Calculation for the working stock solution.

Note

For accuracy, remove the equivalent blocking buffer volume and add the same volume of stock.

7.1 In a **dilution plate**, dispense  200 μL of working stock to row A of designated wells (**Figure 3**). Control wells will also receive the working stock.

Note

- The dilution plate allows for preparing the dilutions ahead of time, while the working plate is for blocking. Also, in case of mistakes, you can toss the dilution plate without having to repeat the whole experiment.
- The dilution plate is an exact same copy of the working plate.

7.2 Add  150 μL of blocking buffer to the remaining wells.

7.3 Perform 1:4 vertical serial dilutions( 50 μL transfers, mix 5x, discard the very last  50 μL). Ensure equal volumes visually.

See **Figure 4** for reference.




Figure 4. Schematic representation of 1:4 vertical serial dilution.

Publication Citation: Created in BioRender. de Vasconcellos Castro, J. (2026)

<https://BioRender.com/vc899id>

7.4 Wash the working plate twice.

7.5 Transfer  100 μL /well from the dilution plate to the working plate.

7.6 Seal the plate and incubate for  01:00:00 at  Room temperature .

1h

OPTION B - Direct ELISA: Primary antibody preparation and incubation (~30 min preparation + 1 hour incubation)

1h 30m

8 Use the same coated and blocked plates as in **Section 6**.

8.1 Prepare HRP-conjugated TA-pAb and TA-mAb at  40 nanomolar (nM) working stocks in blocking buffer.

Note

The layout of the working plate is identical to that in **Section 7**.

Volume Calculations V_2 (in Table 2):

5 wells/plate/each antibody (3 wells for the triplicate + 2 wells for the control)

2 plates total

200 μ L of primary antibody solution/well

$V_2 = 5 \times 2 \times 200 = 2000 \mu\text{L}$ - Make it in excess

$V_2 = 2500 \mu\text{L}$



Goat Anti-Rhesus IgA [GARI]-HRP **Nonhuman Primate Reagent Resource (NHPRR)** Catalog #CAT-00365



Anti-IgA [NHPMab001]-HRP **Nonhuman Primate Reagent Resource (NHPRR)** Catalog #CAT-00392

	Antibody	Initial Concentration (C1)	Final Concentration (C2)	Final Volume (V2)	Initial Volume (V1)
	HRP Goat anti-rhesus IgA [GARI]	0.5 mg/mL, 150 kDa, 3.3 μ M	40 nM	2500 μ L	30 μ L
	HRPMouse anti-rhesus IgA	0.5 mg/mL, 150 kDa, 3.3 μ M	40 nM	2500 μ L	30 μ L

	Antibody	Initial Concentration (C1)	Final Concentration (C2)	Final Volume (V2)	Initial Volume (V1)
	[NHPMa b001]				

Table 2: Calculation for the working stock solution.

Note

For accuracy, remove the equivalent blocking buffer volume and add the same volume of stock.

8.2 Set up the dilution plate as in **sections 7.1, 7.2, and 7.3**.

8.3 Wash the working plate twice.

8.4 Transfer  100 µL /well from the dilution plate to the working plate.

8.5 Seal the plate and incubate for  01:00:00 at  Room temperature .

Note

For Direct ELISA, **skip** the secondary antibody preparation and incubation in **Section 9**, and go directly to Detection & Reading (**Section 10**).

Secondary antibody preparation and incubation (~20 min preparation + 1 hour incubation)

1h 20m

- 9 Prepare each HRP secondary antibody at  1 nanomolar (nM) in blocking buffer:
- mouse anti-goat HRP (detects TA-pAb)
 - goat anti-mouse HRP (detects TA-mAb)
 - Assign secondary for the Negative control Ab species as applicable.

Plate layout (Figure 5):

- Columns 1-3: HRP anti-TA-pAb, mouse anti-goat HRP;
- Columns 4-6: HRP anti-TA-mAb, goat anti-mouse HRP;
- Columns 7-9: HRP anti-negative control Ab, depending on the control antibody of your choice;
- Row A10-11: Blocking buffer (coat only);
- Row B10-11: Blocking buffer (TA-pAb controls);
- Row C10-11: Blocking Buffer (TA-mAb controls);
- Row D10-11: Blocking buffer (Negative control Ab controls);
- Row E10-F11: HRP secondary antibodies controls;
- RowG10-11: Blocking buffer, blank.




Figure 5: Secondary Antibody Plate Layout

Publication Citation: Created in BioRender. de Vasconcellos Castro, J. (2026)
<https://BioRender.com/n5v2prn>

Volume calculations V_2 (in Table 3):

26 well/plate will receive the secondary antibody

2 plates total

100 μL of coating solution/well

$V_2 = 26 \times 2 \times 100 = 5200 \mu\text{L}$ - Make it in excess

$V_2 = 5600 \mu\text{L}$

⚗ Mouse Anti-Goat IgG Fc-HRP (SB115d) **Southern Biotech Catalog #6158-05**

⚗ Peroxidase AffiniPure® F(ab)₂ Fragment Goat Anti-Mouse IgG, Fcy fragment specific **Jackson ImmunoResearch Laboratories, Inc. Catalog #115-036-071**

	Antibody	Initial Concentration (C1)	Final Concentration (C2)	Final Volume (V2)	Initial Volume (V1)
	Mouse anti-goat IgG Fc-HRP (SB115d)	0.75 mg/mL, 150 kDa 5 μM	1 μM	5600 μL	1.12 μL
	HRP F(ab) ₂ Fragment Goat Anti-Mouse IgG	0.6 mg/mL, 100 kDa (it's a F(ab') ₂) = 6 μM	1 μM	5600 μL	0.933 μL

Table 3: Calculation for the secondary antibody solution.

Note

Depending on the Control Antibody of your choice, you will need to adjust the volume of secondary antibody (V₂).

Note

For accuracy, remove the equivalent blocking buffer volume and add the same volume of stock.

- 9.1 Wash working plates twice in wash buffer.
- 9.2 Add  100 μL /well.
- 9.3 Seal the plate and incubate for  01:00:00 at  Room temperature .

1h



Detection & Reading (~15 min)

15m

10 Wash the working plates twice with wash buffer.

10.1 Add  100 μ L /well of TMB substrate.

10.2 Incubate for 5-10 min.

Note

Monitor the color development as this can vary depending on the antibodies.

10.3 Stop the reaction with  100 μ L /well ELISA Stop solution.

10.4 Read the plates at 450 nm.

Data Extraction and Organization

11 Export OD₄₅₀ readings from the plate reader as a ".csv" or ".xlsx" file.

11.1 Option A - Excel

- Organize the data with columns for concentration (nM) and rows for the TAs (raw OD₄₅₀ values);
- Average the triplicates and calculate %CV;

Note

%CV = Coefficient of Variation, measures the variability among replicates. Lower values indicate precision between replicates. For ELISA assays $\leq 15\%$ is generally acceptable.

Calculate it as follows:

$\%CV = (\text{Standard Deviation, SD} / \text{Mean}) \times 100$

- Save the table and import into GraphPad Prism or other analysis software;
- Go to **Section 12**.

11.2 Option B - GraphPad Prism



- Open GraphPad Prism and create a new XY data table ($X = \text{Concentration}$, $Y = OD_{450}$ readings);
- Under X, enter the antibody concentrations and label the column "Concentration (nM)";
- Under Y, paste the corresponding raw OD_{450} values in three replicate subcolumns for each sample;
- Prism will automatically generate a graph, under the "Graphs" tab;
- Adjust the graph layout as desired (color, font, size, legends, etc...);
- Since concentrations are expressed in linear values, set the X-axis to logarithmic scale: Double-click the X-axis → under Scale, select "Log 10".

11.3 Now you are ready to analyze the samples.

Data Analysis & Interpretation

12 In Prism, choose **Analyze → Nonlinear regression (curve fit)**

12.1 Select "**Dose-response - sigmoidal, 4PL, X is the concentration**"

12.2 Under **Settings → Constraints**, leave unconstrained unless you have historical limits already set.

12.3 Outputs per curve:

- R^2
- EC_{50}
- **Top (maximum OD)**
- **Bottom (baseline OD)**
- **Hill slope**

Acceptance Criteria

13 The acceptance criteria will be the following:

	Parameter	Acceptance Criterion	Rationale
	Goodness of fit (R^2)	≥ 0.98	Indicates strong fit of the 4PL model to the experimental data

Parameter	Acceptance Criterion	Rationale
Top (maximum OD)	within 1.0 and 4.0 OD units (depends on the LoD of the reader)	Ensures the signal is within the reader's linear range and avoids saturation
Bottom (baseline)	≤ 0.10 OD	Confirms minimal background and proper washing/blocking
Hill Slope	Between 0.8 and 2.0	Reflects a good dose-response curve steepness. Values outside may indicate hook effect or non specific binding
EC50 (half-maximal concentration)	within 0.5 to 2 x of an established range or a reference lot if available	Ensures the antibody affinity is consistent with its identity
Replicate precision (%CV)	≤ 15 %	Confirms intra-assay repeatability and accurate pipetting
Specificity check (Negative control Ab)	≤ 0.10 OD	Verifies that binding is specific to IgA and not background or nonspecific interactions
Hook effect check	No drop in OD at highest concentrations	Ensures no antigen excess. You can exclude the point and refit the curve

Summary of acceptance criteria

- 14 **The test articles (antibodies) identity is confirmed when the 4PL regression meets the above criteria on both monomeric and dimeric IgA plates.**

Protocol references

Abcam. (2024). *Required controls for ELISA*. Abcam.com. <https://www.abcam.com/en-us/technical-resources/guides/elisa-guide/controls-for-elisa>

BosterBio. (n.d.). *How to design positive and negative controls for ELISA*. <https://www.bosterbio.com/how-to-design-positive-negative-controls-ihc-western-blot-elisa>

de Vasconcellos Castro, J. (2026) <https://BioRender.com/m2svlm2>

de Vasconcellos Castro, J. (2026) <https://BioRender.com/rh7lu1d>

de Vasconcellos Castro, J. (2026) <https://BioRender.com/kvj3gee>

de Vasconcellos Castro, J. (2026) <https://BioRender.com/vc899id>

de Vasconcellos Castro, J. (2026) <https://BioRender.com/n5v2prn>

ELISA Guide; Part 1: Introduction to ELISA, Formats and Signal Amplification. (2023, November 30). <https://www.jacksonimmuno.com/secondary-antibody-resource/immuno-techniques/elisa-guide-part-1/>

Jackson ImmunoResearch. (n.d.). *ELISA Guide Part 3: ELISA Optimization*. <https://www.jacksonimmuno.com/secondary-antibody-resource/immuno-techniques/elisa-guide-part-3-elisa-optimization/>

Acknowledgements

This protocol and reagents were developed by the Nonhuman Primate Reagent Resource (NHPRR; RRID:SCR_012986; nhpreagents.org)