



Apr 07, 2023

Version 1

QIAGEN DNeasy PowerMax Soil Kit V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr4ke3gmk/v1>

QIAGEN¹

¹QIAGEN



Yuanyu Cheng

McGill University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4ke3gmk/v1>

Protocol Citation: QIAGEN 2023. QIAGEN DNeasy PowerMax Soil Kit. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr4ke3gmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Protocol used successfully by Sales et al., 2021 to detect fish sedDNA

Created: December 21, 2022



Last Modified: April 19, 2023

Protocol Integer ID: 74334

Keywords: QIAGEN, DNeasy, PowerMax, SedDNA, Sedimentary DNA, qiagen dneasy powermax soil kit for the isolation, qiagen dneasy powermax soil kit, microbial dna, microbial dna from large quantity, genomic dna from environmental sample, isolated dna, extracted genomic dna, dneasy powermax soil kit, isolating genomic dna, specific for bacteria, low microbial load, genomic dna, successful pcr amplification from sample, bacillus anthracis, bacillus subtili, total dna, dna exposure, dna, pcr amplification, successful pcr amplification, pcr with primer, soil, using pcr, genomic, various soil type, bacteria, fungi, environmental sample, soil type, streptomyce, actinomycete, qiagen, sediment sample

Abstract

For the isolation of microbial DNA from large quantities of soil - great for samples with low microbial load

The DNeasy PowerMax Soil Kit comprises a novel and proprietary method for isolating genomic DNA from environmental samples using Inhibitor Removal Technology® (IRT). With this kit, it is possible to process samples that have proven difficult in the past due to high levels of humic-like substances. The isolated DNA has a high level of purity, which allows for successful PCR amplification from samples. Total DNA isolated from various soil types has been successfully amplified using PCR with primers specific for bacteria (*Bacillus subtilis*, *Bacillus anthracis*), fungi (yeast, mold), and actinomycetes (*Streptomyces*).

Using the DNeasy PowerMax Soil Kit, environmental samples are added to a bead-beating tube with a kit-supplied proprietary buffer for rapid and thorough homogenization. Cell lysis and DNA exposure occur by mechanical and chemical methods. Extracted genomic DNA is captured on a silica membrane in a spin column format. The DNA is washed and eluted from the membrane and is ready for PCR and other downstream applications.

Protocol successfully used by Sales et al., 2021 to characterize fish species diversity and richness in water and sediment samples

Materials

DNeasy PowerMax Soil Kit (10 preparations) already includes:

MB Maxi Spin Columns * 10

PowerMax Bead Pro Tubes * 10

PowerBead Solution  20 mL

Solution C1  6.6 mL * 2

Solution C2  28 mL * 2

Solution C3.  44 mL

Solution C4  330 mL

Solution C5  30 mL * 4

Solution C6  66 mL

Collection Tubes (50ml) 5 * 8

Quick Start Protocol * 1

Additional equipments and reagents to be supplied by user:

Centrifuge capable of spinning 50 ml tubes at 2500 x g using swing-out rotor

Pipettes (1ml and 10 ml)

Vortex-Genie 2 Vortex

Vortex Adapter for 2 (50 ml) tubes (cat. no. 13000-V1-50)

Safety warnings

-  Solution C5 contains ethanol and is flammable.
DO NOT add bleach or acidic solutions directly to the sample preparation waste.



Sample preparation & cell lysis

1h 27m 30s

1 **ADD**  15 mL of PowerBead Solution to a PowerMax Bead Pro Tube

1m

ADD up to  10 g of soil sample to the PowerMax Bead Pro Tube containing PowerBead Solution

VORTEX vigorously for  00:01:00

Note

Refer to manufacturer's Troubleshooting Guide before deciding on the amount of soil to process. However, higher volumes of sediment (10 g) have shown better detection rates for fish sedDNA.

2 **ADD**  1.2 mL of Solution C1 to the PowerMax Bead Pro Tube

30s

VORTEX vigorously for  00:00:30

3 **PLACE** PowerMax Bead Pro Tube on a vortex adapter

40m

VORTEX for  00:10:00 at the highest speed

ALTERNATIVELY, place the tube in a shaking water bath set at  65 °C and shake at maximum speed for  00:30:00

4 **CENTRIFUGE** at  2500 x g for  00:03:00 at  Room temperature

3m

TRANSFER supernatant to a clean collection tube

Inhibitor removal

1h 27m 30s

5 **ADD**  5 mL of Solution C2

10m

INVERT twice to mix

INCUBATE at  2 °C to  8 °C for  00:10:00



6 **CENTRIFUGE** at 2500 x g for 00:04:00 at Room temperature 4m

AVOIDING the pellet, transfer the supernatant to a clean collection tube

7 **ADD** 4 mL of Solution C3 10m

INVERT twice to mix

INCUBATE at 2 °C to 8 °C for 00:10:00

8 **CENTRIFUGE** at 2500 x g for 00:04:00 at Room temperature 4m

AVOIDING the pellet, transfer the supernatant to a clean collection tube

Bind DNA

1h 27m 30s

9 **SHAKE** to mix Solution C4

ADD 30 mL of Solution C4 to supernatant

INVERT twice to mix

10 **FILL** an MB Maxi Spin Column with the solution from step 9 2m

CENTRIFUGE at 2500 x g for 00:02:00 at Room temperature

DISCARD the liquid flow-through

11 **ADD** a second volume of supernatant to the same MB Maxi Spin Column 2m

CENTRIFUGE at 2500 x g for 00:02:00 at Room temperature

DISCARD the liquid flow-through

REPEAT until the entire volume has been processed. This will take up to **4 total spins**

Wash spin column

1h 27m 30s

12 **ADD** 10 mL of Solution C5 3m



CENTRIFUGE at  2500 x g for  00:03:00 at  Room temperature

DISCARD the liquid flow-through

- 13 **CENTRIFUGE** again at  2500 x g for  00:05:00 at  Room temperature to 5m

CAREFULLY place the MB Maxi Spin Column in a new collection tube. Avoid splashing any residual Solution C5 onto the column

Elute the DNA

1h 27m 30s

- 14 **ADD**  5 mL of sterile Solution C6 to the center of MB Maxi Spin Column membrane 3m

CENTRIFUGE at  2500 x g for  00:03:00 at  Room temperature

DISCARD the MB Maxi Spin Column

DNA is now ready for downstream applications

Protocol references

<https://www.qiagen.com/us/products/discovery-and-translational-research/dna-rna-purification/dna-purification/microbial-dna/dneasy-powermax-soil-kit>

Successful fish sedDNA publications using this protocol:

Sales, N. G., Wangenstein, O. S., Carvalho, D. C., Deiner, K., Præbel, K., Coscia, I., McDevitt, A. D., & Mariani, S. (2021). Space-time dynamics in monitoring neotropical fish communities using eDNA metabarcoding. Science of The Total Environment, 754, 142096. <https://doi.org/10.1016/j.scitotenv.2020.142096>