



Apr 02, 2023

## qPCR of $\alpha$ -synuclein, TNF and NF- $\kappa$ B

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzb2z2vx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzb2z2vx1/v1)

Michael J Hurley<sup>1,2</sup>

<sup>1</sup>Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland, USA



Michael J Hurley

UCL

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzb2z2vx1/v1>

**Protocol Citation:** Michael J Hurley 2023. qPCR of  $\alpha$ -synuclein, TNF and NF- $\kappa$ B . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rm7vzb2z2vx1/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** February 27, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 77702

**Keywords:** ASAPCRN, synuclein, transcription pcr,  $\kappa\beta$

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000420

## Abstract

A procedure for quantitative real time reverse transcription PCR of  $\alpha$ -synuclein, TNF and  $NF-\kappa\beta$

## Materials

*$\beta$ -actin*

Forward 5'-GCTGTCCCTGTATGCCTCTG-3'

Reverse 5'-GATGTCACGCACGATTTCCC-3'

*$\alpha$ -synuclein*

Forward 5'-CAGAGGCAGCTGGAAAGACA-3'

Reverse 5'-CACCACTGCTCCTCCAACAT-3'

*TNF*

Forward 5'-GAGCCAGCGCGCCAACGCCCTCCT-3'

Reverse 5'-TGAGGAGCACGTAGTCGGGGCAGC-3'

*$NF-\kappa\beta$*

Forward 5'-ACTGCCGAGCTCAAGATCTGCCGA-3'

Reverse 5'-AAGGAGCCTCGTGCCTCCCAGCCT-3'

## Troubleshooting



## qPCR

- 1 Extract total RNA and protein from cells or tissue using TRI Reagent® solution according to the manufacturer's instructions. Use RNase free glycogen as a carrier for the total RNA.  
Keep protein for use in western blot.
- 2 Dissolve the RNA in nuclease free water.
- 3 Measure the concentration of RNA with a spectrophotometer (e.g. Nanodrop 1000).
- 4 Make cDNA from 2 micrograms of total RNA with a high-capacity cDNA reverse transcription kit (Invitrogen) according to manufacturer's instructions.
- 5 Dilute cDNA from the reverse transcription reaction 1:10 with nuclease free water.
- 6 Amplify 2 microliters of cDNA with gene specific intron-spanning primers (see materials) using PowerUp™ SYBR™ Green Master Mix (Invitrogen) according to the manufacturer's instructions with a StepOne Real-Time PCR system (Applied Biosciences).
- 7 Analyse using the using the comparative  $2^{-\Delta\Delta C_T}$  method.