



Feb 03, 2025

qPCR – Power SYBR Green Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rxw6g1b/v1>

Carolina Lopez¹

¹Washington University



Carolina Lopez

Washington University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rxw6g1b/v1>

Protocol Citation: Carolina Lopez 2025. qPCR – Power SYBR Green Protocol . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5rxw6g1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 29, 2024



Last Modified: February 03, 2025

Protocol Integer ID: 100777

Keywords: protocol for qpcr, green protocol, qpcr, protocol

Abstract

Protocol for qPCR using SYBR Green

Materials

REAGENT	SOURCE	IDENTIFIER
SYBR green	Thermo	Cat # 4367660
Working Primer Mix: 5μM (See directions below) Note: All primers should be cataloged in the Google Drive. Make sure to include at least 2 housekeeping genes.		
cDNA from 250ng-2ug of RNA starting material (500ng-1ug is preferred).		



qPCR – Power SYBR -

2m

1. If needed, prepare **5 μ M Working Primer Dilutions** for each gene of interest
 - a)  20 μ L  100 micromolar (μ M) Forward primer (final concentration: 2.5 μ M)
 - b)  20 μ L  100 micromolar (μ M) Reverse primer (final concentration: 2.5 μ M)
 - c)  760 μ L nuclease-free water

Note

d) In general,  100 micromolar (μ M) primers should be stored in the general lab primer boxes at  -20 $^{\circ}$ C . Once you make your own Working Primer Dilution, you may keep this in your personal  -20 $^{\circ}$ C box for further use.

e) IMPORTANT: Dilute both virus primers and housekeeping genes using this procedure, ALL PRIMERS SHOULD BE THE SAME CONCENTRATION. If the primer concentration is too high, there may be non-specific binding/amplification.

2. Prepare your plate map and PCR tubes. Quick tip:

Organize your cDNA samples in PCR tubes



so you can use the multichannel pipette properly. The pipette only fits every other row/column in the 384 well qPCR plate.

Therefore, you want to split your cDNA samples into groups, alternating the order between groups. For

the SYBR Green/Primer Master Mix, you will need to aliquot your mix into multiple PCR tubes. To determine the number of tubes, take # of samples in plate (rows)/2 and round up if needed.

a) Ex: for 12 samples, will need GAPDH Master Mix in 6 separate PCR tubes to use the multichannel pipette.



3. **Prepare cDNA** (will use  4 μL per well). Need a 1:40 dilution from an RT reaction with 250 ng-2 μg of starting material (500 ng-1 μg is preferred). Prepare the dilutions in  0.2 mL PCR tubes, mix each dilution (pipet or vortex) and spin down on a benchtop centrifuge.

a) Example dilutions to obtain 1:40 (*use nuclease free water):

- i.  1 μL cDNA in 39 μL H₂O
- ii.  2 μL cDNA in 78 μL H₂O
- iii.  3 μL cDNA in 117 μL H₂O
- iv.  4 μL cDNA in 156 μL H₂O

b) Amount of cDNA needed: # of genes x # of repetitions x 4 μL x 1.25

c) Do at least 3 technical replicates per sample (the x1.25 accounts for errors)

4. **Prepare SYBR Green + Primer Master Mix** (will use  6 μL per well). Prepare the dilutions in  0.2 mL PCR tubes or 1.5mL eppendorf tubes, mix each Master Mix well (pipet or vortex) and spin down quickly in the benchtop centrifuge. If needed, can split into multiple tubes for easy multichannel pipetting (see Step 2).



a) SYBR Green (will use  5 μL per well)

i. Amount needed = # of samples x # of repetitions x 5 μL x 1.25

b) Working Primer Dilution 5 μM (will use  1 μL per well)

i. Amount needed = # of samples x # of repetitions x 1 μL x 1.25

5. **Add**  4 μL **of cDNA per well** (use multichannel pipette with  12.5 μL filter pipet tips)

a) Can gently tap plate after to bring cDNA to the bottom of each well

b) Best to visually inspect the multichannel to ensure you are drawing up the correct volume in each tip.

6. **Add**  6 μL **of SYBR Green + Primer Master Mix** per well

a) Mix gently using the pipet, being careful to avoid air bubbles.

b) Best to visually inspect the multichannel to ensure you are drawing up the correct volume in each tip.

7. **Seal plate** with plate cover. It is important the plate is fully sealed or the reaction mix will evaporate during the run. Use the scraper edge to seal each edge of the plate cover VERY well.

8. **Spin plate** at 2000 rpm for  00:02:00 .

9. **Run** the Standard Lopez Lab Protocol – Power SYBR on qPCR machine Settings are set as what is listed below. Double check the settings are correct before each run!

a) **Lopez Lab Protocol – Power SYBR settings**

i. Experiment Properties:

- Block type: 384-Well Block
- Experiment type: Comparative C_T ($\Delta\Delta C_T$)
- Chemistry: SYBR® Green Reagents
- Run Mode: Standard

ii. Experiment Method:

- Volume:  10 μL
- Cover:  105 $^{\circ}\text{C}$

iii. Plate Attributes:

- Passive Reference: ROX
- Reference Sample: Sample 1
- Endogenous Control: Varies from experiment

STAGE	STEP	TEMPERATURE	TIME
Hold	Step 1	50° C	2 minutes
	Step 2	95° C	10 minutes
PCR Stage (runs for 40 cycles)	Step 1	95° C	15 seconds
	Step 2	55° C	25 seconds
	Step 3*	60° C	35 seconds
Melt Curve Stage	Step 1	95° C	15 seconds
	Step 2	60° C	1 minute
	Step 3 (Dissociation)*	95° C	15 seconds

*Image capture/data collection point

10. When your experiment is done, export your data from the machine and analyze.

2 Data Export and Analysis (Lopez-Lab specific)

At this point, it can be helpful to look at your data on the qPCR machine. In particular, the melt and amplification curves are useful. For the amplification curves, you should look to see that your technical replicates are as close to overlapping as possible. For the melting curves, you want to see a single peak for each gene of interest, and the peaks should be clustered together.. If your curves do not look correct, there was a problem during some step of the protocol (ie sample and plate preparation, primer design, machine, etc).

11. Save your raw data. This file is very important and needs to be saved in your RAW DATA folder.

12. To analyze your data, open the Lopez Lab qPCR Data Analysis template.

- Locate the template in the Google Drive qPCR Protocol Folder under the name: "DO NOT EDIT_qPCR analysis template_with fold induction_DO NOT EDIT"
- Do not edit the analysis template in the Google Drive
- Download your own copy of the template. Each time the template is used, the "Save As" function should be used to create a new excel file with your processed data from the experiment. The template should always remain the same.

13. Go to the RESULTS tab on the raw data excel file (the file saved from the qPCR machine). Copy the CT values from each of your samples into the Analysis Template. Make sure your samples and genes of interest are properly labeled and there are no errors in the copy/pasting.



14. If there are any samples that have an undefined CT, assign a value of 40 (because the PCR ran for 40 cycles).

15. At this point, the analysis template should auto-populate. The CT values for each of the housekeeping genes will be averaged to create a housekeeping index (HKI). Then the CT values for each gene of interest will be used to calculate Copy Numbers and Fold Induction.

Note

For the fold induction, the template is set up so that Sample 1 is the control (ie the sample entered in the first row). If this is not your control, change the excel calculations accordingly.

16. For more reading on qPCR data analysis, there are references in the Google Drive (REFERENCE DOCUMENTS → RT-qPCR)

Protocol references

Updated 2023.05.01 RBR