

Aug 09, 2023

qRT-PCR sample preparation

DOI

dx.doi.org/10.17504/protocols.io.4r3l22713l1y/v1

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis

WEHI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l22713l1y/v1>

Protocol Citation: Louise Uoselis 2023. qRT-PCR sample preparation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l22713l1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024







Protocol Integer ID: 86219

Keywords: ASAPCRN, pcr sample preparation, pcr sample preparation protocol, pcr sample preparation for analysis, rotorgeneq machine, qiagen



Abstract

Protocol for qRT-PCR sample preparation for analysis using a RotorGeneQ machine (Qiagen).

Troubleshooting

- 1 Synthesise cDNA libraries from total RNA from each sample using a High-Capacity cDNA Reverse Transcription Kit (Applied Biosystems), using Oligo(dT)₂₀ primers.
- 2 Place the synthesised libraries on ice, and dilute the libraries 1:3 by adding  60 µL of DEPC-treated H₂O to each sample, pipetting up and down gently to mix ~ 8 times.
- 3 Dilute primer stocks by adding  2 µL of the forward primer, and  2 µL of the reverse primer, to  46 µL of DEPC-treated H₂O on ice to make a final combined forward and reverse primer stock of  4 micromolar (µM) forward primer and  4 micromolar (µM) reverse primer for each gene target to be analysed.
- 4 Thaw the 2x QuantiNova SYBR green master mix (Qiagen) on ice.
- 5 Depending on the number of samples being analysed, assemble a sample master mix (not containing the cDNA library) containing the following:

Reagent	Volume (1x)
2x SYBR Green Master Mix	5 uL
4 uM forward and 4 uM reverse primer master mix	1 uL
cDNA library (to be added individually to each tube)	4 uL
Total volume	10 uL

- 6 To qRT-PCR tubes sitting on ice (Qiagen – RotorGeneQ compatible), add  6 µL of the master mix for each target gene to the desired tubes, and add  4 µL of the appropriate cDNA library to each sample.
- 7 Tap the rack holder to ensure all solution is at the bottom of the tube before capping each sample.(NOTE: samples do not need to be mixed prior to analysis)
- 8 Load and run the samples on a RotorGeneQ (Qiagen) machine using the desired parameters for the gene targets you are analysing.

