

Oct 16, 2025

qRT-PCR (SYBR Green)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjebxwgk5/v1>

Armin Bayati¹, Yuhang Zhao¹

¹Massachusetts General Hospital and Harvard Medical School



Jackson Schumacher

Massachusetts General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjebxwgk5/v1>

Protocol Citation: Armin Bayati, Yuhang Zhao 2025. qRT-PCR (SYBR Green). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvjebxwgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2025



Last Modified: October 16, 2025

Protocol Integer ID: 229915

Keywords: ASAPCRN, transcription pcr, using trizol rna extraction, trizol rna extraction, reverse transcription, relative gene expression, qrt, pcr, sybr, rna

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-237603

Abstract

Quantitative reverse-transcription PCR (qRT-PCR) to measure relative gene expression using TRIzol RNA extraction, PrimeScript™ RT Master Mix, and SYBR™ Green on a 96-well real-time instrument (e.g., StepOnePlus™). Analysis by $2^{-\Delta\Delta CT}$ with melt-curve confirmation.

Attachments



qRT-

PCR_protocol_Bay...

37KB

Guidelines

Pause Points:

- RNA stable at -80°C .
- cDNA stable at -20°C .
- Pre-run plate: keep on ice.

Troubleshooting:

- Low RNA yield: ensure clean separation, avoid overdrying.
- Multiple melt peaks: increase annealing temp, redesign primers.
- High variability: mix reagents well, centrifuge plate.
- NTC amplification: replace reagents, check contamination.

Record Keeping:

Record RNA purity, input amounts, primer sequences, plate map, run file, and $\Delta\Delta CT$ calculations.



Materials

Instruments: NanoDrop, PCR thermal cycler, real-time PCR system (96-well), plate centrifuge.

Consumables: Pipettes/tips, RNase-free tubes, 96-well qPCR plates, optical seals.

Reagents: TRIzolTM Reagent; chloroform; isopropanol; 70% ethanol (RNase-free); RNase-free water; PrimeScriptTM RT Master Mix (TAKARA RR036A); SYBRTM Green Universal Master Mix (Applied Biosystems 4309155); forward/reverse primers.

Safety warnings

- ! - Work RNase-free (gloves, RNase-free plastics).
- TRIzolTM, chloroform, and isopropanol are hazardous—use in a fume hood.
- Keep samples on ice unless otherwise specified.

Before start

- Pre-cool a centrifuge to 4 °C.
- Thaw reagents on ice; mix and briefly spin.
- Design primers (amplicon 70–200 bp; T_m ~60 °C when possible).



Step 1 – Total RNA Extraction (TRIzol)

- 1 Lyse/Homogenize: Add 1 mL TRIzol per 50–100 mg tissue or $5\text{--}10 \times 10^6$ cells. Homogenize thoroughly; incubate 5 min at RT.
- 2 Phase Separation: Add 0.2 mL chloroform per 1 mL TRIzol, vortex 10 s, incubate 10 min at RT. Centrifuge $12,000 \times g$, 10 min, 4 °C. Transfer aqueous phase.
- 3 RNA Precipitation: Add 0.5 mL isopropanol per 1 mL TRIzol, invert gently, incubate 10 min RT. Centrifuge $12,000 \times g$, 10 min, 4 °C.
- 4 Wash: Wash pellet 3× with 70% ethanol, centrifuge each time. Air-dry 30 min RT.
- 5 Resuspend & QC: Dissolve in 50 µL RNase-free water, measure purity ($A_{260}/A_{280} > 1.8$). Store at –80 °C.

Step 2 – Reverse Transcription (cDNA Synthesis)

- 6 Reaction (20 µL): 1 µg RNA, 4 µL 5× PrimeScript RT Mix, water to 20 µL. Program: 37 °C 15 min, 85 °C 5 s, 4 °C hold. Dilute to 100 µL, store at –20 °C.

Step 3 – qPCR Reaction Setup

- 7 Reaction mix (10 µL): 5 µL SYBR Mix, 1 µL cDNA, 0.5 µL primers each, 3 µL water. Set up triplicates, NTCs. Seal, centrifuge plate.

Step 4 — Cycling Conditions & Melt Curve

- 8 Option A ($T_m \geq 60$ °C): 50 °C 2 min; 95 °C 2 min; 40× (95 °C 15 s, 60 °C 1 min).
Option B ($T_m < 60$ °C): 50 °C 2 min; 95 °C 2 min; 40× (95 °C 15 s, 55–60 °C 15 s, 72 °C 1 min).
Melt Curve: 95 °C 15 s; 60 °C 1 min; slow ramp to 95 °C.

Step 5 — Data Analysis ($2^{-\Delta\Delta CT}$)

- 9 Export CT values.



10 $\Delta CT = CT(\text{target}) - CT(\text{housekeeping}).$

11 $\Delta\Delta CT = \Delta CT(\text{experimental}) - \Delta CT(\text{control}).$

12 Expression = $2^{-\Delta\Delta CT}.$

13 5. Confirm single melt peak. Plot in Prism.