

Jul 05, 2018

Quantification of circulating microRNA using droplet digital PCR

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.q8edzte

Helle Glud Binderup¹, Jonna Skov Madsen¹, Kim Houliind², Rikke Fredslund Andersen³, Claus Lohman Brasen¹

¹Biochemistry and Immunology, Lillebaelt Hospital, Kolding and Vejle, Denmark;

²Department of Vascular Surgery, Lillebaelt Hospital, Kolding, Denmark;

³Biochemistry and Immunology, Lillebaelt Hospital, Kolding and Vejle, Denmark



Helle Glud Binderup

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.q8edzte

External link: <https://doi.org/10.1371/journal.pone.0201069>

Protocol Citation: Helle Glud Binderup, Jonna Skov Madsen, Kim Houliind, Rikke Fredslund Andersen, Claus Lohman Brasen 2018. Quantification of circulating microRNA using droplet digital PCR. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q8edzte>

Manuscript citation:

Binderup HG, Madsen JS, Heegaard NHH, Houliind K, Andersen RF, Brasen CL (2018) Quantification of microRNA levels in plasma – Impact of preanalytical and analytical conditions. PLoS ONE 13(7): e0201069. doi: [10.1371/journal.pone.0201069](https://doi.org/10.1371/journal.pone.0201069)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 22, 2018

Last Modified: July 05, 2018

Protocol Integer ID: 13286

Keywords: Circulating microRNAs, sample preparation, microRNA purification, reverse transcription, ddPCR, TaqMan assays

Abstract

A protocol for quantification of circulating microRNA using droplet digital PCR with TaqMan assays.

MicroRNAs are purified from platelet-poor plasma using Nucleospin columns

During the purification step samples are spiked with cel-miR-39 as a mean of technical normalization

Normalization are performed by calculating the relative concentration of the target microRNA and the reference.

Guidelines

Blood samples should be obtained using a minimum of venous stasis and with discard of the first 3 mL of blood.

Platelet-poor plasma should be prepared within 2 hours from blood sampling



Materials

MATERIALS

⊗ K2-EDTA containing tubes **Becton Dickinson (BD) Catalog #366643**

⊗ Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**

⊗ Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**

⊗ DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**

STEP MATERIALS

⊗ Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**

⊗ DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**

⊗ Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**

⊗ Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**




























⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**


⊗ DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**

⊗ Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**

Protocol materials

-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  K2-EDTA containing tubes **Becton Dickinson (BD) Catalog #366643**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**
-  Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**
-  ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
-  Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**
-  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**
-  DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**
-  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**
-  TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**
-  ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
-  TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**
-  ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
-  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**
-  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**
-  TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**
-  ddPCR Supermix for probes (no dUTP) **Bio-Rad Laboratories Catalog #1863024**
-  DG32™ Automated Droplet Generator Cartridges **Bio-Rad Laboratories Catalog #1864108**



 Automated Droplet Generation Oil for Probes **Bio-Rad Laboratories Catalog #1864110**

Before start

We recommend to use a few samples to test if the volume of spike-in added during microRNA purification is adequate before purifying all your samples.

The concentration of the spike-in and the target miRNA should be within the same range, otherwise adjust the volume of spike-in added during miRNA purification.



Preparation of platelet-poor plasma (PPP)

- 1 Dual centrifugation
 1. use 10 ml of EDTA anticoagulated whole blood
 2. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
 3. transfer plasma phase to new tube, leaving approximately 1 mL of plasma on top of the buffy coat
 4. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
 5. transfer plasma phase to cryo tubes, leaving approximately 1 ml of plasma in the bottom of the tube
 6. store at -80 °C

Equipment

new equipment

NAME

Hettich centrifuge

BRAND

4706-01

SKU

Rotina 420R

SPECIFICATIONS

Note















Alternatively, a prolonged single centrifugation may be used

1. use 5 ml of EDTA anticoagulated whole blood
2. centrifuge at 3000 g for 30 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to cryo-tubes, leaving approximately 0.5 ml of plasma on top of the buffy coat

MicroRNA purification

- 2  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

- 3 Follow the instructions given by the manufacture (se notes before you start):

NucleoSpin® miRNA Plasma		
1 Prepare sample		300 µL plasma or serum* 90 µL MLP Vortex 5 s RT, 3 min
2 Precipitate protein	 	30 µL MPP Vortex 5 s RT, 1 min 11,000 x g, 3 min
3 Transfer supernatant		Transfer clear supernatant to Collection Tube (2 mL, lid)
4 Adjust binding conditions	 	400 µL isopropanol Vortex 5 s
5 Bind RNA and DNA	 	Load sample on NucleoSpin® miRNA Column RT, 2 min 11,000 x g, 30 s
6 Optional: Digest DNA	 	Optional: 1st 700 µL MW2 11,000 x g, 30 s 2nd 250 µL MW2 11,000 x g, 2 min 50 µL rDNase in Reaction Buffer for rDNase RT, 15 min
7 Wash and dry	 	1st 100 µL MW1 11,000 x g, 30 s 2nd 700 µL MW2 11,000 x g, 30 s 3rd 250 µL MW2 11,000 x g, 2 min
8 Elute RNA	 	30 µL RNase-free H ₂ O RT, 1 min 11,000 x g, 1 min

Note

STEP 1: use 300 µL platelet-poor plasma

STEP 4: after addition of isopropanol, add also 5 µL spike-in (cel-miR-39, 2.75×10^{-12} M)

STEP 6: perform the optional DNA digest

Reverse transcription

4



TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**



20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

5

Prepare RT master mix:

Component	Master mix volume per 15 µL reaction*
100 mM dNTPs (with dTTP)	0.15 µL
Multiscribe RT enzyme (50 U/µL)	1 µL
10x RT buffer	1.5 µL
RNase inhibitor (20 U/µL)	0.19 µL
Nuclease free water	to 15 µL in total**

*add 10-20% excess volume

**add water to a total reaction volume of 15 µL (including microRNA and RT-primers, see step 7)

Mix gently and place on ice

Note

Remember to include a RT-negative sample (no template)

It is also a good idea to include a RT-positive sample (a microRNA-sample included in all runs)

6

For each RT reaction, combine RT master mix with 2 µL of purified microRNA in a 0.2 µL polypropylene reaction tube

(the volume of RT master mix is dependent on the number of RT-primers, see step 7)

Keep on ice

7

Prepare RT primer mix:

Add for each RT reaction 0.75 µL 20x RT primer from each microRNA assay set

Example with two microRNA assays:

Component	Volume per 15 μ L reaction*
20x RT specific primer #1	0.75 μ L
20x RT specific primer #2	0.75 μ L
Total volume	1.5 μ L**



*add 10-20% excess volume

**when using two microRNA assays the volume of water needed in step 5 is 8.66 μ L for each reaction, and the volume of master mix added in step 6 is 11.5 μ L

Mix gently

- 8 Add the appropriate volume of RT primer mix to the reaction tubes prepared in step 6
Mix gently and incubate on ice for 5 minutes
- 9 Transfer the reaction tubes to a thermocycler
Incubate at 16 °C for 30 min, at 42 °C for 30 min and at 85 °C for 5 min.
- 10 Continue immediately to the PCR amplification or store the RT-reaction in -20 °C.

Droplet digital PCR

- 11  20X TaqMan MicroRNA Assay Thermo Fisher Scientific Catalog #4440887
 ddPCR Supermix for probes (no dUTP) Bio-Rad Laboratories Catalog #1863024

- 12 Prepare ddPCR reaction mix:

Component	Volume per 20.3 μ L reaction*
2x ddPCR supermix	24 μ L
20x TaqMan microRNA assay	2.4 μ L
Nuclease free water	19 μ L
Total volume	45.4 μL


*add 10-20% excess volume


Mix gently

- 13 1. Dilute each of the RT-reactions (cDNAs) 1:10 with nuclease free water (e.g. 2 µL cDNA + 18 µL water)
2. Combine 45.4 µL of ddPCR reaction mix with 2.6 µL of diluted RT-reaction (cDNA) in a 96 well plate

Droplet generation

14

 DG32™ Automated Droplet Generator Cartridges **Bio-Rad**
Laboratories Catalog #1864108

 Automated Droplet Generation Oil for Probes **Bio-Rad**
Laboratories Catalog #1864110

Equipment

new equipment

NAME

AutoDG

BRAND

1864101

SKU

Automated Droplet Generator from BioRad

SPECIFICATIONS

15

1. Transfer 22 µL of the reaction mixture into each of two wells on an empty 96 well plate
2. Heat sealing
3. Transfer plate to the AutoDG
4. The AutoDG transfers 20 µL from each well to a GD32 cartidge for droplet generation
5. Heat sealing of output plate

PCR

- 16 Transfer plate to termocycler
Incubate in 95 °C for 5 minutes
Proceed with 44 cycles of 95 °C for 15 sec and 60 °C for 60 sec.



Final incubation in 98 °C for 10 min

Cool to 4 °C

Droplet reading

17

Software

QuantaSoft Software

NAME

<http://www.bio-rad.com/en-dk/sku/1864011-quantasoft-software-regulatory-edition?ID=1864011>

SOURCE
LINK

18 ddPCR analysis is performed with QX100 Droplet Reader and QuantaSoft Software

Normalization

19 Normalization is performed by calculating the relative concentration of target miRNA and cel-miR-39