

Jul 05, 2018

# Quantification of circulating microRNA using TaqMan Low Density Array (TLDA)

 [PLOS One](#)

DOI

[dx.doi.org/10.17504/protocols.io.q62dzge](https://dx.doi.org/10.17504/protocols.io.q62dzge)

Helle Glud Binderup<sup>1</sup>, Jonna Skov Madsen<sup>1</sup>, Kim Houliind<sup>2</sup>, Rikke Fredslund Andersen<sup>1</sup>, Claus Lohman Brasen<sup>1</sup>

<sup>1</sup>Biochemistry and Immunology, Lillebaelt Hospital, Kolding and Vejle, Denmark;

<sup>2</sup>Department of Vascular Surgery, Lillebaelt Hospital, Kolding, Denmark



Helle Glud Binderup

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.q62dzge](https://dx.doi.org/10.17504/protocols.io.q62dzge)

External link: <https://doi.org/10.1371/journal.pone.0201069>

**Protocol Citation:** Helle Glud Binderup, Jonna Skov Madsen, Kim Houliind, Rikke Fredslund Andersen, Claus Lohman Brasen 2018. Quantification of circulating microRNA using TaqMan Low Density Array (TLDA). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q62dzge>

## Manuscript citation:

Binderup HG, Madsen JS, Heegaard NHH, Houliind K, Andersen RF, Brasen CL (2018) Quantification of microRNA levels in plasma – Impact of preanalytical and analytical conditions. PLoS ONE 13(7): e0201069. doi: [10.1371/journal.pone.0201069](https://doi.org/10.1371/journal.pone.0201069)

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 21, 2018

**Last Modified:** July 05, 2018

**Protocol Integer ID:** 13242



**Keywords:** Circulating microRNAs, sample preparation, microRNA purification, reverse transcription, TLDA, TaqMan assays

## Abstract

A protocol for quantification of circulating microRNA using TaqMan Low Density Array.

MicroRNAs are purified from platelet-poor plasma using Nucleospin columns

During the purification step samples are spiked with cel-miR-39 as a mean of technical normalization

Normalization is performed using the  $\Delta$ Ct-method

## Guidelines

Blood samples should be obtained using a minimum of venous stasis and with discard of the first 3 mL of blood.

Platelet-poor plasma should be prepared within 2 hours from blood sampling



## Materials

### MATERIALS

⊗ K2-EDTA containing tubes **Becton Dickinson (BD) Catalog #366643**

⊗ Nucleospin®miRNA Plasma **Macherey-Nagel Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

### STEP MATERIALS

⊗ Nucleospin®miRNA Plasma **Macherey-Nagel Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ Nucleospin®miRNA Plasma **Macherey-Nagel Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**


⊗ 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**



Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**


## Protocol materials


 Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

 K2-EDTA containing tubes **Becton Dickinson (BD) Catalog #366643**


 Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**


 TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**


 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

 TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**


 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**


 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**


 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**


 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**


 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**


 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**


 TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

 Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

 TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**

 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**

 TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

 Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**



⊗ Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

⊗ TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

⊗ TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**

⊗ 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

⊗ 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**

⊗ Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

## Before start

We recommend to test if the volume of spike-in added during microRNA purification is adequate before purifying all your samples.

When testing the preamplification product, the Ct-values of the spike-in and the target miRNA should be within the same range, otherwise adjust the volume of spike-in added during miRNA purification.



## Preparation of platelet-poor plasma (PPP)

- 1 Dual centrifugation
  1. use 10 ml of EDTA anticoagulated whole blood
  2. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
  3. transfer plasma phase to new tube, leaving approximately 1 mL of plasma on top of the buffy coat
  4. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
  5. transfer plasma phase to cryo tubes, leaving approximately 1 ml of plasma in the bottom of the tube
  6. store at -80 °C

### Equipment

new equipment

NAME

Hettich centrifuge

BRAND

4706-01

SKU

Rotina 420R

SPECIFICATIONS

### Note















Alternatively, a prolonged single centrifugation may be used

1. use 5 ml of EDTA anticoagulated whole blood
2. centrifuge at 3000 g for 30 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to cryo-tubes, leaving approximately 0.5 ml of plasma on top of the buffy coat

## MicroRNA purification

- 2  Nucleospin®miRNA Plasma **Macherey-Nagal Catalog #740971.50**

- 3 Follow the instructions given by the manufacture (se notes before you start):

<b>NucleoSpin® miRNA Plasma</b>		
<b>1 Prepare sample</b>		300 µL plasma or serum* 90 µL MLP Vortex 5 s RT, 3 min
<b>2 Precipitate protein</b>	 	30 µL MPP Vortex 5 s RT, 1 min 11,000 x g, 3 min
<b>3 Transfer supernatant</b>		Transfer clear supernatant to Collection Tube (2 mL, lid)
<b>4 Adjust binding conditions</b>	 	400 µL isopropanol Vortex 5 s
<b>5 Bind RNA and DNA</b>	 	Load sample on NucleoSpin® miRNA Column RT, 2 min 11,000 x g, 30 s
<b>6 Optional: Digest DNA</b>	 	<b>Optional:</b> <b>1<sup>st</sup></b> 700 µL MW2 11,000 x g, 30 s <b>2<sup>nd</sup></b> 250 µL MW2 11,000 x g, 2 min 50 µL rDNase in Reaction Buffer for rDNase RT, 15 min
<b>7 Wash and dry</b>	 	<b>1<sup>st</sup></b> 100 µL MW1 11,000 x g, 30 s <b>2<sup>nd</sup></b> 700 µL MW2 11,000 x g, 30 s <b>3<sup>rd</sup></b> 250 µL MW2 11,000 x g, 2 min
<b>8 Elute RNA</b>	 	30 µL RNase-free H <sub>2</sub> O RT, 1 min 11,000 x g, 1 min



**Note**

STEP 1: use 300  $\mu\text{L}$  platelet-poor plasma

STEP 4: after addition of isopropanol, add also 5  $\mu\text{L}$  spike-in (cel-miR-39,  $2.75 \times 10^{-12}$  M)

STEP 6: perform the optional DNA digest

**Reverse transcription**

4



TaqMan®MicroRNA Reverse Transcription Kit **Applied Biosystems (ThermoFisher Scientific) Catalog #4366597**



Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

5 Prepare RT master mix:

Component	Master mix volume per 10 $\mu\text{L}$ reaction*
Customs RT primer pool (10X)	1 $\mu\text{L}$
100 mM dNTPs (with dTTP)	0.27 $\mu\text{L}$
Multiscribe RT enzyme (50 U/ $\mu\text{L}$ )	2 $\mu\text{L}$
10x RT buffer	1 $\mu\text{L}$
MgCl <sub>2</sub> (25 mM)	1 $\mu\text{L}$
RNase inhibitor (20 U/ $\mu\text{L}$ )	0.1 $\mu\text{L}$
Nuclease free water	1.63 $\mu\text{L}$

\*add 10-20% excess volume

Mix gently and place on ice

**Note**

Remember to include a RT-negative sample (no template)

It is also a good idea to include a RT-positive sample (a microRNA-sample included in all runs)


6 For each RT reaction, combine 7  $\mu\text{L}$  RT master mix with 3  $\mu\text{L}$  of purified microRNA in a 0.2  $\mu\text{L}$  polypropylene reaction tube


Keep on ice



- 7 Mix gently and incubate on ice for 5 minutes
- 8 Transfer the reaction tubes to a thermocycler  
Incubate in 40 cycles of 16 °C for 2 min, at 42 °C for 1 min and 50 °C for 1 sec.  
Finish with 85 °C for 5 min and cool to 4 °C.
- 9 Continue immediately to the PCR amplification or store the RT-reaction in -20 °C.

## Preamplification of cDNA

- 10  Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

 TaqMan™ PreAmp Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4391128**

- 11 Prepare preamplification reaction mix:

Component	Volume per 25 µL reaction*
2x TaqMan PreAmp master mix	12.5 µL
Customs Primers (10X)	2.5 µL
Nuclease free water	5 µL
<b>Total volume</b>	<b>20 µL</b>

\*add 10-20% excess volume

Mix gently


- 12 Combine 20 µL of qPCR reaction mix with 5 µL of RT-reaction (cDNA) in a 96 well plate  
Seal plate
- 13 Incubate in 95 °C for 10 min, 55 °C for 2 min and 72 °C for 2 min  
Proceed with 14 cycles of 95 °C for 15 sec and 60 °C for 4 min.  
Finish with 99.9 °C for 10 min and cool to 4 °C

Continue immediately or store at 4 °C for up to 12 hours or at -20 °C for up to 1 week.

## Test preamplification product

- 14 Use a specific TaqMan assay for the spike-in (Cel-miR-39, assay 000200) and one of the target miRNAs included in the Array to test the preamplification product before loading it on the array.

 20X TaqMan MicroRNA Assay **Thermo Fisher Scientific Catalog #4440887**

 2X TaqMan Universal PCR Master Mix **Applied Biosystems (ThermoFisher Scientific) Catalog #4318157**

- 15 Prepare qPCR reaction mix (run each sample in doublets):

Component	Volume per 20.3 $\mu$ L reaction*
2x TaqMan Universal PCR master mix	10 $\mu$ L
20x TaqMan microRNA assay	1 $\mu$ L
Nuclease free water	8 $\mu$ L
<b>Total volume</b>	<b>19 <math>\mu</math>L</b>

\*add 10-20% excess volume

Mix gently

- 16 Combine 19  $\mu$ L of qPCR reaction mix with 1.3  $\mu$ L of RT-reaction (cDNA) in a 96 well plate  
Seal plate



17

**Software****Sequence Detection System (SDS)**

<https://www.thermofisher.com/search/results?query=sequence+analysis+software&persona=Catalog&navId=10949&refinementAction=>

18

Transfer plate to ABI Prism 7900HT  
Use the SDS software to set up the run  
Incubate in 50 °C for 2 min and 95 °C for 10 min  
Proceed with 40 cycles of 95 °C for 15 sec and 60 °C for 60 sec.

19

Evaluate the Ct-values obtained

**Loading of the array**

20



Custom TaqMan®Array MicroRNA Cards **Applied Biosystems (ThermoFisher Scientific) Catalog #4449135**

21

Prepare samples:

Component	Volume per 120 µL reaction
2x TaqMan Universal PCR II mix No AmpErase UNG	60 µL
PreAmplification product	1,2 µL
Nuclease free water	58,8 µL



<b>Total volume</b>	<b>120 <math>\mu</math>L</b>
---------------------	------------------------------

- 22 Load 100  $\mu$ L sample-dilution to each port on the array  
Centrifuge The array at 1200 rpm for 2 min  
Seal the array

## Real-time PCR

23

### Software

#### Sequence Detection System (SDS)

<https://www.thermofisher.com/search/results?query=sequence+analysis+software&persona=Catalog&navId=10949&refinementAction=>

- 24 Transfer array to ABI Prism 7900HT  
Use the SDS software to set up the run  
Incubate in 50 °C for 2 min and 94.5 °C for 10 min  
Proceed with 40 cycles of 97 °C for 30 sec and 59.7 °C for 60 sec.

## Normalization

- 25 Normalization is performed using the  $\Delta$ Ct-method ( $2^{-\Delta\text{Ct}}$ )  
 $\Delta\text{Ct} = \text{mean Ct}_{\text{target miRNA}} - \text{mean Ct}_{\text{cel-miR-39}}$  (mean of triplets)