

Nov 17, 2025

Quantification of human mDA/progenitor cells in grafted animal brain tissue



Forked from [Stereological estimation of human TH+ cells in grafted animal brain tissue](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnrwopl5d/v1>

Roberto Garcia Swinburn¹, Ernest Arenas¹

¹Karolinska Institute Stockholm

SOX6 mDA differentiation



Roberto Garcia Swinburn

Karolinska Institute Stockholm

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnrwopl5d/v1>

Protocol Citation: Roberto Garcia Swinburn, Ernest Arenas 2025. Quantification of human mDA/progenitor cells in grafted animal brain tissue. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.14egnrwopl5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 16, 2025

Last Modified: November 17, 2025

Protocol Integer ID: 227399

Keywords: grafted animal brain tissue, th positive cells within grafted animal, animal brain tissue, grafted animal, progenitor cell, human nuclei, positive cell

Abstract

We used this Qupath-based quantification for the percentage of human nuclei and TH positive cells within grafted animals with several other markers of interest.

Protocol materials

✕ OCT (Optimal Cutting Temperature compound) **Sakura Finetek Catalog #4583**

✕ Sheep Tyrosine Hydroxylase **PeI-Freez Catalog #P60101-0**

✕ Anti-Nuclei Antibody clone 235-1 **Merck Millipore (EMD Millipore) Catalog #MAB1281**

✕ DAPI **Merck MilliporeSigma (Sigma-Aldrich) Catalog #10236276001**

Imaging

- 1 Slice the fixed frozen brain in 5 series at 20 µm using a Cryostat at  -20 °C . The brain can be embedded in  OCT (Optimal Cutting Temperature compound) **Sakura Finetek Catalog #4583** to stabilize better. Each slice could be either in cold PBS1X or a cryoprotective solution for storage.

We ended with 15-25 slices of striatum with grafts.
- 2 Stain slices against TH, human nuclei and another marker of interest for mDA. We used  Sheep Tyrosine Hydroxylase **Pel-Freez Catalog #P60101-0** ,  Anti-Nuclei Antibody clone 235-1 **Merck Millipore (EMD Millipore) Catalog #MAB1281** and  DAPI **Merck MilliporeSigma (Sigma-Aldrich) Catalog #10236276001** for the TH and human nuclei. The remaining channel was used for a marker of interest
- 3 Mount slices on glass slides.
- 4 Image human cells at 20X magnification. TH cells, human nuclei and DAPI-stained nuclei should be discernible.
We used Zeiss LSM880 confocal microscopes



Analysis

- 5 Creat a New Project on Qupath.

Software

QuPath

NAME

Pete Bankhead

DEVELOPER

<http://qupath.github.io>

SOURCE LINK

- 6 Load images onto QuPath project
- 7 Create a Full Image annotation.
- 8 Use the Cell Detection tool (Analyze>Cell Detection>Cell Detection) to detect Human Nuclei.
Adjust the detection channel to the human nuclei and threshold to optimize detection.
- 9 Use the Classify tool (Classify>Object Classification>Create composite classifier). Adjust channel filter to TH channel, and threshold for optimized classification.
Above threshold should be TH, below should display nothing.
We recommend activating Live Preview to optimize Threshold.

Save the classifier.
- 10 Repeat step 9 with another channel with associated marker
- 11 Create a composite classifier (Classify>Object Classification>Create composite classifier) with your relevant markers.
- 12 Once you have saved the Cell Detection parameters and the Classification, re-use on all images.
- 13 Export the Image measurements in CSV
- 14 Calculate the percentage of double positive cells