



May 14, 2024

# Quantification of IHC Using an Inverted Confocal Microscope and NIS-Elements Program-Killinger Lab 2024

DOI

[dx.doi.org/10.17504/protocols.io.5qpvokrbzl4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvokrbzl4o/v1)

Solji Choi<sup>1</sup>

<sup>1</sup>Rush University Medical Center



Bryan Killinger

Rush University

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvokrbzl4o/v1>

**Protocol Citation:** Solji Choi 2024. Quantification of IHC Using an Inverted Confocal Microscope and NIS-Elements Program-Killinger Lab 2024. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvokrbzl4o/v1>

### Manuscript citation:

Choi SG, Tittle T, Garcia-Prada D, Kordower JH, Melki R, Killinger BA. Alpha-synuclein aggregates are phosphatase resistant. bioRxiv [Preprint]. 2024 Apr 9:2023.11.20.567854. doi: 10.1101/2023.11.20.567854. PMID: 38645137; PMCID: PMC11030248.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 14, 2024

**Last Modified:** May 14, 2024

**Protocol Integer ID:** 99764

**Keywords:** quantification of ihc, ihc, using ni, stained section, killinger lab, inverted confocal microscope, confocal microscope, dab, killinger lab 2024 protocol

**Funders Acknowledgements:**

**Michael J. Fox Foundation**

Grant ID: ASAP-021030

**Michael J. Fox Foundation**

Grant ID: ASAP-024442

## Abstract

Protocol to quantify DAB stained sections using NIS-Elements.

## Troubleshooting

- 1 Capture brightfield images with an inverted confocal microscope with a 20X objective (Nikon A1R).
- 2 Conduct annotation of each tissue section within a bounding box of 2000×2000 pixels for mouse tissues and 2863×2454 pixels for human tissues.
- 3 Use a manual RGB-based color thresholding algorithm for mouse tissues to exclude pure methyl green nuclei staining from the measurement. For human tissues, use NIS-elements (version 5.10.01, <https://www.microscope.healthcare.nikon.com/products/software/nis-elements>, RRID:SCR\_014329) auto-thresholding algorithm.
- 4 Export the percentage area of the thresholded signal (object area fraction provided in the program) and normalize to the average value of non-CIAP treated tissues.