

Sep 02, 2022

QuantiGene multiplex assay

DOI

dx.doi.org/10.17504/protocols.io.kqdg39ew7g25/v1

Amy R Hicks¹

¹UCL



Amy R Hicks

UCL

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg39ew7g25/v1>

Protocol Citation: Amy R Hicks 2022. QuantiGene multiplex assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg39ew7g25/v1>

Manuscript citation:

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6834638/>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 02, 2022

Last Modified: June 01, 2024

Protocol Integer ID: 69492

Keywords: Gene expression, Mid throughput, ASAPCRN, captured rna signal, quantigene multiplex, rna signal, rna, gene, set of gene, series of incubation, lysed sample, luminex, lab of professor gill bate, incubation, assay

Abstract

This protocol outlines a method of simultaneously measuring the expression of up to 80 genes in any one sample using the QuantiGene multiplex assay from ThermoFisher. This protocol was adapted from the lab of Professor Gill Bates. Cells are grown in 96-well plates, lysed, snap frozen and stored at -80°C . Lysed samples are then incubated overnight with magnetic capture beads as well as a probe panel to detect a custom designed set of genes. A series of incubations and washes amplify the captured RNA signals, which are subsequently measured using a fluorescent probe on a Magpix (Luminex).



Materials

Normal cell culture materials

- Media

QuantiGene reagents

- Lysis mixture
- Proteinase K
- Blocking reagent
- Pre-amplifier solution
- Amplifier solution
- Label probe solution
- SAPE
- SAPE diluent
- Wash buffer component 1
- Wash buffer component 2
- SAPE wash buffer
- Magnetic separation plate
- Plate seals
- Hybridization plate
- Pressure seals
- Probe set
- Capture beads







Lab equipment

- 200ul multichannel pipette
- Reagent reservoir
- Vortemp shaking incubator
- Magpix plate reader
- Handheld magnetic plate washer
- RNase-free water
- Dry ice

Troubleshooting





Sample preparation

- 1 Pre-warm lysis mixture at **37°C** for **30 mins**, followed by gentle swirling.  00:30:00 30m
- 2 Prepare **working lysis mixture** by adding  10 μL of **proteinase K** to each  1 mL of **lysis mixture** required.
- 3 Add **1/2 volume** of working lysis mixture to cells in **culture media** from a reagent reservoir using the multichannel pipette (e.g. for a 96 well plate containing  100 μL of media per well, add  50 μL of working lysis mixture).
- 4 **Mix** by pipetting up and down 3-4 times, discard the tips before continuing with the consecutive wells (more is more).
- 5 **Snap freeze** on a bed of **dry ice** and store in the **-80°C** until required.
- 6 When required, incubate the cell culture plate in the **Vortemp** pre-warmed to **50°C** for  01:00:00 without shaking. 1h
- 7 **Verify** cell lysis using the cell culture microscope.
- 7.1 For new QuantiGene plexes, appropriate dilutions must be ascertained via a serial dilution experiment- refer to [Papadopoulos et al, 2019](#) and pages **29–30** of the **QuantiGene Plex Gene Expression Assay User Guide** for more details.


Storage of lysates

- 8 Samples should be stored long term in the **-80°C** freezer. Samples do not need to be thawed on ice and are **stable** at **RT**.

Assay day 1

- 9  0 μL Pre-warm **lysis mixture** at **37°C** for  00:30:00 followed by gentle swirling. 30m



- 10 Arrange sufficient **pre-vortexed** and **diluted** samples (for **experiment**) or **serial dilutions** and **reference RNA** (for plex **optimisation**) as per your plate plan and keep at **RT**- each  40 μL sample should be run in **duplicate**, include a **background control** by making sufficient diluted lysis mixture (1 volume lysis mixture plus 2 volumes of RNase-free water).
- 11 Handle the reagents listed below as follows:
- 11.1 **Probe Set & Blocking Reagent** (kept at **-20°C** in QuantiGene reagents box): **thaw** and **vortex** briefly to mix, then **centrifuge probe set** briefly to collect contents at the bottom of the tube.
- 11.2 **Proteinase K** (kept at **-20°C** in QuantiGene reagents box): keep on **ice**.
- 11.3 **Capture Beads** (kept at **4°C** in QuantiGene capture bead box): take out of storage right before use and **protect** from **light**.
- 12 Prepare an appropriate volume of **working bead mix** by combining the following reagents in the order listed (**this is for 2- to 64-plex assays**), scale according to the number of wells on your QuantiGene plate(s), keep working bead mix at RT and protected from light.




12.1

	A	B	C	D
	Order	Reagent	1 well (μL)	96 wells (+14 for extra) (μL)
	1	Nuclease-free water	2.6	286
	2	Lysis Mixture	3.3	363
	3	Blocking Reagent	1	110
	4	Proteinase K	0.1	11
	5	Capture Beads (vortex for 30 seconds before adding)	0.5	55
	6	Probe Set	2.5	275



	A	B	C	D
	TOTAL:	10	1100	






Working bead mix

- 13 Vortex **Working Bead Mix** for **10 seconds** and then carefully pipette **10µl** into each well of a **magnetic separation plate**, avoiding bubbles, add  **40 µL** of each sample (including background controls) as per your plate plan into the magnetic separation plate (load each sample with a new pipette tip).
- 14 **Seal** magnetic separation plate with a pressure seal. Use the backing of the pressure seal to firmly and evenly apply pressure across the whole seal and lastly run your finger along each edge of the plate to seal.
- 15 Place the magnetic separation plate in the Vortemp shaking incubator for  **18:00:00** to  **22:00:00** at **54°C** at **600rpm**.

1d 16h













Assay day 2

5h 17m 15s


- 16 Turn on **Magpix** and computer to allow lasers time to **warm up**.
- 17 Warm **pre-amplifier solution**, **amplifier solution**, **label probe solution** and **SAPE diluent** at **37°C** at least  **00:30:00** prior to use.
- 18 Prepare **1 x wash buffer** by adding  **3 mL** **wash buffer component 1** and  **50 mL** **wash buffer component 2** and topping up to  **1 L** with **nuclease-free water** from the Milli-Q or Hyclone water.
- 19 **Remove** the magnetic **separation plate** from the shaking incubator and adjust temperature to **51°C** at **600rpm**.
- 20 **Centrifuge** magnetic separation plate at **240 × g** for  **00:01:00** at RT.

1m



- 21 In the **fume hood**, insert magnetic separation plate into **handheld magnetic plate washer** and ensure it is securely locked, allow  00:01:00 to allow magnetic beads to accumulate on bottom of each well. 1m
- 22 Keep plate inserted in handheld magnetic plate washer at all times for this step: add  100 μL of **1 x wash buffer**, wait  00:00:15 to allow the magnetic beads to accumulate at the bottom of each well, remove solution by quickly **inverting** over a waste container and **gently blot** on several layers of paper towel to remove residual solution. 15s
- 23 Repeat previous step **two more times**.
- 24 Add  50 μL of **pre-amplifier solution** to each well and **seal** with a **foil plate seal**, return to Vortemp and incubate for  01:00:00 at **51°C** at **600rpm** (the minimum time for incubation is  00:45:00 and maximum is  02:00:00). 3h 45m
- 25 **Repeat** steps **21-24** for **amplifier** solution in place of **pre-amplifier** solution.
- 26 **Repeat** steps **21-24** for **label probe** solution in place of **pre-amplifier** solution.
- 27 Prepare **SAPE working reagent** by mixing  3 μL of **SAPE** to  1 mL **SAPE diluent** (scaled to the number of the **wells** on your plates + **10%**), **vortex** to mix and keep at **RT protected** from **light**.
- 28 **Repeat** steps **21-24** then add  50 μL of **SAPE working reagent** to each well and **seal** with a foil plate deal, return to Vortemp and incubate for  00:30:00 at **51°C** at **600 rpm** before turning off Vortemp. **Do not exceed**  00:30:00 **incubation**. 1h
- 29 **Repeat** steps **21-24** with **SAPE wash buffer** in place of regular wash buffer.



- 30 To prepare the plate for analysis, add  130 μL of **SAPE wash buffer** to each well, **seal** the plate with a foil plate seal, **tape** sealed magnetic separation plate down onto a **shaker** and shake **vigorously** (~800rpm) at **RT**, immediately run plate on **Magpix**.
- 31 If analysing more than one plate, keep consecutive plates at room temperature protected from light until required and then prepare as in step 32 and shake vigorously (~800rpm) at RT before analysing on the Magpix instrument.
- 32 Plates can be stored **long-term** in at **4°C**, for reanalysis of stored plates, repeat steps 29-30.