

Jun 26, 2025

Version 2

Quantitative Bloch Decay for Spin-1/2 Nuclei V.2

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge5r7yg47/v2>

Alexander L. Paterson¹

¹National Magnetic Resonance Facility at Madison (NMRFAM), University of Wisconsin-Madison, Madison, WI, United States



NMRFAM Facility

UW-Madison

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5r7yg47/v2>

Protocol Citation: Alexander L. Paterson 2025. Quantitative Bloch Decay for Spin-1/2 Nuclei. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge5r7yg47/v2> Version created by **NMRFAM Facility**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 26, 2025

Last Modified: June 26, 2025

Protocol Integer ID: 221044

Keywords: Quantitative MAS Bloch Decay: Spin-1/2, 1d nmr spectra, acquisition of 1d nmr spectra, quantitative bloch decay, t1 relaxation of each resonance, spin nuclei, d1 relaxation period, 1d mas spectrum, nucleus with quantitative intensity, t1 relaxation, nuclei, nucleus, spin, spectrum, resonance

Funders Acknowledgements:

National Science Foundation

Grant ID: 1946970

Abstract

Purpose

To acquire a 1D MAS spectrum of a spin-1/2 nucleus with quantitative intensities.

Scope

The acquisition of 1D NMR spectra of spin-1/2 nuclei with quantitative relative intensities requires attention to be paid to the T_1 relaxation of each resonance. This SOP should only be followed when quantitative relative intensities are required. If only lineshapes or chemical shifts are required, the **d1** relaxation period can be reduced in duration.

Guidelines

If there is a significant rotor or probe background, an echo sequence may be more appropriate.

The quantitative intensities described here are relative intensities, not absolute intensities. Determining absolute intensities, sometimes known as spin counting, is outside the scope of this SOP.

Materials

Definitions:

1. T_1 : Longitudinal relaxation time constant

Safety warnings

-  High-power decoupling can cause damage to a probe if applied incorrectly. Ensure that the decoupling power and duration are set within safe limits.



Before start

User should be familiar with the power limitations and duty cycle of the probe being used.

User should be familiar with the maximum safe spinning speed of the probe and rotor being used.

Accurate T_1 relaxation times and calibrated pulse lengths must have been obtained prior to starting.

The expected amount of time to completion is highly sample-dependent and cannot be accurately estimated ahead of time.

Procedure

- 1 Load pulse sequence **hpdec**.
- 2 Set the relaxation time **d1** to $5 \times T_1$, where T_1 is the largest value previously measured.
- 3 Set the pulse length **p1** and pulse power **plw1** to a previously optimized 90° pulse.
- 4 Set the ^1H decoupling program, **CPDPRG2**, pulse length, **PCPD2**, and pulse power, **plw12**, to a previously optimized ^1H decoupling condition.
 - 4.1 If no ^1H decoupling is required, set the pulse power **plw12** to 0 W.
- 5 Set the acquisition time **aq** to a sufficiently long value such that the FID is not truncated, but less than 50 ms if decoupling is used.

Safety information

Applying high-power decoupling for more than 50 ms can cause physical damage to the probe.

Safety information

Ensure that the relaxation time **d1** is long enough to respect the duty cycle of the probe. This typically requires **d1** to be at least 20x longer than **aq**.

- 6 Adjust **ns** to effectively make use of available experimental time while achieving sufficient signal-to-noise ratios.

Note

If a good number for **ns** cannot be estimated, spectra can be acquired in blocks and added via **fidsum**.





- 7 Process the spectrum while paying careful attention to background subtraction and baseline correction if necessary.
- 8 Obtain the relative intensities of the resonances via integration if well-resolved, or deconvolution if not.

Protocol references

Y. Nishiyama and N. T. Duong, "Practical guides for ^1H detected solid-state NMR under fast MAS for small molecules," Journal of Magnetic Resonance Open, vol. 10-11, p. 100062, Jun. 2022.

<https://doi.org/10.1016/j.jmro.2022.100062>

Protocol

NAME

Saturation Recovery

CREATED BY

NMRFAM Facility

Preview

Example Data at 1.1 GHz: usnan.org/ark:/83454/e1691d38bb-2212-46a2-b1d4-890dc7c444e5

Example Data at 900 MHz: usnan.org/ark:/83454/e17c6d4174-b34a-49fe-bffd-9ad348f39b29

Example Data at 600 MHz: usnan.org/ark:/83454/e11c1749a0-e544-4c54-aad4-adfd78c8e349