

Jun 04, 2019

# Quantitative RT-PCR using the TaqMan® Gene Expression Assays on StepOnePlus™ Real-Time PCR System

 [PLOS One](#)

DOI

[dx.doi.org/10.17504/protocols.io.2v8ge9w](https://dx.doi.org/10.17504/protocols.io.2v8ge9w)

Lisa-Maria Rosenthal<sup>1</sup>, Giang Tong<sup>1</sup>, Katharina Schmitt<sup>1</sup>

<sup>1</sup>Deutsches Herzzentrum Berlin



Lisa-Maria Rosenthal

Deutsches Herzzentrum Berlin

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.2v8ge9w](https://dx.doi.org/10.17504/protocols.io.2v8ge9w)

External link: <https://doi.org/10.1371/journal.pone.0226005>

**Protocol Citation:** Lisa-Maria Rosenthal, Giang Tong, Katharina Schmitt 2019. Quantitative RT-PCR using the TaqMan® Gene Expression Assays on StepOnePlus™ Real-Time PCR System. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.2v8ge9w>

**Manuscript citation:**

Rosenthal L, Leithner C, Tong G, Streitberger KJ, Krech J, Storm C, Schmitt KRL (2019) RBM3 and CIRP expressions in targeted temperature management treated cardiac arrest patients—A prospective single center study. PLoS ONE 14(12): e0226005. doi: [10.1371/journal.pone.0226005](https://doi.org/10.1371/journal.pone.0226005)

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** May 15, 2019

**Last Modified:** June 04, 2019

**Protocol Integer ID:** 23200



## Abstract

### **Quantitative RT-PCR using the TaqMan® Gene Expression Assays on StepOnePlus™ Real-Time PCR System**

## Materials

### MATERIALS

 TaqMan® Gene Expression Assay Thermo Fisher Scientific Catalog #4331182

 StepOnePlus™ Real-Time PCR System Thermo Fisher Scientific Catalog #4376600



1

- Allow the TaqMan® Gene Expression Assay components to thaw on ice

## 2 Primers

Primers			
Gene	Gene Assay ID	Gene	Gene Assay ID
RBM3	Hs00943160_g1	CIRP	Hs00989762_g1
IL-6	Hs00174131_m1	MCP-1 (Ccl2)	Hs00234140_m1
iNOS	Hs01075529_m1	GAPDH	Hs02786624_g1

P

- ## 3
- Prepare the TaqMan® master mix on ice.

Com pone nt	Volume ( $\mu$ L)
Gene Assay Prime rs	0.5
TaqM an Gene Expre ssion PCR Mast er Mix	5.0
cDNA	0.5
Nucle ase- free H <sub>2</sub> O	4.0
Total per reacti on	10.0

- 4 ■ Run the StepOnePlus™ Real-Time PCR System under the following conditions.

Step	Times and Temperatures			
2. PCR Step	Initial Steps		PCR (40 Cycles)	
	AmpErase® UNG Activation	AmpliTaq Gold® DNA Polymerase Activation	Melt	Anneal/Extend
	HOLD	HOLD	CYCLE	
	2 min @ 50 °C	10 min @ 95 °C	15 sec @ 95 °C	1 min @ 60 °C

