


Jun 22, 2023

🌐 Quick and dirty sequencing microbial genome extraction v1

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1qjx7gr2/v1>



Sung won Lim¹, Sebastian Cocioba¹

¹Binomica labs



Sung won Lim

Binomica Labs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1qjx7gr2/v1>

External link: <https://naturepoker.wordpress.com/2023/02/26/halobacteria-mutant-sequencing-3-years-in/>

Protocol Citation: Sung won Lim, Sebastian Cocioba 2023. Quick and dirty sequencing microbial genome extraction v1. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov1qjx7gr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

The protocol works and is currently being improved. Different microbes will require optimization.

Created: June 13, 2023

Last Modified: June 22, 2023

Protocol Integer ID: 83293

Keywords: dirty microbial genome extraction protocol, sequencing microbial genome extraction v1, microbial genome extraction v1, difficult microbes in the future, lab archaeal evolution study, difficult microbe, deinococcus, rapid sequencing kit, certain microbe, such as deinococcus, ont deinococcus, rad004 rapid sequencing kit, read capable dna, using edwards buffer, microbe, adding initial lysozyme incubation step, sequencing project, initial lysozyme incubation step in the very beginning, binomica lab, edwards buffer, capable dna

Abstract

Quick and dirty microbial genome extraction protocol using Edwards buffer.

I've been using this method for in-lab Archaeal evolution study for the last four years or so and it's been working fantastically well, especially for short read sequencing.

Certain microbes (such as *Deinococcus*) can be difficult to process properly using just this protocol - standard caveats and fixes tend to improve the output, such as adding initial lysozyme incubation step in the very beginning, sometimes combined with a freeze-thaw cycle for especially difficult samples (I'll update and upload a separate version of this protocol for difficult microbes in the future).

This particular protocol is aimed at getting as much intact, long-read capable DNA out from a microbe as cheaply as possible, as shown in the included gel picture.

Originally written up for my **lab note** at <https://naturepoker.wordpress.com/2023/02/26/halobacteria-mutant-sequencing-3-years-in/>

Ideation and initial testing performed with Sebastian S. Cocioba at Binomica Labs for *ONT Deinococcus radiophilus* genome sequencing project, using RAD004 rapid sequencing kit.

https://www.ncbi.nlm.nih.gov/assembly/GCF_020889625.1

Image Attribution

Image taken by author - *Halobacterium* mutant strain genome extraction



- 1 Spin down 1 ml of sample for 1 minute at max speed and decant
- 2 Resuspend vigorously with 100ul Edward's buffer
- 3 Transfer carefully to PCR tube - mixture will be viscous
- 4 Add 2ul of RNase A and mix vigorously, vortex for 10 seconds
- 5 Incubate at 37C for 15 minutes
- 6 Add 2ul of Proteinase K and mix vigorously, vortex for 10 seconds
- 7 Incubate at 55C for 1 hour, and deactivate via incubation at 95C for 10 minutes
- 8 Transfer to 1.5ml eppendorf tube
- 9 Add 10% 3M (pH 5.4) sodium acetate, and 1:1 volume of 100% isopropyl alcohol
- 10 Invert tube 10 times - precipitates should begin to form
- 11 Spin down at max speed for 5 minutes
- 12 Decant the supernatant carefully
- 13 Add 1ml of 70% EtOH and resuspend the pellet



- 14 Spin down at max speed for 5 minutes
- 15 Repeat 70% EtOH resuspension and washing step at least 2 more times
- 16 Decant completely and dry the pellet for 5 minutes - do not let the pellet overdry
- 17 Resuspend in storage buffer of choice or dH₂O to wanted volume
- 18 Incubate in a 37C shaker for 1 hour - the extract is likely to be not dissolved fully
- 19 Incubate in 4C overnight - depending on yield, the extract will be extremely viscous
- 20 Check extraction quality using both standard gel electrophoresis and UV-Vis method such as Nanodrop.