

May 29, 2024

## Quick protocol for protein extraction from adherent fish-derived fibroblasts

DOI

[dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1)

Joao M Moreno<sup>1,2</sup>, Vitor C Sousa<sup>1</sup>, Romana Santos<sup>2</sup>

<sup>1</sup>cE3c – Centre for Ecology, Evolution and Environmental Changes & CHANGE – Global Change and Sustainability Institute, Faculdade de Ciências da Universidade de Lisboa, Lisboa, Portugal;

<sup>2</sup>MARE – Centro de Ciências do Mar e do Ambiente (MARE) & ARNET—Aquatic Research Network, Faculdade de Ciências da Universidade de Lisboa, Lisboa, Portugal



**Joao M Moreno**

cE3c – Centre for Ecology, Evolution and Environmental Chang...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1>

**Protocol Citation:** Joao M Moreno, Vitor C Sousa, Romana Santos 2024. Quick protocol for protein extraction from adherent fish-derived fibroblasts. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 16, 2024

**Last Modified:** May 29, 2024

**Protocol Integer ID:** 99955

**Keywords:** Protein extraction, Protein purification, Cell culture, total proteins from adherent fish, using fish cell culture, extracting total protein, fish cell culture, quick protocol for protein extraction, protein extraction, comprehensive solubilization of cellular protein, total soluble protein, cellular protein, proteomic, studying proteomic, cells from the culture flask, adherent fish, derived fibroblast, quality protein, fish, mass spectrometry, extraction

**Funders Acknowledgements:**

Fundação para a Ciência e a Tecnologia

Grant ID: SFRH/BD/143199/2019

Fundação para a Ciência e a Tecnologia

Grant ID: PTDC/BIA-EVL/4345/2021

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://dx.doi.org/10.17504/protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://dx.doi.org/10.17504/protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

We present a rapid and efficient protocol for extracting total proteins from adherent fish-derived fibroblasts, specifically optimized for applications in Western blotting and mass spectrometry. The method involves directly dissociating cells from the culture flask into a lysis buffer (e.g. RIPA buffer) followed by centrifugation to collect the total soluble proteins. The use of buffers like RIPA ensures comprehensive solubilization of cellular proteins while preserving their integrity and functionality. This straightforward and reproducible protocol yields high-quality protein extracts suitable for various downstream analytical techniques. Its simplicity and reliability make it an invaluable tool for researchers studying proteomics using fish cell culture.



## Materials

### Solutions

- 1X Phosphate Buffered Saline (PBS)
- RIPA lysis buffer: 50 mM Tris-HCl (pH 7.4), 150 mM NaCl, 1mM EDTA, 1% NP-40, 1% sodium deoxycholate, 0.1% SDS, sterile-filtered.

### Materials

- 2 mL or 1.5 mL microcentrifuge tubes
- Cell scrappers
- Pipettes and pipette tips
- Serological pipettes (variable volume)
- Pasteur pipettes


### Equipments

- Refrigerated microcentrifuge
- Flow hood chamber

## Troubleshooting



- 1 Carefully remove all culture media from the flask and add enough ice-cold 1X PBS to wash the cells
- 2 Carefully remove the ice-cold 1x PBS and add ice-cold lysis buffer (RIPA buffer) according to the estimated number of cells:  

 1 mL for  $10^7$  cells (roughly a T75 flask).
- 3 Use a cell scraper to dissociate the cells from the bottom of the flask.
- 4 Resuspend the cells in the lysis buffer and transfer the suspension to a microcentrifuge tube.
- 5 Agitate for 20 minutes at 4°C.
- 6 Centrifuge at 13,000 x g for 20 min at 4°C.
- 7 Carefully transfer the supernatant containing the soluble protein to a new tube and keep on the ice. Discard the pellet.

Note: The protein solution can be kept in the freezer for longer storage periods until further use.