

Oct 30, 2023

Quick-Start Protocol for DNeasy® PowerSoil® Pro Kit

DOI

dx.doi.org/10.17504/protocols.io.j8nlkww45l5r/v1



Carlos Carlos Goller¹

¹NC State University

BIT-Protocols

Tech. support phone: +91 95134-135 email: ccgoller@ncsu.edu



Nidhi Grover

North Carolina State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkww45l5r/v1>

Protocol Citation: Carlos Carlos Goller 2023. Quick-Start Protocol for DNeasy® PowerSoil® Pro Kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlkww45l5r/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 18, 2022

Last Modified: February 28, 2024

Protocol Integer ID: 71497

Keywords: PowerSoil, metagenomics, soil, start protocol, protocol, kit

Abstract

This protocol will allow users to operate the DNeasy® PowerSoil® Pro Kit.

Guidelines

Solution CD2 should be stored at 2–8°C upon arrival. All other reagents and kit components should be stored at room temperature (15–25°C).

Further information

- DNeasy® PowerSoil® Pro Kit Handbook: www.qiagen.com/HB-2495
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Materials

The DNeasy® PowerSoil® Pro Kit. All requirements included within.

Troubleshooting

Safety warnings

- ! Use caution when operating the vortex. Ensure tubes are balanced before turning on the instrument. Dispose of materials in appropriate waste containers.



Before start


- Ensure that the PowerBead Pro Tubes rotate freely in the centrifuge without rubbing.
- If Solution CD3 has precipitated, heat at 60°C until the precipitate dissolves.
- Perform all centrifugation steps at room temperature (15–25°C).



Using the DNeasy® PowerSoil® Pro Kit

18m 10s

- 1 Spin the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom. Add up to  250 mg of soil and  800 µL of Solution CD1. Vortex briefly to mix.



- 2 Secure the PowerBead Pro Tube horizontally on a Vortex Adapter for 1.5–2 ml tubes (cat. no. 13000-V1-24). Vortex at maximum speed for  00:10:00

10m

Note

If using the Vortex Adapter for more than 12 preps simultaneously, increase the vortexing time by 5-10 min

For more information about other bead beating methods, see the “Protocol: Detailed” section of DNeasy® PowerSoil® Pro Kit Handbook.

- 3 Centrifuge the PowerBead Pro Tube at  15000 x g for  00:01:00

1m







- 4 Transfer the supernatant to a clean 2 ml Microcentrifuge Tube (provided)

Note

Expect 500–600 µl. The supernatant may still contain some soil particles

- 5 Add  200 µL of Solution CD2 and vortex for  00:00:05

5s

- 6 Centrifuge at  15000 x g for  00:01:00 at  Room temperature . Avoiding the pellet, transfer up to  700 µL of supernatant to a clean 2 ml Microcentrifuge Tube (provided)



1m






Note

Expect 500–600 µl.



7 Add  600 μL of Solution CD3 and vortex for  00:00:05 .

5s

8 Load  650 μL of the lysate onto an MB Spin Column and centrifuge at  15000 x g for  00:01:00

1m






9 Discard the flow-through and repeat step 8 to ensure that all of the lysate has passed through the MB Spin Column.

10 0. Carefully place the MB Spin Column into a clean 2 ml Collection Tube (provided).

Safety information




Avoid splashing any flow-through onto the MB Spin Column

11 Add  500 μL of Solution EA to the MB Spin Column. Centrifuge at  15000 x g for  00:01:00

1m



12 . Discard the flow-through and place the MB Spin Column back into the same 2 ml Collection Tube.

13 Add  500 μL of Solution C5 to the MB Spin Column. Centrifuge at  15000 x g for  00:01:00

1m





14 Discard the flow-through and place the MB Spin Column into a new 2 ml Collection Tube (provided).



15 Centrifuge at up to  16000 x g for  00:02:00 . Carefully place the MB Spin Column into a new 1.5 ml Elution Tube (provided)

2m



16 . Add  50 mL -  100 mL of Solution C6 to the center of the white filter membrane.



- 17 Centrifuge at  15000 x g for  00:01:00 . Discard the MB Spin Column. The DNA is now ready for downstream applications

1m

**Note**

We recommend storing the DNA frozen (–30 to –15°C or –90 to –65°C) as Solution C6 does not contain EDTA. To concentrate DNA, please refer to the Troubleshooting Guide.