

Oct 14, 2022

Radiolabeled polyamine uptake in cells

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1>

Marine Houdou¹, Nathalie Jacobs¹, Peter Vangheluwe¹

¹KU Leuven



Nathalie Jacobs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1>

External link: <https://doi.org/10.3390/biom13020337>

Protocol Citation: Marine Houdou, Nathalie Jacobs, Peter Vangheluwe 2022. Radiolabeled polyamine uptake in cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1>

Manuscript citation:

Houdou M, Jacobs N, Coene J, Azfar M, Vanhoutte R, Haute CVd, Eggermont J, Daniëls V, Verhelst SHL, Vangheluwe P, Novel Green Fluorescent Polyamines to Analyze ATP13A2 and ATP13A3 Activity in the Mammalian Polyamine Transport System. *Biomolecules* 13(2). doi: [10.3390/biom13020337](https://doi.org/10.3390/biom13020337)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 04, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 70812

Keywords: ASAPCRN, radiolabeled polyamine uptake capacity, polyamine uptake capacity in cell, polyamine uptake in cell, polyamine uptake, dpm

Abstract

This protocol provides a technique to determine the radiolabeled polyamine uptake capacity in cells, via the acquisition of disintegrations per minute (DPM) using a Liquid Scintillation Counter.

Guidelines

Proper guidelines for working with radiolabeled materials should be followed at all times.

Safety warnings

 Radiation hazards



- 1 Cells are seeded in 12-well plates, such that 70-80% confluency is reached on the day of the assay. Seed out 2 'treatment' wells, 1 'quick wash' well and 1 'untreated' well per cell line and treatment dose.

Note

	A	B	C	D
	treatment	treatment	quick wash	untreated

- 2 Remove the culturing medium in 'treatment wells' and add  500 μL of medium, containing the desired concentration of radiolabeled polyamines. Leave regular culturing medium in the 'quick wash' and 'untreated' wells.

- 3 Incubate  37 $^{\circ}\text{C}$  00:30:00

30m

- 4 4 minutes before reaching the 30-minute mark, perform a Quick-wash step as follows: add  500 μL of medium, containing the desired concentration of radiolabeled polyamines in the 'quick wash' wells.  00:00:00 immediately remove the treatment medium, and wash with  500 μL of DPBS (-/-) containing  50 micromolar (μM) of the respective unlabeled (=cold) polyamine at  4 $^{\circ}\text{C}$. Continue with 2 more washing steps with  800 μL of DPBS (-/-)  4 $^{\circ}\text{C}$.

- 5 Proceed with washing the other wells ('treatment' and 'untreated'). Remove the treatment medium, and wash with cold  500 μL of DPBS (-/-) containing  50 micromolar (μM) of the respective cold polyamine at  4 $^{\circ}\text{C}$. Continue with 2 more washing steps with  800 μL of DPBS (-/-) at  4 $^{\circ}\text{C}$.

- 6 After the last wash, add  200 μL 0.1% SDS-DPBS (-/-) per well to lyse the cells.



7 Incubate  00:10:00 at  Room temperature .

10m

8 Scrape the cells and pipette the lysates into scintillation vials, that are filled with  7 mL of scintillation fluid (Ecolite (+), MP). Rinse each well with  200 μ L DPBS (-/-) and pipette into the respective scintillation vial.

9 Mix scintillation vials well prior acquisition of disintegrations per minute (DPM) in the Liquid Scintillation Counter.