



Jul 09, 2024

Radiolabeled spermine uptake in cells

 Forked from [Radiolabeled polyamine uptake in cells](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkorm5v5r/v1>

Stephanie Vrijzen^{1,2}, Marine Houdou³, Nathalie Jacobs³, Peter Vangheluwe^{1,2}

¹Laboratory of Cellular Transport Systems, Department of Cellular and Molecular Medicine, KU Leuven, B-3000 Leuven, Belgium;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network Chevy Chase, MD 20815, USA;

³KU Leuven



Stephanie Vrijzen

Laboratory of Cellular Transport Systems

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkorm5v5r/v1>

External link: <https://doi.org/10.3390/biom13020337>

Protocol Citation: Stephanie Vrijzen, Marine Houdou, Nathalie Jacobs, Peter Vangheluwe 2024. Radiolabeled spermine uptake in cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkorm5v5r/v1>

**Manuscript citation:**

Houdou M, Jacobs N, Coene J, Azfar M, Vanhoutte R, Haute CVd, Eggermont J, Daniëls V, Verhelst SHL, Vangheluwe P, Novel Green Fluorescent Polyamines to Analyze ATP13A2 and ATP13A3 Activity in the Mammalian Polyamine Transport System. Biomolecules 13(2). doi: [10.3390/biom13020337](https://doi.org/10.3390/biom13020337)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 05, 2023

Last Modified: July 09, 2024

Protocol Integer ID: 91849

Keywords: ASAPCRN, radiolabeled spermine uptake, spermine uptake in cell, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000458

Abstract

This protocol provides a technique to determine radiolabeled spermine uptake in cells, via the acquisition of counts per minute (CPM) using a Liquid Scintillation Counter.

Guidelines

Proper guidelines for working with radiolabeled materials should be followed at all times.

Materials

- ¹⁴C-labeled spermine: 3139-50 µCi, ARC
- Cold spermine: 85590-5G, Merck
- RIPA Lysis and Extraction Buffer: 89900, Invitrogen
- EcoLite Liquid Scintillation Cocktail: 01882475-CF, MP Biomedicals

Safety warnings

 Radiation hazards



- 1 Cells are seeded in 12-well plates, such that 70-80% confluency is reached on the day of the assay. Seed out 2 'treatment' wells and 1 'background' well per cell line.
- 2 Remove the culture medium from all wells.
- 3 In the treatment wells: add  300 μL of medium, containing 5 μM radiolabeled spermine.
In the background wells: add  300 μL of medium, containing 5 μM radiolabeled spermine and 100 μM unlabeled (cold) spermine.
- 4 Incubate  37 $^{\circ}\text{C}$  00:30:00 30m
- 5 Aspirate the medium.
- 6 Wash wells 2x with  1 mL PBS (-/-).
- 7 Remove the last wash, then add 200 μL RIPA buffer to the wells.
- 8 Incubate  00:10:00 at  Room temperature 10m
- 9 Scrape the cells and pipette the lysates into scintillation vials, that are filled with  7 mL of scintillation fluid (Ecolite (+), MP). Rinse each well with  200 μL PBS (-/-) and pipette into the respective scintillation vial.
- 10 Mix scintillation vials well prior acquisition of counts per minute (CPM) in the Liquid Scintillation Counter.