

Jun 07, 2023



RADSeq protocol EFGL



Forked from BestRAD protocol

DOI

dx.doi.org/10.17504/protocols.io.kqdg39bjqg25/v1



¹Eagle Fish Genetics Lab

Mimulus



EagleFish GeneticsLab

Eagle Fish Genetics Lab



Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: https://dx.doi.org/10.17504/protocols.io.kqdg39bjqg25/v1

Protocol Citation: EagleFish GeneticsLab 2023. RADSeq protocol EFGL. protocols.io

https://dx.doi.org/10.17504/protocols.io.kqdg39bjqg25/v1

Manuscript citation:





License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 16, 2023

Last Modified: June 07, 2023

Protocol Integer ID: 77120

Keywords: Massively parallel sequencing, MPS, restriction-site associated DNA, RAD, sequence capture, genotyping, population genetics, Eagle Genetics Lab, EFGL, Radseq, radseq protocol efgl modified, radseq protocol efgl modified from protocol, protocol

Abstract

Modified from protocol of Sean O'Rourke and Mike Miller published in:

Citation

Omar A. Ali, Sean M. O'Rourke, Stephen J. Amish, Mariah H. Meek, Gordon Luikart, Carson Jeffres and Michael R. Miller

(2016). RAD Capture (Rapture): Flexible and Efficient Sequence-Based Genotyping. GENETICS.

https://doi.org/10.1534/genetics.115.183665

LINK

Attachments





BestRAD protocol v1.... BestRAD_library_prep...

460KB

75KB



Materials

MATERIALS

NEBNext Multiplex Oligos for Illumina (Index Primers Set 1) - 24 rxns New England Biolabs Catalog #E7335S

Pstl - 10,000 units New England Biolabs Catalog #R0140S

🔀 Bfal - 500 units New England Biolabs Catalog #R0568S

X NEBNext Ultra II DNA Library Prep Kit for Illumina - 24 rxns New England Biolabs Catalog #E7645S

Restriction enzymes:

Pstl: NEB Bfal-HF: NEB CutSmart buffer

BestRAD plate adaptors:

BestRAD adaptors allow the addition of inline barcodes, and the isolation of RAD tags through purification by and enzymatic liberation of fragments from streptavidin beads. Well-specific Hamming barcodes (septamers in this case) are specified in an Excel spreadsheet and are not actually ordered as N's. Top oligos are 5'-biotinylated and contain the 3' Pstl/Sbfl overhang. Bottom oligos are 5'-phosphorylated to promote ligation.

Top oligo

Anatomy: biotin spacer SbfI barcode overhang Sequence: /5Biosg/GTACGT CCTGCAGG NNNNNNN TGCA

Bottom oligo

Sequence: /5Phos/NNNNNN CCTGCAGG ACGTAC

Ordering Specs

Integrated DNA technologies

- Standard plate oligos
- 25 nmol standard desalting
- drv
- \$1859.52 total for a set of 48



Streptavidin purification:

Dynabead M280 strepatividin, 2mL (GrizMart, ~ \$500)

2X Binding and Wash Buffer:

- 1. 10 mM TrisHCI (pH 7.5)
- 2. 1 mM EDTA pH 8.0
- 3. 2 M NaCl
- 4. Concentrated (or dry) stocks should be available in chemical cabinet.

NextGen library prep:

NEBNext Ultra II

Kit

- includes reagents for end-repair, A-tail, ligation
- NEB E7645S 24 rxns, GrizMart, Fisher
- **~** \$590

Oligos (12-plex)

- Indexed oligos containing Illumina sequencing primer sequences and required for annealing to flow cell. Added via PCR to NEBNext adaptor-ligated fragments.
- 12 barcoded i7 indexing primers
- 1 universal (i5) oligo
- NEBNext adaptor w/ USER enzyme
- NEB E7335S
- Grizmart
- **~**\$110

Universal primer with molecular barcode

- Modeled after i5 index primers from NEB #E7600 (p21 of manual)
- N's specify the equimolar addition of dATP, dTTP, dGTP, and dCTP during synthesis
 - Not truly random. Some GC-bias in addition
 - Should be sufficiently diverse to detect PCR duplicates
- Used in place of the universal i5 oligo in NEB #E7335S
- IDT, HPLC-cleaned \$90.50
- Primer with molecular barcode in bold:

5'-AAT GAT ACG GCG ACC ACC GAG ATC TAC ACN NNN NNN NAC ACT CTT TCC CTA CAC GAC GCT CTT CCG ATC*T-3'

Compare to NEBNext i506 primer with indexing barcode in bold:



5´-AAT GAT ACG GCG ACC ACC GAG ATC TAC ACT AAT CTT AAC ACT CTT TCC CTA CAC GAC GCT CTT CCG ATC*T-3

Reaction purification and size selection:

- Ampure XP magnetic beads
- Polyethylene glycol (PEG-8000)
- NaCl
- Nice summary as of Dec 2017: http://core-genomics.blogspot.com/2012/04/how-do-spribeads-work.html

Troubleshooting

Safety warnings

• Please refer to the SDS (Safety Data Sheet) for safety warings and hazard information.

Annealing TOP/BOTTOM BestRAD adapters



1 **Preparing BestRAD adapters** (Skip section if BestRAD adapters have been previously annealed)

Note

We have used 3 restriction enzymes with 6 base pair recognition sites (Pstl, BamHl, and HindIII) and 1 enzyme with an 8 bp recognition site (Sbfl) in our protocol. Adapters should be compatible with the base pair overhang that each restriction enzyme leaves. Both Pstl and Sbfl leave the same overhangs so the same adapters can be used with either enzyme. Adapters are ordered from Integrated DNA Technologies as single stranded top adapters in one plate and single stranded bottom adapters in another plate. They have to be rehydrated, diluted to the proper concentration, and annealed prior to beginning the BestRAD protocol. See attached documents below for the sequences of our top and bottom adapters that we use. Note that the file format is for importing into IDTdna's Custom DNA Oligo, single-stranded DNA, 96 well plate form (https://www.idtdna.com/pages/products/custom-dna-rna/dna-oligos/custom-dna-oligos). We've ordered adapters a variety of ways, but the best/easiest is to order the

adapters at 25nmole scale, standard desalting purification, PCR plate type, dry ship option,

normalized yield and nmol quantity of 10. All top adapters are biotinylated.

We typically make a 1μ M working concentration of annealed adapters if the restriction enzyme we are using is a 6 base-pair cutter. We make a 50nM working concentration of annealed adapters if the restriction enzyme is an 8 base-pair cutter. A 6 bp cutter will cut more frequently so we need a higher concentration of it.





If you already have a plate of diluted, annealed adapters, skip to Step 5

Equipment and supplies needed for steps 1-4:

thermalcycler (We have several different types of thermalcyclers. Here's an example of one we use often)

Equipment

Mastercycler Pro

NAME

thermalcycler

TYPE

Eppendorf

BRAND

No Longer Manufactured

SKU



Equipment

EPPENDORF SCIENTIFIC CENTRIFUGE 5804

NAME

plate centrifuge

TYPE

Eppendorf

BRAND

02-262-8153PM; discontinued

SKU





Equipment	
Repeater-M4	NAME
pipette	TYPE
Eppendorf	BRAND
14-287-150	SKU

Equipment NAME ALPS 25 V TYPE manual heated plate sealer BRAND Thermo Scientific SKU AB-0384/110; discontinued

p10 multi-channel pipette and tips

benchtop centrifuge

1 mL Eppendorf combi-tips (Fisher Scientific; 13-683-703)

heat seals (Fisher Scientific; AB-0745)

unskirted 96 well PCR trays (Fisher Scientific; AB-0700)

Reagents needed:

Bottom adapter plate

Top adapter plate

Duplex Buffer IDT Catalog #11-01-03-01

Main Invitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36

1.1 Making a **100μM** stock of bottom and top adapters



Assuming the bottom and top adapters were ordered dry and at a 25 nmole scale and normalized to 10nmol, first re-hydrate to **100µM** by adding 4 100 µL of nuclease-free water to each well with repeater pipette and combitip

Heat seal the plate

Expected result

Now you should have a plate of Δ 100 μL top adapters at 100μM concentration and a plate of 4 100 µL bottom adapters at **100µM** concentration.

1.2 Making a **10uM** stock of bottom and top adapters

Note

 $C1 = 100 \mu M$

 $C2 = 10 \mu M$

 $V2 = 100 \mu l$

C1V1 = C2V2

V1 = (C2V2)/C1

V1 = (10 * 100)/100

V1 = 10

1.3 Mix the **100µM** adapter plates by gently vortexing, then briefly spin the plates down



- 1.4 Use a p10 multi-channel pipette to transfer Δ 10 μL of each top adapter to a new plate.
- 1.5 Using a repeater pipette and combitip, add \perp 90 μ L of nuclease-free water to the same wells of the new plate.

Heat seal the plate



Even though you need to use the plates immediately in the next steps, you still need to heat seal them after every dilution so you can vortex the plate and spin it down

1.6 Repeat for the bottom adapters from step 1.4.

Heat seal the plate

Expected result

Now you should have a plate of Δ 100 μL top adapters at **10μM** concentration and a plate of \perp 100 µL bottom adapters at **10µM** concentration.

- 2 Choose between the following options:
 - If using 6 base pair cutter (Pstl, BamHl, Hindlll), need to make **1µM** dilution of BestRAD adapters
 - If using 8 base pair cutter (SbfI), need to make 50nM dilution of BestRAD adapters

STEP CASE

6 base pair cutter, 1µM adapters From 22 to 212 steps

3 Making a **1.0µM** working concentration of **COMBINED** bottom and top adapters

Note

 $C1 = 10 \mu M$

 $C2 = 1.0 \mu M$

 $V2 = 50 \mu l$

C1V1 = C2V2

V1 = (C2V2)/C1

V1 = (1 * 50)/10

V1 = 5

3.1 Mix the **10µM** adapter plates by gently vortexing, then briefly spin the plates down



- 3.2 Using a multi-channel pipette, add $\underline{\underline{}}$ 2.5 μL of the **10\mu M** top adapters to a new plate. Discard tips.
- 3.3 Using a multi-channel pipette, add \triangle 2.5 μ L of the corresponding **10\muM** bottom adapters (the same well position as the top adapter plate) to the same well of the new plate. Discard tips.
- 3.4 Mix the duplex buffer by vortexing, then briefly spin down
- 3.5 Using a repeat pipette and a 1mL combi-tip, add \perp 45 μ L of duplex buffer to all wells.

Heat seal the plate

Expected result

Now you have a single plate with $\Delta 50 \mu$ of bottom and top adapters combined at **1µM** concentration.

4 Annealing adapters

3m

Run the plate of **COMBINED** adapters ($_____$ 50 μ L) on a thermalcycler using the following program (In our lab, the program can be found using the menu options General \rightarrow RAD \rightarrow **duplexAssay**)

CHECK THE PROGRAM TO MAKE SURE IT'S CORRECT


- 1. § 95 °C for (5) 00:02:00
- 2. Reduce temperature by \$\mathbb{8} 1 \cdot \Cong for \cdot \Omega 00:01:00
- 3. repeat until temperature is down to \$\ \ 4 \ ^C\$



Expected result

Now you have a single plate with $\[\] 50 \ \mu \]$ of bottom and top adapters combined at either $1\mu M$ concentration if you were following the 6 bp cutter steps or 50nM concentration if you were following the 8 bp cutter steps. These BestRad adapters include the barcodes for each sample

Check to see if there's sufficient stock made of **2X** B+W buffer (needed for step 15, size selection) and stored in 50mL conical tube at room temperature. If not, make it



- 5.1 Add the following to a 50mL conical tube:
 - 1. \perp 500 μ L 1M Tris-HCl pH 7.5
 - 2. \perp 100 μ L 0.5M EDTA pH 8.0
 - 3. <u>▲</u> 5.844 g NaCl
 - 4. Add enough water to obtain 🚨 50 mL final volume

Note

Final concentration is 2X: 10mM Tris-HCl pH 7.5, 1mM EDTA pH 8.0, 2M NaCl

6 Choose from the following options:
Quantification of <48 samples
Quantification of >48 samples

CTED CACE			
OLEE CHOE			

<48 Samples 202 steps

DNA quantification <48 samples

7 Quantify DNA with Qubit

Note

Note that for sample sizes less than 48ish, we quantify each sample individually using a Qubit fluorometer. Otherwise, we use a plate reader.



Note **Equipment and supplies needed for this step: Equipment** NAME **Qubit Fluorometer** TYPE Fluorometer to quantify DNA **BRAND** Invitrogen SKU Q33226 Qubit assay tubes (Fisher Scientific; Q32856) Eppendorf 2mL Snap-Cap Microcentrifuge Safe-Loc Tubes (Fisher Scientific; 05-402-7) unskirted 96 well PCR trays (Fisher Scientific; AB-0700) p10 pipette and tips p200 pipette and tips Reagents needed: DNA tray(s) Qubit dsDNA HS assay kit Fisher Scientific Catalog #Q32854 Main Invitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36 if dilution necessary

- 7.1 Enter sample data info into "sampleData" tab of BestRAD_library_prep_labData v1.2.xlsx (attached in the DESCRIPTION tab of this protocol. It has example data for 20 samples.)
- 7.2 Qubit all samples with Qubit dsDNA HS assay kit.

 Follow the manufacturer's protocol to make a working solution (1:200) of dye https://tools.thermofisher.com/content/sfs/manuals/Qubit_dsDNA_HS_Assay_UG.pdf
- 7.3 Make sure you set the sample volume and the output to be read in $ng/\mu l$, not ng/mL



Take **2** readings per sample and enter them into columns F and G in "qubitData" tab in BestRAD_library_prep_labData v1.2.xlsx. The spreadsheet has a column that will take the average. The target is **10ng/μl**.

- 7.4 If DNA concentration is off the Qubit scale,
 - 1. perform a 1:5 dilution (\perp 1 μ L of DNA, \perp 4 μ L of water) in a new PCR plate
 - 2. Re-Qubit diluted samples
- 7.5 If DNA concentration is too low (less than 9ng/µl), reach out to lead biologist for replacement samples, or re-extract DNA then re-Qubit new samples

DNA Normalization

8 Normalize DNA to 10ng/μl

Note

Equipment and supplies needed for this step:

manual heated plate sealer
plate seals (Fisher Scientific; AB-0745)
unskirted 96 well PCR trays (Fisher Scientific; AB-0700)
p10 pipette and tips
p200 pipette and tips
vortex
plate centrifuge
Results from "Norm" tab of BestRAD_library_prep_labData v1.2.xlsx
Reagents needed:
DNA tray(s)


Invitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36

The table in "Norm" tab of $BestRAD_library_prep_labData\ v1.2.xlsx$ has been formulated to normalize all samples to $10ng/\mu l$ at a $30~\mu L$ final volume, and you'll need $10~\mu L$ for each sample in next step.

- 8.1 Label a new PCR plate: Normalized DNA for [insert RADseq project name], date, initials
- 8.2 Pipette the amount of DNA indicated in column F of "Norm" tab of BestRAD_library_prep_labData v1.2.xlsx from the DNA tray to corresponding well position in new Normalization plate



- 8.3 Pipette the amount of nuclease-free water indicated in column G of "Norm" tab to corresponding well position in new Normalization plate
- 8.4 Heat seal plate, Mix by gently vortexing, then briefly spin the plate down



Day 1 - Digestion

9 Digest the normalized DNA with a restriction enzyme



Equipment and supplies needed for this step:

thermal cycler
manual heated plate sealer
plate seals (Fisher Scientific; AB-0745)
p10 pipette and multi-channel pipette and tips
p200 pipette and tips
Repeater-M4 pipette
0.1mL combi-tip (Fisher Scientific; 13-683-700)
unskirted 96 well PCR trays (Fisher Scientific; AB-0700)
1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific; 13-698-791)
vortex
plate centrifuge
timer

Reagents needed:

10µl of normalized DNA per sample

- Invitrogen Nuclease-Free water **Fisher Scientific Catalog** #43-879-36
- CutSmart Buffer 5.0 ml New England Biolabs Catalog #B7204S
- 1 of the following restriction enzymes
- Pstl 50,000 units New England Biolabs Catalog #R0140L
- SbfI-HF 2,500 units New England Biolabs Catalog #R3642L
- BamHI 10,000 units New England Biolabs Catalog #R0136S
- HindIII 10,000 units New England Biolabs Catalog #R0104S

NOTE: All restriction enzymes' concentrations are 20,000U/µl

- 9.1 Thaw reagents to prep for master mix
 - Whichever restriction enzyme chosen for the project (Pstl, Sbfl, BamHl, or HindIII)
 - CutSmart Buffer 10X
 - Normalized DNA if it's been in the freezer
- 9.2 Add reagents from table below to a 1.5mL Lo-Bind tube to make a master mix. Add 5-10 to your sample size to account for pipette error. For example, if you have 20 samples, multiply the volumes per sample by 25.



ample
1
1

Recipe above assumes 20units/µl of enzyme and 2.4 units per reaction

- 9.3 Vortex and briefly spin down
- 9.4 Label PCR plate "Digestion"
- 9.5 Using a p10 pipette or 0.1mL combitip, add 🚨 2 μL of master mix to each well
- 9.6 Using a multi-channel p10 pipette, add \perp 10 μ L of normalized DNA to corresponding wells
- 9.7 Heat seal PCR plate, vortex gently and quickly spin down.
- 9.8 Run the following thermalcycler program (In our lab, the program is "1RE-digest" and can be found using the menu options General → RAD)

1h 30m

- 1. 37 °C for (5) 01:00:00
- 2. \$\\$5 \circ\$ for \(\chi_{\chi}\) 00:30:00

CHECK THE PROGRAM TO MAKE SURE IT'S CORRECT

9.9 Set a timer for (5) 01:15:00

1h 30m

Prep next step, about 600:15:00 prior to end of thermalcycler program




Expected result

Tray labeled "Digestion" with \perp 12 μ L of \sim \perp 100 ng of DNA per sample that has been digested with a restriction enzyme

Day 1 - Ligation

10 Ligate adapters onto digested DNA



- 10.1 1. Thaw BestRAD adaptors (prepared in steps 1-4) and reagents for Ligation Master Mix (listed above)
- 10.2 gently vortex and spin down all reagents (
- 10.3 Make Ligation master mix in a 1.5mL Lo-Bind tube by combining reagents below. Multiply by the number of samples you have + 5-10 extra to account for pipette error.



Ligation Master Mix	Per sample		
Water	1.28 µl		
CutSmart Buffer 10X	0.4 µl		
ATP	0.16 µl		
T4 ligase	0.16 µl		
Total	2 µl		

- 10.4 Remove Digestion PCR plate from thermocycler (Step 9.9), quick spin
- 10.5 With a multi-channel p10 pipette, add Δ 2 μL of annealed adaptors to each well in Digestion PCR plate, corresponding to the same well setup in the "library prep" tab of BestRAD_library_prep_labData v1.2.xlsx. Discard tip after each use.
- 10.6 Using a p10 pipette, add $\stackrel{\blacksquare}{\Delta}$ 2 μ L of Ligation master mix to each well. Discard tip after each use.
- 10.7 Heat seal PCR plate, gently vortex, quick spin
- 10.8 Overnight

16h 30m

Run thermalcycler program: (In our lab, the program is called "**2 Ligation**" and can be found using the menu options : General \rightarrow RAD)

- 1. \$\colon 20 \colon \colon 16:00:00
- 2. \$\mathbb{8}\$ 80 °C for \(\bigotimes \) 00:30:00

Note


Do not let plate sit at 4°C hold for too long, recommended to proceed to next step as soon as program ends.



Expected result △ 16 μL of digested DNA per sample with BestRad adapters attached

Day 2 - Sonication

11 **Prepping the Sonicator**



11.1 Fetch ~ 🗸 2 L of de-ionized water from the Fish Health Lab. Take containers with you





11.2 Connect all corresponding tubes to the machine as labeled and screw on the filter




11.3 Empty 🚨 1 L of de-ionized water into sonicator and start the machine.





11.4 When the filter is fully saturated, and water cycle is complete, add more de-ionized water until the level is just below the top white spout.



11.5 Proceed to next step, and let machine cool to \$\ 4 \circ\$



Day 2 - Pool samples

12

Note

Equipment and supplies needed for this step:

p200 pipette and tips

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific; 13-698-791)

2.0mL Lo-Bind microcentrifuge tube (Fisher Scientific; 13-698-792)

Reagents needed:

"Digestion" plate of digested samples with BestRad adapters attached

12.1 **Pooling samples**

- Using a p200 single channel pipette, pool all 🚨 16 μL of each well from Digestion PCR plate into a new 1.5mL Lo-Bind tube
- If library contains more than 93 samples, pool into a 2.0mL Lo-Bind tube instead

Expected result

A library tube with pooled samples, volume = \perp 16 μ L X number of samples

Day 2 - Sonicate, precipitate, ligate library

13 purifying pooled library



Equipment and supplies needed for this step:

EquipmentDynaMag SpinNAMEMagnet for microcentrifuge tubesTYPEInvitrogenBRAND12320DSKU

p200 pipette and tips

p1000 pipette and tips

1.5mL Lo-Bind microcentrifuge tubes (Fisher Scientific; 13-698-791)

Bioruptor Plus TPX microtubes (Fisher Scientific; C30010010-50)

5.0mL Axygen Polypropylene Snaplock non-sterile tube (Fisher Scientific; 14-568-100) OR

15mL conical tube (Fisher Scientific; 14-959-53A)

vortex

benchtop centrifuge

Reagents needed:

tube with pooled library

Ampure XP beads **Beckman Coulter Catalog** #A63881

Minvitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36

100% ethanol (laboratory grade, non-denatured!)

Note

Ampure XP beads have been aliquoted into 2 mL tubes to prevent contamination of the stock. In general, whenever they are used in the protocol, they should be vortexed thoroughly to ensure mixing



if library contains >18 samples

- Rule of thumb: Each tube should have 200-300µl of pooled samples to allow adequate room for 80% ethanol wash step
- At the end, the combined total low TE volume of all tubes must be \perp 210 μ L

Using the table below, calculate the respective volumes required, or enter in number samples in "Pooling_calculator" tab in library prep spreadsheet to have volumes calculated

		N = number of samples			
No. Samples	Tubes needed	Vol. of pooled samples	Vol. of AMPure XP Beads	Vol. of 80% EtOH	Vol. of Low TE in each tube
0-18	1	16*N	16*N	4 *(16*N)	210
19-37	2	(16*N) ÷ 2	(16*N) ÷ 2	4 * ((16*N) ÷ 2)	105
38-56	3	(16*N) ÷ 3	(16*N) ÷ 3	4 * ((16*N) ÷ 3)	70
57-74	4	(16*N) ÷ 4	(16*N) ÷ 4	4 * ((16*N) ÷ 4)	52.5
75-93	5	(16*N) ÷ 5	(16*N) ÷ 5	4 * ((16*N) ÷ 5)	42
94-100	6	(16*N) ÷ 6	(16*N) ÷ 6	4 * ((16*N) ÷ 6)	35

- 13.1 Determine how many 1.5mL Lo-Bind tubes you need based on table above.
- 13.2 Aliquot the pooled library into the tubes equally
- 13.3 Add an equal volume of AMPure XP beads to pooled library amount in a 1:1 ratio
- 13.4 Gentle vortex, quick spin.
- 13.5 Incubate tube(s) at room temperature for 👏 00:10:00

10m

13.6 During the 10 min incubation, prep 80% ethanol in 5mL tube(s) or 15ml conical



```
e.g. to make \perp 5 \text{ mL} of 80%, mix \perp 4 \text{ mL} of 100% ethanol and \perp 1 \text{ mL} of
nuclease-free water
For 100 samples: need \perp 12 mL 100% ethanol, \perp 3 mL water
For 74 samples: need 4 9.6 mL 100% ethanol, 4 2.4 mL water
For 56 samples: need 4 8 mL 100% ethanol, 4 2 mL water
For 37 samples: need 4 5.6 mL 100% ethanol, 4 1.4 mL water
The volumes above take into account the extra ~2mL of 80% ethanol needed in the rest of
the protocol
```

- 13.7 Place tube(s) on DynaMag Spin, and allow the magnetic beads to stick to the side of the tube(s) for (5) 00:05:00
- 13.8 While the tube(s) are still on the magnet, using a p200 pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard supernatant
- 13.9 While tube(s) are still on the magnet, add 2X total volume of 80% ethanol to the beads and incubate for 60 00:00:30

Note

```
e.g. If 48 samples: \Delta 256 \muL of library in 3 tubes: 2X total volume = ( \Delta 256 \muL
pooled library + \Delta 256 \muL beads)*2 = \Delta 1024 \muL 80% ethanol per tube
```

- 13.10 Discard supernatant as in step 13.8
- 13.11 repeat step 13.9 one more time

5m

30s



- 13.12 Discard supernatant as in step 13.8
- 13.13 Leave tube(s) on magnet and allow any residual ethanol to evaporate by letting stand for (c) 00:05:00 uncovered

5m

13.14 Take tube(s) off the magnet

> Add low TE to each tube, refer to table at the beginning of this section for how much Pipette up and down a few times to completely resuspend the bead pellet

13.15 Incubate tubes(s) at room temperature for 00:05:00.

5m

13.16 Place tube(s) in magnetic stand for 00:05:00

5m

Note


THE SUPERNATANT HAS THE DNA IN IT NOW, YOU WANT TO KEEP THE SUPERNATANT

13.17 Using a p200 pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove the supernatant and put into a TPX tube.

Repeat for all library tubes on the magnet. Deposit all supernatants into same TPX tube

13.18 Pipette 4 210 µL of water into another TPX tube

LABEL BOTH TPX TUBES TO AVOID MIXUP!!



14 Sonication



Equipment and supplies needed for this step:

Sonicator (prepped in step 7)

Reagents needed:

TPX tube with library

TPX tube with water as balance

14.1 Sonicate library for 600:04:30 at 20% capacity, and at 4 °C

4m 30s

14.2 Remove tubes after sonication has completed

OK to sit in room temperature for next step.

Discard TPX vial containing water

Turn off sonicator and coolant pump



 \perp 210 μ L of pooled, sonicated library in tube.

Note

DNA fragments should be between 200-500bp long. Can check this on gel or Bioanalyzer here if troubleshooting needed

Some fragments now lack BestRad adapters after sonication. The next step gets rid of these

Day 2 - Bead Selection

15 **Prep buffers for bead binding in Step 17**



Equipment and supplies needed for this step:

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791) p200 pipette and tips

p1000 pipette and tips

heat block

Reagents needed:

2X B+W buffer (directions on how to prepare in step 5)

- CutSmart Buffer 5.0 ml New England Biolabs Catalog #B7204S
- Mainvitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36
- 15.1 Thaw 10X CutSmart Buffer at room temperature
- 15.2 Check heat block to make sure it's on and set to \$\ \bigsep\$ 56 °C
- 15.3 To a 1.5mL Lo-Bind tube add:
 - Δ 25 μL 10X CutSmart Buffer
 - Δ 225 μL water
- 15.4 To a 1.5mL Lo-Bind tube add:
 - Δ 250 μL 2X B+W buffer
 - 🚣 250 μL water
- 15.5 To a 1.5mL Lo-Bind tube add:
 - Δ 175 μL 2X B+W buffer
 - 🚨 175 μL water
- 15.6 Place tube from step 15.5 in heat block set to \$\ \ \ 56 \circ\$



Expected result

1 tube of Δ 250 μL 1X CutSmart Buffer

1 tube of \perp 500 μ L 1X B +W Buffer

1 tube of warm 4 350 μL 1X B +W Buffer in heat block

16 Prep Dynabead M280 streptavidin magnetic beads

Note

Equipment and supplies needed for this step:

p200 pipette and tips

vortex

benchtop centrifuge

DynaMag Spin

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791)

Reagents needed:

2X B+W buffer made in Step 5

Dynabeads M-280 Streptavidin Fisher Scientific Catalog #11-206-D

16.1 Vortex stock bottle of Dynabeads and transfer Δ 20 μL to new 1.5mL Lo-Bind tube

Note


Dynabeads are temperature sensitive, place beads back in fridge immediately after use

16.2 Place tube on the DynaMag Spin for 00:00:10

10s

- 16.3 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the 🚨 20 μL supernatant
- 16.4 Remove tube from magnet and add \perp 100 μ L **2X** B+W buffer





Day 2 - binding DNA to beads

22m

17 Bind DNA to beads

Note

Equipment and supplies needed for this step:

p1000 pipette and tips p200 pipette and tips benchtop centrifuge DynaMag Spin

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791)

Reagents needed:

sonicated library from Step 14.2 prepped beads from 16.9

2X B+W buffer made in Step 5 1X B+W buffer from Step 15.4

warm 1X B+W buffer from Step 15.6



- 17.1 Transfer all $\underline{\underline{A}}$ 210 μL of sonicated library to freshly prepped $\underline{\underline{A}}$ 200 μL Dynabeads tube from Step 16.9
- 17.2 Incubate for 00:20:00 at room temperature while gently vortexing every 2 minutes for 10 seconds (e.g 20:00 minute mark, vortex till 19:50. At 18 minutes, vortex till 17:50, etc)
- 17.3 Quick spin then place tube on the DynaMag Spin for 00:00:30
- 17.4 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the $410 \,\mu$ L supernatant
- 17.5 Take tube off magnet, re-suspend beads in \triangle 150 μ L of **1X** B+W buffer from Step 15.4 Mix by pipetting
- 17.6 Place tube on the DynaMag Spin for 00:00:30

17.7 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the $4 \times 150 \,\mu$ L supernatant

17.8

Repeat steps 17.5 - 17.7 2 more times

17.9 Take tube off magnet, re-suspend beads in Step 15.6 Take tube off magnet, re-suspend beads in Step 15.6

Mix by pipetting

Note

For following steps, keep warm **1X** B+W buffer in heat block at \$\mathbb{g}\$ 56 °C when not in use

20m

30s

30s



17.10 Place tube on the DynaMag Spin for 00:00:30

- 30s
- 17.11 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the \perp 150 μ L supernatant
- 17.12 Repeat steps 17.9 17.11 1 more time
- 17.13 Take tube off magnet, re-suspend beads in \perp 200 μ L of **2X** B+W buffer

Mix by pipetting

Expected result

A tube with Δ 200 μ L of DNA fragments with BestRad adapters bound to the beads.

18 Liberate DNA from beads

Note

Equipment and supplies needed for this step:

thermal cycler
p10 pipette and tips
p200 pipette and tips
p1000 pipette and tips
plate centrifuge
DynaMag Spin
magnet stand for plate

unskirted PCR plate (Fisher Scientific; AB-0700)

manual heated plate sealer

plate seal (Fisher Scientific, AB-0745)

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791)

Reagents needed:

DNA and beads from Step 17.13

1X CutSmart buffer from Step 15.3

Sbfl-HF - 2,500 units New England Biolabs Catalog #R3642L



18.1 Re-suspend beads in \perp 100 μ L of 1X CutSmart buffer from Step 15.3

Mix by pipetting

18.2 Place tube on the DynaMag Spin for 00:00:30

30s

- 18.3 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the $400 \, \mu$ supernatant
- 18.4 Repeat steps 18.1 18.3. Note that this time around the volume of the supernatant to remove will be $400 \, \mu L$
- 18.5 Take tube off magnet, re-suspend beads with bound DNA in Δ 40 μ L 1X CutSmart Mix by pipetting
- 18.6 Transfer all Δ 40 μ L beads mix into well C3 of PCR plate

Note

In our lab, we don't use any PCR tubes, only plates. This step can be done in a tube that would fit a thermalcycler.

- 18.7 Pipette 🚨 2 µL Sbfl-HF (NEB R3642L) into the well
- 18.8 Use a p200 pipette set to 20µl to mix the solution

DO NOT VORTEX or SPIN DOWN PLATE

- 18.9 Heat seal PCR plate
- 18.10 Run thermalcycler program: 3 37hold (General \rightarrow RAD \rightarrow 3 37hold)
 - \$\cdot 37 °C for \cdot 01:00:00

1h



18.11 After program has completed, give the plate a quick spin

18.12 Place plate on magnet in plate format for 00:00:30

30s

Note

In our lab, we rigged a magnet to work on a plate. However, there are commercially available plate magnets. If you performed the incubation in a PCR tube, transfer the liquid to a tube that will fit the tube magnet stand.

18.13 transfer all $\stackrel{\blacksquare}{4}$ 42 μ L supernatant to new 1.5mL Lo-Bind tube

Expected result

 Δ 42 μ L of DNA fragments with BestRad adapters attached in tube

19 Size select DNA fragments

Note

Equipment and supplies needed for this step:

DynaMag Spin unskirted PCR plate (Fisher Scientific; AB-0700) 1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791) **Reagents needed:**

80% ethanol from Step 13.6

p200 pipette and tips

DNA and beads from Step 18.13

Now TE Fisher Scientific Catalog #AAJ75793AP

Ø Agencourt AMPure XP Beckman Coulter Catalog #A63880



19.1 Add 🚨 42 µL of AMPure XP beads to 🚨 42 µL DNA tube Mix by pipetting 19.2 Incubate 4 84 µL mixture at room temperature for 60 00:10:00 10m 19.3 Place tube on the DynaMag Spin for 00:05:00 5m 19.4 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the $\parallel \Delta \parallel$ 84 $\mu \parallel$ supernatant 19.5 While still on magnet, add \perp 168 μ L of 80% ethanol (made the same day) to beads 30s and incubate for 600:00:30 19.6 While the tube is still on the magnet, using a pipette and tip, insert the tip into the tube away from the bead pellet and slowly remove and discard the 🚨 168 👊 supernatant 19.7 Repeat wash steps 19.5 - 19.6 19.8 Leave tube on magnet and allow any residual ethanol to evaporate by letting stand for 5m (5) 00:05:00 uncapped 19.9 Resuspend beads with 4 55 µL low TE Mix by pipetting 19.10 Incubate at room temperature for 00:05:00 5m 19.11 Place on DynaMag Spin for 00:05:00 5m Note NOW THE SUPERNATANT HAS THE DNA IN IT



19.12 Transfer Δ 55 µL supernatant to PCR plate (well C3 is preferred)

Note

If stopping is necessary, transfer supernatant to a new 1.5mL Lo-Bind tube and store in **₽** -20 °C

Expected result

Δ 55 μL purified DNA fragments with BestRad adapters attached

Day 2 - Ligation

20 **Ligating NEBNext adapters to DNA fragment**

Note

Equipment and supplies needed for this step:

thermal cycler p10 pipette and tips p200 pipette and tips manual heated plate sealer plate seals (Fisher Scientific, AB-0745) plate centrifuge vortex

1.5mL Lo-Bind microcentrifuge tube (Fisher Scientific, 13-698-791) timer

Reagents needed:

NEBNext Ultra II DNA Library Prep Kit for Illumina - 96 rxns New England **Biolabs Catalog #**E7645L

2 1M Tris-HCl pH7.5 Fisher Scientific Catalog #AAJ60636K2 purified DNA from Step 19.12



The protocol calls for 10mM Tris-HCl pH 7.5, so make a dilution if needed The NEBNext kit has 2 components: End Prep reagents and Ligation reagents

- 20.1 Thaw reagents from NEBNext End prep kit
 - Ultra II End Prep Enzyme Mix
 - Ultra II EP Buffer

Mix by pipetting (set p200 to 40µl)

- 20.3 Heat seal the plate and quick spin 😝
- 20.4 Run thermalcycler program: 4 End Prep (General \rightarrow RAD \rightarrow End Prep)
 - \$ 20 °C for ♠ 00:30:00
 - \$\cdot 65 °C for \cdot 00:30:00

Set a timer for 45 min

- 20.5 While program is running, drain sonicator.
 - 1. Disconnect the tube labeled "coolant supply" from the sonicator and remove the white nozzle end completely.
 - 2. Turn on the pump, and tip the sonicator backwards to drain as much water as possible, then turn pump off.
 - 3. Use paper towels to absorb all remaining water in sonicator
 - 4. Leave lid open for a day or 2 to air dry. Unscrew the filter, empty out any remaining water and let it sit to air dry.
- 20.6 When timer goes off, thaw the following reagents for ligation:
 - Ultra Prep II Ligation master mix
 - NEBNext Adaptor
 - Ligation Enhancer
- 20.7 Dilute NEBNext adapter 1:10 in a new 1.5mL Lo-Bind tube

$$\bot$$
 1 μL of NEBNext Adaptor + \bot 9 μL 10mM Tris HCl pH 7.5

1h



The manufacturer's protocol says to dilute the NEBnext adapter IF the DNA input is ≤ 100 ng. We don't check this and assume the adapter always needs diluting. Also, the protocol says to mix the 10mM Tris HCl pH 7.5 with 10mM NaCl, but we haven't been doing this (May, 2023).

20.8 Make Ligation master mix in a 1.5mL Lo-Bind tube by combining reagents below. Multiply by the number of samples you have + 5-10 extra to account for pipette error.

Adaptor ligation	Vol. per sample			
Ultra Prep II Ligation master mix	15 μΙ			
*NEBNext Adaptor (diluted)	2.5 μΙ			
Ligation Enhancer	1 μΙ			
Total	18.5 µl			

Note

The NEBNext adapter is a double stranded hairpin that will ligate to both sides of the DNA fragments. It contains the priming sites for the PCR primers in Step _

- 20.9 Add master mix to well of plate with DNA (the volume should be 483 µL) Mix by pipetting (set p200 to 40µl and mix)
- 20.10 Heat seal the plate, quick spin
- 20.11 Run thermal cycler program: 5 Adaptor Ligation (General → RAD → 5 Adapter Ligation \$ 20 °C for (5) 00:15:00

20.12 After program has ended add 🚨 3 μL of NEBNext USER enzyme

Mix by pipetting as above

15m



The NEBNext adapter has a uracil base in it. The USER enzyme cleaves the hairpin adapter at the uracil so now the DNA fragments are linear and have the PCR priming sites on both sides needed in Step 22

- 20.13 Heat seal the plate, quick spin
- 20.14 Run thermalcycler program: 6 USER enzyme (General \rightarrow RAD \rightarrow 6 USER enzyme) \$ 37 °C for (5) 00:15:00

15m



A well in a PCR plate with \perp 86 μ L of DNA fragments that have BestRad adapters on one side with the barcodes and NEBNext adapters on both sides with the PCR priming sites

Day 2 - Ampure Size Selection and Clean-Up

20m

21 Size selection and cleanup of adaptor-ligated DNA

Note

Equipment and supplies needed for this step:

p10 pipette and tips

p200 pipette and tips

1.5mL Lo-Bind microcentrifuge tubes (Fisher Scientific, 13-698-791)

DynaMag Spin

unskirted 96 well PCR tray (Fisher Scientific; AB-0700)

Reagents needed:

- Invitrogen Nuclease-Free water **Fisher Scientific Catalog #**43-879-36
- Agencourt AMPure XP **Beckman Coulter Catalog** #A63880
- low TE Fisher Scientific Catalog #AAJ75793AP

80% ethanol



- 21.3 Incubate at room temperature for 00:05:00
- Place tube on DynaMag Spin, and allow the magnetic beads to stick to the side of the tube for 00:05:00
- 21.5 Transfer all \perp 143.5 μ L supernatant to new 1.5mL Lo-Bind tube
- 21.6 Add $\underline{\underline{}}$ 26 μL AMPure XP beads

mix by pipetting

21.7 Incubate at room temperature for (5) 00:05:00

21.8 Place tube on DynaMag Spin, and allow the magnetic beads to stick to the side of the

- tube for 👏 00:05:00
- 21.9 Remove and discard supernatant (~ Δ 170 μL)
- 21.10 While still on magnet, add \triangle 200 μ L of fresh 80% ethanol to beads and incubate for 00:00:30
- 21.11 Repeat 21.9 21.10
- 21.12 Leave tube on magnet and allow any residual ethanol to evaporate by letting stand for 00:05:00 uncovered


5m

30s


5m

5m





22 PCR amplification



22.1 Thaw reagents from NEBNext Multiplex Oligo kit (orange box) labeled: "RADseq PCR Amp"

22.2 Select an Index (i7) from the NEBNext Multiplex Oligo kit to use and record it in the "Final Library" tab in *BestRAD_library_prep_labData v1.2.xlsx*.

Note

The i7 MUST BE unique per library. If you are preparing more that 1 library, use different i7s.

22.3 Add the following reagents to the $\ \ \underline{\ \ }$ 15 μL of DNA:

PCR Amplification	Volume		
Ultra II Q5 Hot Start HiFi PCR Master Mix	25 µl		
Index Primer/i7 Primer	5 µl		
Universal PCR Primer/i5 Primer	5 µl		

Total volume is Δ 50 μL

- 22.4 Heat seal the plate, gently vortex, quick spin
- 22.5 Run thermal cycler program: General \rightarrow RAD \rightarrow 7 NEBNext PCR 11x

Runs for about 00:28:00

Cycle Step	Temp	Time	Cycles	
Initial Denaturation	98°C	30 s	1	
Denaturation	98°C	10 s	122	
Annealing/Extension	65°C	75 s	11	
Final Extension	65°C	5 min	1	
Hold	4°C	- 00	1	

Note

The BestRAD protocol says to perform 4-12 cycles of PCR. Currently we use 11.

28m




22.6 Transfer Δ 50 μL PCR product to new 1.5mL Lo-Bind tube

Note

Do NOT stop here. Complete PCR purification below

23 PCR purification



- Add Δ 45 μL AMPure XP beads to the Δ 50 μL of DNA from Step 22.6 mix by pipetting
 Incubate at room temperature for 00:05:00
- Place tube on DynaMag Spin, and allow the magnetic beads to stick to the side of the tube for 00:05:00
- 23.4 Remove and discard supernatant (~ Δ 95 μL)



23.5 While still on magnet, add 🚨 200 μL of fresh 80% ethanol to beads and incubate for **(5)** 00:00:30 23.6 Repeat 21.9 - 21.10 23.7 Leave tube on magnet and allow any residual ethanol to evaporate by letting stand for (5) 00:05:00 uncovered 23.8 Take the tube off the magnet and re-suspend the beads in 4 33 µL low TE Mix by pipetting 23.9 Incubate for 00:02:00 at room temperature 23.10 Put tube on DynaMag Spin for 00:05:00 23.11 Transfer 4 28 µL into a new 1.5mL Lo-Bind tube Label the tube with species "Radseq", library number, date Store at
\$\mathbb{4} -20 °C if stopping for the day **Expected result** Δ 28 μL of final, undiluted library

Day 3 - quantify library

1h 3m

24 **Qubit Rad library**



Equipment and supplies needed for this step:

Qubit fluorometer Qubit assay tubes (Fisher Scientific, Q32856) 1.5mL or 2.0mL microcentrifuge tube p10 pipette and tips p200 pipette and tips

vortex

benchtop centrifuge

Reagents needed:

DNA library from Step 23.11

Qubit dsDNA HS assay kit Fisher Scientific Catalog #Q32854

- 24.1 Thaw library to room temperature
- 24.2 Make 1:200 diluted dye mix in a microcentrifuge tube

↓ 199 μL Qubit buffer ↓ 1 μL Qubit dye

24.3 Pipette \perp 198 μ L of diluted dye mix from Step 24.2 and \perp 2 μ L of the library into a qubit assay tube

Vortex thoroughly, quick spin

24.4 Incubate at room temperature for (5) 00:03:00

3m

24.5 Put qubit assay tube in Qubit Fluorometer and take 2 readings per sample and enter them into BestRAD_library_prep_labData v1.1 spreadsheet in "Final Library" tab. The spreadsheet has a column that will take the average. The target is 4-20 ng/µl.



Don't use any tube other than the qubit assay tubes in the Qubit.


Make sure you set the sample volume and the output to be read in $ng/\mu l$, not ng/mL

Expected result

 \triangle 26 μL of your Radseq library at a concentration >= **4 ng/μl**.

25 E-gel of library





25.1 Thaw 1 Kb ladder, loading dye, and RAD library if needed

Gentle vortex, quick spin

- 25.2 Cut strip of parafilm, tape to benchtop if necessary
- 25.3 Onto parafilm, pipette \perp 15 μ L of 1 Kb ladder. Add \perp 5 μ L of nuclease free water to the ladder. Mix by pipetting



25.4 Onto parafilm, pipette \perp 5 μ L of the RAD library. Add \perp 3 μ L of loading dye and Δ 12 μL of nuclease free water to the library

RETURN LIBRARY TO FREEZER

- 25.5 Unwrap e-gel, remove comb, remove camera (top) from the Power Snap unit and load egel into machine
- 25.6 Pipette all 4 20 µL of diluted ladder from parafilm into lane 1 of e-gel
- 25.7 Pipette all A 20 µL of diluted RAD library from parafilm into lane 2 of e-gel
- 25.8 Pipette Δ 20 μL of water into remaining empty wells of e-gel
- 25.9 Close lid, turn Power Snap on, select "E-gel 0.8-2%" program, push start, put camera back on Power Snap unit

Set timer for 15 min


Note

The program's duration is 26 min but the gel only needs to be run for 10-15min


- 25.10 When timer goes off, tap "view gel" → adjust brightness to max → Take picture → Tap "Export" → insert thumb drive into usb port → Tap "Export"
- 25.11 Paste e-gel picture in the BestRAD_library_prep_labData v1.1 spreadsheet under "Final Library" tab

Crop picture to include only lane1 and lane2

Use this image to determine approximate size of library and enter it into "Final Library" tab



size of library should be ~300-500bp



26 **qPCR library quantification**



Equipment needed

p10 pipette and tips p200 pipette and tips repeater M-4 pipette

1.5mL Lo-Bind microcentrifuge tubes (Fisher Scientific, 13-698-791)

Applied Biosystems MicroAmp Optical 7500 Plate (Fisher Scientific 43-168-13)

MicroAmp Optical Adhesive Film (Fisher Scientific 43-119-71)

vortex

benchtop centrifuge

0.1 mL Eppendorf combi-tips (Fisher Scientific; 13-683-700)

KimWipes

thumb drive

Equipment

7500 Real Time PCR system

NAME

qPCR

TYPE

Thermo Scientific

BRAND

discontinued

SKU

Reagents needed

- NEBNext Library Quant Kit for Illumina 500 rxns New England Biolabs Catalog #E7630L
- Minvitrogen Nuclease-Free water Fisher Scientific Catalog #43-879-36
- TE pH 8.0 (1X TE Solution) **IDT Technologies Catalog #11-01-02-05**

26.1 Thaw reagents in NEBNext library quantification kit and the RAD library

Gently vortex the library, quick spin



26.2 Check small bottle in fridge for 1X dilution of qPCR dilution buffer from NEBNext Library Quant kit. If low, make more:

🕹 150 μL 10X Buffer

Δ 1350 μL nuclease-free water



26.3 Make sure the master mix in the qPCR kit has been prepped. Look for the "Prep checkbox". If it hasn't, follow instructions below

1.2. Prepare NEBNext Library Quant Master Mix (with primers)

NEB #E7630S: Add 100 µl NEBNext Library Quant Primer Mix to the tube of NEBNext Library Quant Master Mix (1.5 ml). Mix by vortexing for 10 seconds. Write the date on the master mix tube to indicate that primer mix has been added.

NEB #E7630L: Add 500 µl NEBNext Library Quant Primer Mix to the bottle of NEBNext Library Quant Master Mix (7.5 ml). Mix by vortexing for 10 seconds. Write the date on the master mix bottle to indicate that primer mix has been added.

Note: If using ROX for normalization, add ROX to the NEBNext Library Quant Master Mix. Add 20 µl (NEB #E7630S) or 100 µl (NEB #E7630L) ROX and vortex. See Table 1 in Section 2 for ROX selection.

26.4 Make serial dilutions in 1.5 Lo-Bind tubes as follows:

Note

The accuracy of the qPCR depends on the accuracy and precision of the preparation of the serial dilutions below. In the following steps, each time a dilution is made, it should be thoroughly mixed but you also want to avoid bubbles forming. Annoying, yes.

Note

If the Qubit score from Step 24.5 was > 25 ng/ μ l, use dilutions starting from 1:100,000 or 1:200,000 in qPCR

If the Qubit score was > 30 ng/µl, make addition dilutions up to 1:2,000,000 dilution

If the Qubit score was high, you may not want to run the more concentrated dilutions (1 and 2).

- 1. 1:1000 dilution: Δ 1 μ L RAD library + Δ 999 μ L 1X Buffer
- 2. 1:10,000 dilution: \bot 10 μ L of 1:1000 dilution + \bot 90 μ L 1X Buffer
- 3. 1:100,000 dilution: $\underline{\mbox{\sc L}}$ 10 $\mu\mbox{\sc L}$ of 1:10,000 dilution + $\underline{\mbox{\sc L}}$ 90 $\mu\mbox{\sc L}$ 1X Buffer
- 4. 1:200,000 dilution: Δ 50 μ L of 1:100,000 dilution + Δ 50 μ L 1X Buffer
- 6. 1:1,000,000 dilution: $\underline{\underline{A}}$ 10 $\mu \underline{L}$ of 1:100,000 dilution + $\underline{\underline{A}}$ 90 $\mu \underline{L}$ 1X Buffer
- 7. 1:2,000,000 dilution: $\underline{\bot}$ 50 μL of 1:1,000,000 dilution + $\underline{\bot}$ 50 μL 1X Buffer



	1	2	3	- 4	- 5	6	7	8	9	10	11	12
A	STANDARD 2 10pM	STANDARD 2 10pM	STANDARD 2 10pM	LIBRARY 1:100,000	LIBRARY 1:100,000	LIBRARY 1:100,000						
В	STANDARD 3	STANDARD 3 1pM	STANDARD 3 1pM	LIBRARY 1:200,000	LIBRARY 1:200,000	LIBRARY 1:200,000						
c	STANDARD 4 0.1pM	STANDARD 4 0.1pM	STANDARD 4 0.1pM	LIBRASIY 1:500,000	LIBRARY 1:500,000	LIBRARY 1:500,000						
D	STANDARD 5 0.01pM	STANDARD 5 0.01pM	STANDARD 5 0.01pM	LIBRARY 1:1,000,000	LIBRARY 1:1,000,000	LIBRARY 1:1,000,000						
E	NTC	NTC	NTC									
F												
G												
н												

Example of plate layout - all standards and library dilutions are in triplicate

26.6 Using a p10 pipette, add $\stackrel{\text{$\square$}}{=}$ 2 $_{\mu L}$ of standards (use standards 2-5) from qPCR kit, library dilutions and no template control (NTC; 1X Buffer) to each well according to plate layout. Discard tip after each use.

Note

Add to side of well. To avoid bubbles, DO NOT pump!

26.7 Use MicroAmp Optical Adhesive Film to seal plate (this isn't heat sealed). Use a KimWipe to press down on the seal to avoid smudges

Quick spin **B DO NOT VORTEX**

Check underside of plate for dirt/lint, check wells for bubbles

- 26.8 Turn on 7500 qPCR system. Push to eject tray, load plate, push to close
- 26.9 1. Open 7500 v2.3 software
 - 2. User name: GUEST
 - 3. Click "Template" → Open "E7630 slow_v2" (screen will be on "Plate setup")
- 26.10 1. In left "Setup" column: Click on "Experiment Properties"




- 2. Enter in Experiment name (e.g. CarpRADseq_library1_09292022)
- 3. Leave everything else as is
- 26.11 1. In left "Setup" column: Click on "Plate Setup"
 - 2. Click on "Assign Targets and Samples" tab
 - 3. Highlight all wells that you're not using a Right click and select "Clear"
- 26.12 Click on "Save as" and save to project folder on computer Start Run (~) 01:00:00)

1h

26.13 Export data.xls onto flash drive, and enter in results on NEB website to get the undiluted concentration of your RAD library

https://nebiocalculator.neb.com/#!/qPCRlibQnt



Example of Nebiocalculator – Efficiency needs to be as close to 100% as possible (+/- 10%), and $R^2 = 1$

Note

Fragment size is the approximate band size obtained from E-gel

26.14 Dilute RAD library to 4nM using V1 = (C2*V2/V1) using 1X TE



Note C1 = concentration of library according to qPCR C2 = 4nMV2 = Δ 50 µL (This may change if you have a very high or low concentrated library. 50µl is appropriate if your library was ~150-250nM) Keep original and diluted libraries at <a>♣ -20 °C **Expected result**

Refer to "BESTRAD on NextSeq Protocol" for how to load the library

Citations

Omar A. Ali, Sean M. O'Rourke, Stephen J. Amish, Mariah H. Meek, Gordon Luikart, Carson Jeffres and Michael R. Miller. RAD Capture (Rapture): Flexible and Efficient Sequence-Based Genotyping https://doi.org/10.1534/genetics.115.183665